

400-6688-9966

与世界同步，创国际一流，铸就生命空间的辉煌。



[Http://www.cnnthp.cn/](http://www.cnnthp.cn/)

Twoway 通卫®

南通华普工贸有限公司

地址：江苏省南通市狼山镇84-3号
电话：0086-513-85720986 0086-513-88118101
传真：0086-513-85721396 0086-513-85714586
E-mail：nthp@msn.com
Http：//nthuapu.en.alibaba.com/

Nantong Hope Industry&Trade Co.,ltd

Add:NO.84-3 Langshan Town,Nantong City,
Jiangsu Province,P.R.China.
Tel:0086-513-85720986 0086-513-88118101
Fax:0086-513-85721396 0086-513-85714586
E-mail:nthp@msn.com
Http://nthuapu.en.alibaba.com/

Twoway 通卫®

Nantong Hope Industry&Trade Co.,ltd
南通华普工贸有限公司



Production equipment / 生产设备



Company introduction / 公司简介

Nantong Hope Industry and Trade Co., Ltd. was founded in August of 2003 and located in Nanong, one of the first groups of coastal open cities in China. With registered capital of RMB 500,000 yuan and four years of continuous effort, it has set up a complete production system from weaving to finishing. Now it has two weaving branches, more than 100 weaving machines and nearly 100 employees. The products are mainly sold to hospitals, medical care institutes and foreign trade companies. Since the acquirement of self-management of export at the end of 2004, it is gradually exploring the international markets.

The company attained the quality system certifications of ISO9001 and ISO13485 in 2004, which better guarantees the quality of products.

The company mainly produces and sells the dressing of secondary sanitary materials, including medical absorbent gauze, cotton tampon, gauze piece, belly belt, bandage, non-woven fabrics, cotton stick, cotton printing cloth, embroidered cleaning cloth, mouth protector and handkerchief, etc.

The operation principle for the company is to provide best products at most competitive price to the customers. With the business conception of "credit brings market; quality pushes development", the company will exert itself to provide the best service and most satisfactory products to customers.

南通华普工贸有限公司位于全国第一批沿海开放城市之一南通，公司成立于2003年8月，注册资金50万元。经过近4年来的不懈努力，公司现在已经发展拥有从织造到成品的完整生产体系。拥有2家织造分厂，100多台织机，员工近100人。产品主要由原来销售到各个大中小型医院、医疗机构、各外贸公司。于2004年底获得自营出口权，之后逐步开拓出面向国际市场的销售。

公司在2004年获得ISO90001和ISO13485质量体系认证，为产品的质量提供了良好的保证体系。

公司主要生产销售的产品为二类卫生材料敷料，包括各种规格的医用脱脂纱布、棉球、纱布片、腹部巾、纱布绷带、无纺布片、棉签等，以及各类全棉类印花、绣花抹布、口罩、手帕等。

公司经营宗旨是以最优的性价比为最精明的客户服务，相信在我们检查的“诚信创造销售，质量谋求发展”的经营理念下，一定会为客户提供优质的服务和满意放心的产品。



Absorbent Gauze / Zigzag Gauze 脱脂纱布

Material : 100% Cotton , degreased and bleached

Cotton Yarn : 40s, 32s, 21s

Mesh : 12x8,19x9, 20x12,19x15,24x20,28x24 or as per your request

Size : 44"x20yds, 65cmx100m, 90cm(36")x100y, 36"x50y, 90cmx50m, special size pls contact us

Layer : 1ply, 2ply, 4ply

Type : with x-ray or without can be made

Color : white

Packing : individual poly bag or paper bag + woven bag or outer cartons

OEM : customer's design is welcome

Application : hospital, clinic, first aid, other wound dressing or care

Absorbent Gauze Roll/ 纱布卷

Material : 100% Cotton , degreased and bleached

Cotton Yarn : 40s, 32s, 21s

Mesh : 12x8,19x9, 20x12,19x15,24x20,28x24 or as per your request

Size(Width) : 90cm(36"), 120cm, 130cm, special size pls contact us

Size(Length) : 10m, 25m, 50m, 90m,100yds(91m), 1000m, 2000m, 3000m or as per your request

Layer : 1ply, 2ply, 4ply

Type : with x-ray or without can be made

Color : white(mostly), green, blue

Packing : non-sterile, 1roll/blue paper or white paper packing, a label, a number of rolls/carton


OEM : customer's design is welcome

Application : hospital, clinic, first aid, other wound dressing or care






Gauze pad
Gauze sponge



Vaseline Gauze



Sterile gauze swab / 消毒纱布片

Material : 100% cotton, degreased and bleached

Typ : esfolded or unfolded edge, with or without detectable ray

Yarns : 40s, 32s, 21s

Mesh sizes : 12x8,19x9, 20x12,19x15,24x20,28x24 or as per your request

Sizes : 5x5(2"), 7.5x7.5, 10x10, or customized

Layers : 4plys, 8plys, 12plys, 16plys, or customized

Sterilization : EO or Gamma ray

Packing : *1pc, 2pcs, 5pcs, 10pcs per paper-poly pouch or paper-paper pouch 10pcs~100pcs per mid boxor can be as your request"

Application : hospital, clinic, first aid, other wound dressing or care

non-sterile gauze swab/ 不消毒纱布片

Material : 100% cotton, degreased and bleached

Types : folded or unfolded edge, with or without detectable ray

Yarns : 40s, 32s, 21s

Mesh sizes : 12x8,19x9, 20x12,19x15,24x20,28x24 or as per your request

Sizes : 5x5(2"), 7.5x7.5, 10x10, or customized

Layers : 4plys, 8plys, 12plys, 16plys, or customized

Packing : *50pcs, 100pcs, 200pcs per paper pack or poly bagor can be as your request"

Application : hospital, clinic, first aid, other wound dressing or care



Material : 100% cotton, degreased and bleached
Types : sterile or non-sterile, pre-washed or non-washed
Yarns : 40s, 32s, 21s
Mesh sizes(cm) : 19*15, 20*12, 24*20, 26*18, 28*24, 30*20 or customized
Sizes(cm) : 10*40, 20*40, 10*45, 30*30, 37*45, 40*40, 45*45, or customized
Layers : 4ply, 6plys, 8plys, 12plys, 16plys, or customized
Sterilization : EO or Gamma ray
Packing : "non-sterile -poly bag packing
 sterile - 1pc/pack, 2pcs/pack, 5pcs/pack, 10pcs/pack
 or can be as your request"
Appication : hospital, clinic, first aid, other wound dressing or care



Gauze lap sponge/ 纱布巾
 Abdominal swab/ 腹部垫

Gauze lap sponge
 Abdominal swab





Triangular Bandage/ 三角绷带



Gauze Bandage 纱布绷带

Material : 100% Cotton , degreased and bleached

Cotton Yarn : 40s, 32s, 21s

Mesh : 12x8,19x9, 20x12,19x15,24x20,28x24 or as per your request

Size(Width) : 5cm,7.5cm,10cm,15cm, 20cm, special size pls contact us

Size(Length) : 3y, 3m, 4m, 4.5m, 5m, 5y, 10y, 10m or as per your request

Layer : 1ply

Type : Head bandage or normal one

Color : white

Packing : *Mostly 12rolls/dozen, 80~100dozens/carton or single packing

Paper packing, PE bag packing, mid-box packing can be provided*

OEM : customer's design is welcome

Application : hospital, clinic, first aid, other wound dressing or care

Material : 100% Cotton or Non-woven(viscose and fiber)

Cotton Yarn : 40s

Mesh : 50x48, 50x30, 44x30, 40x30, or as per your request

Size : 36"x36"x50", 40"x40"x56, special size pls contact us

Color : raw white or white

Packing : 1pc/waterproof bag, 240pcs/ctn or as per your request

OEM : customer's design is welcome

Application : Ideal for splinting, outer layer wound bandaging, first aid etc.



Absorbent Gauze Ball/ 纱布球

Material : 100% Cotton , degreased and bleached

Cotton Yarn : 40s, 32s, 21s

Mesh : 12x8,19x9, 20x12,19x15,24x20,28x24 or as per your request

Open Size : 15x16cm, 20x20cm, 24x25cm, 30x30cm, 40x40cm, special size pls contact us

Type : with x-ray or without can be made

Color : white

Packing : sterile pouch packing or non-sterile poly bag packing

OEM : customer's design is welcome

Application : hospital, clinic, first aid, other wound dressing or care

Features and advantages

- 1.100% cotton
- 2.Highly soft,Good absorbency
- 3.Soft,feel comfortable
- 4.Used to absorbing blood and other liquid
- 5.OEM service available
- 6.Application for personal cleaning care, remove cosmetic, making up or other cleaning use



Asorbent cotton roll/ 脱脂棉卷

Material : 100% cotton

Color : white

Type : with and without paper interleaved

Weight : 25g, 50g, 100g, 150g, 250g,454g, 500g,1000g

Packing : individual paper or polybag packing



Cotton Swab/ Cotton Bud 棉签

Material : 100% Cotton , degreased and bleached

Type : single head or double head

Stick : wooden, plastic, bamboo, paper or as per your request

Size : 1- dia of stick:2.5mm-4.0mm

2- dia of head:4-15mm

3- length of the stick:75mm-400mm

4- length of the cotton tip:15-50mm

5- also can be made as per your request

Packing : sterile or non-sterile

OEM : customer's design is welcome

Application : hospital, clinic, first aid, other wound dressing or care

Absorbent Cotton Ball/ 棉球

Material : 100% Cotton, degreased and bleached

Specification : 0.1g, 0.2g, 0.3g, 0.5g, 1g, 2g, other size pls contact us

Color : White or Colored

Type : sterile or nonsterile

Packing : (1) sterile: 5pcs or 10pcs/pouch

(2) non-sterile: 100pcs~1000pcs/poly bag or as per your request

OEM : customer's design is welcome

Application : hospital, clinic, first aid, other wound dressing or care

Absorbent Gauze Face Mask/ 纱布口罩



Material : 100% Cotton , degreased and bleached

Cotton Yarn : 40s

Mesh : 50x50, 30x30mesh or as per your request

Size : 17.5x9.5cm,15x20cm,16.5x9cm,9.5x13.5cm,10.5x8cm,17.5x12cm,11.5x8.5cm, special size pls contact us

Type : ear-loop or tie on

Color : white or printed color

Packing : individual polybag packing or 10~50pcs polybag packing or as per your request

OEM : customer's design is welcome

Application : medical & health, surgery room, laboratory, industry, household etc.

Cleaning cloth 清洁布

Stitch : Pearl stitch

Shape : square or as your request


Size : 20cm x 20cm, 25cm x 25cm, 30cm x 30cm
or as per your request

Fabric style : plain or printed

Packing : 2pc/pack , OPP bag

Features and advantages

- 1.100% cotton gauze fabric
- 2.Breathability and good absorbency
- 3.Different dimension and prints
- 4.OEM service available
- 5.Easy wash,quick dry
- 6.Application for household cleaning etc.



Material : 100% cotton gauze

Stitch : Pearl stitch

Shape : square

Size : 20cm x 20cm, 25cm x 25cm, 30cm x 30cm or as per your request

Fabric style : plain or printed

Packing : 1pc/pack ,100pc/pack PE bag

Handkerchief/ 手帕

