



公司官微



**信阳颐和非织布有限责任公司**

XINYANG YIHE NON-WOVEN CO.,LTD.



**信阳颐和非织布有限责任公司**  
XINYANG YIHE NON-WOVEN CO.,LTD.

厂址：河南省信阳市息县工业集聚区海虹大道1号

ADD: NO.1 Haihong Road, Xixian Industrial Park, Xinyang, Henan, China.

电话：86-371-5857018 传真：86-376-5857019

销售处：河南省郑州市金水区九如路荣盛华府15幢2单元102

Sales Office: #15,2-102, Rongshenghuafu, Jiuru Road, Jinshui District, Zhengzhou, Henan, China

电话：86-371-63247177/63731999

网址：www.yhnonwovenfabric.com



产品研发  
Product Development



售后服务  
After-sale Services



无纺布基材  
Nonwoven Base Material

## 一站式服务 One Stop Services



销售  
Sales



无纺布后整理  
Off-line/On-line Treatment



成品加工  
Finished Product Processing

### 甄选优质原材料体验完美品质

SELECT SUPERIOR QUALITY RAW MATERIAL AND EXPERIENCE PERFECT QUALITY

信阳颐和拥有一体化的供应链，能够更好地掌握每个环节的成本并进行质量控制，生产出更具竞争力的产品。

With a complete and organized supply chain, we are dedicated to providing total Nonwoven solutions to customers, We are positioned to be a **one-stop service provider**.





信阳颐和非织布有限责任公司  
XINYANG YIHE NON-WOVEN CO.,LTD.

## 公司介绍 COMPANY PROFILE

信阳颐和非织布有限责任公司成立于2008年，旗下有分公司河南颐和健康用品有限公司，合计员工300人。每年产出超过10万吨无纺布及大量的无纺布制品。

公司已取得ISO9001、ISO13485质量管理体系，CE和FDA注册，国内医疗器械生产许可证和注册，是一家集无纺布原料生产、工艺研发、制品加工和自营出口为一体的医用和个人防护用品综合服务公司。

信阳颐和公司纺粘无纺布年产能达十万多吨，广泛应用于三大领域：

卫生材料用：亲水、拒水、超柔软无纺布

医疗耗材用：三抗高透气无纺布、抗病毒复合无纺布

工业用：抗老化、防静电、阻燃抗菌等彩色无纺布

我们拥有符合GMP标准的十万级净化车间，配备有专门的理化实验室和EO灭菌设备，可以为客户提供高品质的隔离衣、手术衣、防护服、实验服、口罩、鞋套、条帽、床单、手术洞巾和手术包等一次性医疗防护用品。秉承着“稳健运作、诚信经营、成果共享”的创业理念，信阳颐和公司旨在为客户提供低成本、高品质、高效率的专业化服务。

Xinyang Yihe Non-woven Co.,Ltd was established in 2008,has subsidiary Henan Yihe Health Products Co.,Ltd.With over 300 skilled employees,producing over 100,000 tons nonwoven fabrics and numerous finished goods.

Our company has obtained ISO9001,ISO13485,CE and FDA certificates and also obtained domestic medical device production license and registration,which is one stop supply platform of non-woven raw material production, process research and development,product processing and self export as a whole.

Xinyang Yihe Non-woven currently produces more than 100,000 tons of nonwoven fabrics on an annual basis that are widely used in three industries:

Hygiene market: Nonwoven products for hygiene market include hydrophilic, hydrophobic,super-soft,silky-soft nonwovens.

Medical market: We provide SSMMS/SMMS nonwoven fabric with excellent barrier properties for the medical market,in addition,we offer post-treatments such as alcohol and blood repellency as well as antibacterial laminated composites to meet the market demands for performance,protection,comfort and competitiveness.

Industrial market: We have a wide range of color selections for protective and industrial applications,other post-processing methods include antistatic,flame retardant, antimicrobial,and anti-aging treatments.

We have equipped with 100,000 grade cleanroom which is complied with the requirement of good manufacturing Practice(GMP),Coupled with high standard laboratory and an ethylene oxide(EO) sterilization chamber.

Our company produces and delivers multiple categorie of clinically-superior products including isolation gown,surgical gown,coveralls, lab coat,face mask,shoe cover,bouffant cap,bedsheet,surgical drape and surgical pack etc.Following the spirit of "sound operation, honest management and sharing of results", these values have been the essential elements of our corporate culture and of our very existence.Our goal is to provide our partner the most efficient,cost-effective and high quality infection prevention solutions.




信阳颐和非织布有限责任公司  
XINYANG YIHE NON-WOVEN CO.,LTD.



信阳颐和非织布有限责任公司  
XINYANG YIHE NON-WOVEN CO.,LTD.

# COMPANY CERTIFICATES 公司认证





## 无纺布生产线 PRODUCTION LINE

公司拥有3条高速纺熔生产线及其无纺布后处理设备，采用100%聚丙烯制造的纺粘无纺布，年产能可达十万多吨。

Our company has 3 high-speed spunbond production lines and its non-woven post-processing equipment. Spunbond non-woven fabrics made from 100% polypropylene can be produced at an annual capacity of more than 100,000 tons.

SS纺粘无纺布(幅宽2.0M)  
Spunbond Non woven Fabric (WIDTH: 2.0M)

SS纺粘无纺布(幅宽2.4M)  
SS Non woven Fabric ( WIDTH: 2.4M )


SSS 纺粘无纺布 ( 幅宽3.2M )  
SSS Non woven Fabric( WIDTH: 3.2M )

熔喷无纺布Melt-Blown (幅宽3.2M)

SMS(纺粘/熔喷/纺粘) ( 幅宽3.2M)  
SMS Non woven Fabric(WIDTH: 3.2M )

SSMMS(纺粘/纺粘/熔喷/熔喷/纺粘) ( 幅宽3.2M)  
SSMMS Non woven Fabric ( WIDTH: 3.2M )

无纺布后处理设备 ( 幅宽3.2米 )  
Non-woven reprocessing equipment ( WIDTH: 3.2M)



## SMS/SMMS/SSMMS 无纺布 SMS/SMMS/SSMMS Nonwoven Fabric

### 布料用途：

公司生产的SMS/SMMS/SSMMS 无纺布可用于制作手术洞巾、手术衣、洁净服等医用产品，适用于蒸汽或EO(环氧乙烷)灭菌，并与产品标签和灭菌指示标签相兼容。

### EN13795/YY0560

公司生产的SMS/SMMS/SSMMS 医用无纺布符合EN13795:2019 和YY0506:2019标准的性能要求。

### 生物相容性

公司生产的SMS/SMMS/SSMMS 无纺布通过生物相容性测试，测试结果显示，产品无毒、无致敏性、不刺激皮肤。

### Application:

SMS/SMMS/SSMMS nonwoven fabric has been successfully used in disposable surgical drapes, surgical gowns and clean air suits applications. Steam and EO sterilization are applicable to the SMS/SMMS/SSMMS fabric which is also compatible with product labels and sterilization indicators.


### EN13795/YY0560

Our SMS/SMMS/SSMMS nonwoven fabric complies with N13795:2019 and YY0506:2019 standards for disposable surgical drapes, surgical gowns and clean air suits.

### Biocompatibility

SMS/SMMS/SSMMS nonwoven fabric has been thoroughly evaluated and tested for biocompatibility and the fabric shows no toxicity and no skin irritation and doesn't cause any allergic reaction.





### 纺粘无纺布

#### Spunbond Nonwoven Fabric

公司生产的纺粘无纺布，拥有单层或多层的组合结构，透气性好、强度高、耐磨损、耐老化，且质感柔软。此外还可以通过功能整理，赋予纺粘无纺布亲水、抗静电、抗菌等特殊功能。广泛应用于口罩、医用吸水面层、纸尿裤、卫生巾、擦拭布、无纺布袋等。

Spunbound fabric with multiple or single layer has characteristics of excellent breathability, brilliant mechanical properties, abrasion resistance and UV resistance. In addition spunbond nonwovens can be given a special feature such as hydrophilic, anti-static and anti-bacterial through functional finishing. It's widely used for disposable face mask, surgical drape absorbent top sheet, baby diaper, sanitary napkin, wipes and shopping bag etc.







Baby Diaper



Medical gown/mask



Nonwoven bag



### 无纺布后处理

#### Nonwoven Post-processing Treatment

#### 三抗+抗静电处理

我们可以对SMS/SMMS/SSMMS无纺布进行三抗和抗静电处理，处理后的无纺布具有抗酒精、抗血液、抗油、抗静电的功能。处理后的无纺布符合EN13795:2019和YY/T0506:2019标准的要求。可应用于一次性手术衣及手术铺单。

#### 抗静电处理

主要应用于工业、医疗、防护用品。

#### 吸水处理

主要应用于医用耗材面层吸水材料，如手术床单，一次性湿巾上。

#### 阻燃处理

广泛应用于家具产品和航空产品。

#### 抗紫外线处理

主要应用于农业覆盖布，汽车罩等面料的防晒，抗老化用。

#### Repellent Treatment

Alcohol repellent and anti-static treatments can be applied to the SMS/SMMS/SSMMS non woven fabric simultaneously. The treated fabric features excellent alcohol, blood and oil repellent properties, and it complies with the performance requirements of EN13795:2019 and YY/T0506:2019. The treated fabric has been widely used for disposable surgical gowns and drapes.

#### Anti-static treatment

Mainly used for industrial, medical and protective products.

#### Absorbent treatment

Applied to the top sheet for disposable medical products like surgical drapes, wet wipes.

#### Flame Retardant Treatment

Applied to home furnishing and aviation non woven products.

#### UV Resistance treatment

Primarily used for agriculture & landscaping, land cover, car cover and packaging materials.





## 熔喷无纺布

Meltblown Nonwoven Fabric

公司可以生产不同防护的高效低阻的熔喷无纺布，进口熔喷头及先进的生产工艺，使我司熔喷布可以分别满足国内外不同标准对防护、医疗和卫生等口罩的要求。我司熔喷布性能稳定，带电效果持久，手感柔软。广泛应用于口罩、空气过滤、水过滤、真空吸尘器粉尘过滤等。

Our company can produce meltblown fabric with different protection efficiency and low resistance. The imported melt nozzle and advanced production technology make our meltblown fabric meet the domestic and oversea requirements of different standards for protective, medical and sanitary masks. Our meltblown fabric has stable performance, lasting electrified effect and soft handle. It's widely used for filtration of face mask, air filter, water filter, vacuum cleaner.



### 熔喷布等级

Grade of Meltblown fabric

BFE95, BFE99  
PFE90, PFE95, PFE99

我司不同等级的无纺布适用于以下标准：  
Our meltblown fabric with different grades are applicable to the following standards:

YY/T 0969  
YY 0469  
GB2626 KN90, KN95, KF94  
GB/T 32610  
EN14683 Type I, Type II, Type IIR

# 手术衣系列 SURGICAL GOWNS

## 一次性手术衣

Disposable Surgical Gown

### 产品信息

产品名称：一次性手术衣  
材料：SMS/SMMS/SSMMS  
克重：30-50克  
颜色：卫生蓝，卫生绿  
规格尺寸：S-3XL  
袖口：橡筋或针织袖  
款式：领口系带或魔术贴  
腰部2根或4根带子  
工艺：超声波车缝  
消毒方式：EO灭菌

### Product Information

Name: Surgical gown, operating gown  
Material: SMS/SMMS/SSMMS  
Weight: 30gsm-50gsm  
Color: Medical Blue, Green, or customized  
Size: S-M-L-XL-XXL-XXXL  
Sleeve: Elastic or Knitted cuff  
Style: Collar-Tie/Velcro,  
waist 2 or 4 ties,  
Sewing: Ultrasonic Welding  
Sterilization: Non-sterile or EO Sterile is applicable

Product Name	Specification	Packaging
Disposable Surgical Gown	35gsm, Non-Sterile	10pcs/Bag, 100pcs/Case
Disposable Surgical Gown	40gsm, Non-Sterile	10pcs/Bag, 100pcs/Case
Disposable Surgical Gown	45gsm, Sterile	1pcs/Pouch Bag, 50pcs/Case
Disposable Surgical Gown	50gsm, Sterile	1pcs/Pouch Bag, 50pcs/Case



## 加强型手术衣

Reinforced Surgical Gown


### 产品信息

产品名称：加强型手术衣  
材料：SMS/SMMS/SSMMS, PE复合材料  
克重：35-50克  
颜色：卫生蓝，卫生绿  
规格尺寸：S-3XL  
袖口：针织袖  
款式：领口系带或魔术贴  
腰部2根或4根带子  
工艺：超声波车缝  
消毒方式：EO灭菌

### Product Information

Name: Surgical gown, operating gown  
Material: SMS/SMMS/SSMMS, PE laminated fabric.  
Weight: 35gsm-50gsm  
Color: Medical Blue, Green, or customized  
Size: S-M-L-XL-XXL-XXXL  
Sleeve: Knitted cuff  
Style: Collar-Tie/Velcro,  
waist 2 or 4 ties,  
Sewing: Ultrasonic Welding  
Sterilization: Non-sterile or EO Sterile is applicable

Product Name	Specification	Packaging
Poly-Reinf.Surgical Gown	35gsm, Sterile	1pcs/Pouch Bag, 50pcs/Case
Poly-Reinf.Surgical Gown	40gsm, Sterile	1pcs/Pouch Bag, 50pcs/Case
Poly-Reinf.Surgical Gown	45gsm, Sterile	1pcs/Pouch Bag, 50pcs/Case
Poly-Reinf.Surgical Gown	35gsm, Sterile	1pcs/Pouch Bag, 50pcs/Case





## 一次性隔离衣

**Disposable Isolation Gown**  
产品备案号：豫信械备20210004号

### 产品信息

产品名称：一次性隔离衣  
材料：SPP, SMS  
克重：16-50克  
颜色：蓝色，黄色，绿色，白色  
规格尺寸：S-3XL  
袖口：橡筋或针织袖  
款式：领口系带或魔术贴  
腰部2根带子  
工艺：车缝或超声波

### Product Information

Name: Disposable Isolation Gown  
Material: SPP, SMS  
Weight: 16-50gsm  
Color: Blue, Yellow, Green, White or Customized  
Size: S-3XL  
Sleeve: Elastic or Knitted cuff  
Style: Collar-Tie/Velcro,  
waist 2 ties,  
Sewing: Sewing or Ultrasonic Welding

Product Name	Specification	Packaging
Disposable Isolation Gown	SPP 18gsm	10pcs/Bag, 100pcs/Case
Disposable Isolation Gown	SPP 25gsm	10pcs/Bag, 100pcs/Case
Disposable Isolation Gown	SPP 30gsm	10pcs/Bag, 100pcs/Case
Disposable Isolation Gown	SPP 35gsm	10pcs/Bag, 100pcs/Case

## 一次性医用防护服

**Disposable Medical Coverall**  
产品注册证号：豫械注准20202141053

### 产品信息

产品名称：一次性医用防护服  
材料：透气膜无纺布  
克重：50-65克  
颜色：白色  
规格尺寸：S-3XL  
袖口：橡筋袖  
款式：拉链双盖扁贴条、  
带帽、无鞋套  
消毒方式：非灭菌或EO灭菌

### Product Information

Name: Disposable Medical Coverall  
Material: Microporous Fabric  
Weight: 50-65gsm  
Color: White  
Size: S-M-L-XL-XXL-XXXL  
Sleeve: Elastic cuff  
Style: elastic cuff/ankle, zipper, seam sealing tape  
With hood, without boot  
Sterilization: Non-sterile or EO Sterile is applicable

Product Name	Specification	Packaging
Disposable Coverall	Microporous 50gsm	1pcs/Bag, 50pcs/Case
Disposable Coverall	Microporous 60gsm	1pcs/Bag, 50pcs/Case
Disposable Coverall	Microporous 65gsm	1pcs/Bag, 50pcs/Case



## 一次性覆膜隔离衣

**Impervious Isolation gown**  
产品备案号：豫信械备20210004号

### 产品信息

产品名称：一次性覆膜隔离衣  
材料：淋膜无纺布  
克重：25-50克  
颜色：蓝色，黄色，白色  
规格尺寸：S-3XL  
袖口：橡筋或针织袖  
款式：领口系带或魔术贴  
腰部2根带子  
工艺：车缝或超声波

### Product Information

Name: Impervious Isolation gown  
Material: PP+PE Lamination Fabric  
Weight: 25-50gsm  
Color: Blue, Yellow, White  
Size: S-3XL  
Sleeve: Elastic or Knitted cuff  
Style: Collar-Tie/Velcro,  
waist 2 ties,  
Sewing: Sewing or Ultrasonic Welding

Product Name	Specification	Packaging
Impervious Isolation Gown	PP+PE 25gsm	10pcs/Bag, 100pcs/Case
Impervious Isolation Gown	PP+PE 35gsm	10pcs/Bag, 100pcs/Case
Impervious Isolation Gown	PP+PE 40gsm	10pcs/Bag, 100pcs/Case

## 一次性普通防护服

**Disposable Coverall**

### 产品信息

产品名称：一次性普通防护服  
材料：PP/SMS  
克重：25-50克  
颜色：白色，蓝色  
规格尺寸：S-3XL  
袖口：橡筋袖  
款式：  
A连体带帽，橡筋袖，不带脚套  
B连体带帽和脚套，橡筋袖  
C翻领连体，橡筋袖，不带帽和脚套  
D翻领连体，不带帽和脚套

### Product Information

Name: Disposable Coverall  
Material: PP/SMS  
Weight: 25-50gsm  
Color: White, Blue  
Size: S-M-L-XL-XXL-XXXL  
Sleeve: Elastic cuff  
Style A: with hood, elastic cuff/ankle, zipper, without shoe cover.  
Style B: with hood and shoe cover, elastic cuff, zipper  
Style C: with collar, zipper, elastic cuff/ankle,  
without hood and shoe cover.  
Style D: with collar, zipper, opening cuff/ankle,  
without hood and shoe cover.,

Product Name	Specification	Packaging
Disposable Coverall	SPP 30gsm	1pcs/Bag, 50pcs/Case
Disposable Coverall	SPP 35gsm	1pcs/Bag, 50pcs/Case
Disposable Coverall	SPP/SMS 45gsm	1pcs/Bag, 50pcs/Case





## 一次性实验服

### Disposable Lab Coat

#### 产品信息

产品名称：一次性实验服  
 材料：PP/SMS  
 克重：25-50克  
 颜色：白色，蓝色，深蓝，紫色，粉色  
 规格尺寸：S-3XL  
 款式：单领或双领，翻领或针织领  
 橡筋袖或针织袖  
 扣子或魔术贴  
 口袋可有或无

#### Product Information

Name: Disposable Lab Coat  
 Material: PP/SMS  
 Weight: 25-50gsm  
 Color: White, Blue, Navy Blue, Purple, Pink.  
 Size: S-3XL  
 Style: Collar: Single collar/Double collar  
 Knitted collar/V collar  
 Cuff: Elastic/Knitted Wrists / Opening  
 With Button or Zipper or Velcro  
 Pocket: 1, 2, 3 or no pockets

Product Name	Specification	Packaging
Disposable Lab Coat	SPP 25gsm	1pcs/Bag, 50pcs/Case
Disposable Lab Coat	SPP 35gsm	1pcs/Bag, 50pcs/Case
Disposable Lab Coat	SPP 40gsm	1pcs/Bag, 50pcs/Case

## 一次性挂耳口罩

### Disposable Face Mask

产品注册证号：豫械注准20202141052

#### 产品信息

产品名称：一次性挂耳口罩  
 材料：PP无纺布+熔喷布+PP无纺布  
 规格尺寸：17.5X9.5CM  
 款式：弹力挂耳  
 颜色：蓝色，绿色，白色，粉色，黑色  
 过滤效率：BFE≥98%/≥95%  
 标准：YY/T0969-2013,  
 EN14683:2019 Type I/II/IIR

#### Product Information

Name: Disposable Face Mask  
 Material: PP+M+PP.  
 Size: 17.5\*9cm, 14.5x9cm for children  
 Type: Earloop  
 Color: Blue, Green, White, Pink, Black.  
 BFE: BFE≥98%/≥95%  
 Standard: YY/T0969-2013  
 EN14683: 2019 Type I/II/IIR

Product Name	Specification	Packaging
Disposable Face Mask	3ply, 17.5X9.5cm	50pcs/box, 2000pcs/case
Disposable Face Mask	3ply, 14.5X9.5cm	50pcs/box, 2000pcs/case



## 一次性盥洗服

### Disposable Scrub Suit

#### 产品信息

产品名称：一次性盥洗服（上衣+裤子）  
 材料：PP/SMS  
 克重：35-50克  
 颜色：蓝色，深蓝，绿色，紫色  
 规格尺寸：S-3XL  
 款式：V领或圆领，  
 腰部抽绳或橡筋  
 口袋可有或无

#### Product Information

Name: Disposable Scrub Suit  
 Material: PP/SMS  
 Weight: 35-50gsm  
 Color: Blue, Navy Blue, Green, Purple.  
 Size: S-3XL  
 Style: V shape /round shape collar  
 With short Sleeve and Pockets (1/2/3/4)  
 Drawstring waist/ elastic waist  
 Pants with 1/2 rear pocket

Product Name	Specification	Packaging
Disposable Scrub Suit	SPP/SMS 35gsm	1set/Bag, 50pcs/Case
Disposable Scrub Suit	SPP/SMS 40gsm	1set/Bag, 50pcs/Case
Disposable Scrub Suit	SPP/SMS 45gsm	1set/Bag, 50pcs/Case

## 一次性系带口罩

### Disposable Tie-on Face Mask

#### 产品信息

产品名称：一次性系带口罩  
 材料：PP无纺布+熔喷布+PP无纺布  
 规格尺寸：17.5X9.5CM  
 款式：系带  
 颜色：蓝色，绿色，白色  
 过滤效率：BFE≥98%/≥95%  
 标准：YY/T0969-2013,  
 EN14683:2019 Type I/II/IIR

#### Product Information

Name: Disposable Tie-on Face Mask  
 Material: PP+M+PP.  
 Size: 17.5\*9cm, 14.5x9cm for children  
 Type: Tie-on  
 Color: Blue, Green, White, Pink, Black.  
 BFE: BFE≥98%/≥95%  
 Standard: YY/T0969-2013  
 EN14683: 2019 Type I/II/IIR

Product Name	Specification	Packaging
Disposable Face Mask	3ply, 17.5X9.5cm	50pcs/box, 2000pcs/case
Disposable Face Mask	3ply, 14.5X9.5cm	50pcs/box, 2000pcs/case



## 一次性帽子

### Disposable Cap

#### 产品信息

款式：条帽、圆帽、医生帽  
 材料：PP/SMS/水刺无纺布  
 克重：9-30克  
 颜色：白色，蓝色，绿色，粉色  
 尺寸：18'' ,21'' 24''  
 工艺：机制、手工

#### Product Information

Name: Disposable Cap  
 Style: Clip Cap,Bouffant Cap,surgeons' cap  
 Material: PP/SMS/ spunlace fabric.  
 Weight: 9-30gsm  
 Color: White,blue, green,pink  
 Size: 18'',21''24''  
 Elastic: Single or Double elastic.  
 Feature: soft, breathable, fluid-resistant



Product Name	Specification	Packaging
Clip/Bouffant cap	18'', Blue/Green/White/Pink	100pcs/ bag, 2000pcs/ case
Clip/Bouffant cap	21'', Blue/Green/White/Pink	100pcs/ bag, 2000pcs/ case
Clip/Bouffant cap	24'', Blue/Green/White/Pink	100pcs/ bag, 2000pcs/ case
Surgeons' cap	SPP or spunlace	100pcs/ bag, 2000pcs/ case

## 一次性床罩

### Disposable table cover

#### 产品信息

名称：一次性床罩  
 材料：PP/SMS  
 克重：20-45克  
 颜色：白色，蓝色，粉色  
 款式：全筋或半筋  
 尺寸：:80\*180cm,80\*200m,120\*220cm

#### Product Information

Name: Disposable table cover  
 Material: PP/SMS  
 Weight: 20-45gsm  
 Color: white,blue,pink  
 Style: elastic for ending.  
 Size:80\*180cm,80\*200m,120\*220cm  
 Feature: Comfortable, breathable,  
 dust-proof,fluid-resistant.



Product Name	Specification	Packaging
Table cover	SPP20-45gsm	10pcs/ bag, 100pcs/ case
Table cover	SMS20-45gsm	10pcs/ bag, 100pcs/ case

## 一次性鞋套

### Disposable Shoe Cover

#### 产品信息

名称：一次性鞋套  
 款式：普通、防滑  
 材料：PP/SMS/PE/CPE  
 克重：25-40克  
 颜色：白色，蓝色，绿色，粉色  
 尺寸：40\*15cm, 41\*16cm

#### Product Information

Name: Disposable Shoe cover  
 Style:Normal or Anti-Skid  
 Material: PP/SMS /PE/CPE plastic material  
 Weight: 20-40gsm  
 Color: white, blue, green,pink  
 Size:40\*15cm, 41\*16cm  
 Feature: Comfortable, Non-slip, dust-proof,fluid-resistant.



Product Name	Specification	Packaging
Standard shoe cover	SPP 28-40gsm	100pcs/ bag, 2000pcs/ case
Non-skid shoe cover	SPP 35-40gsm	100pcs/ bag, 2000pcs/ case
CPE shoe cover	2.3g CPE plastic material	100pcs/ bag, 5000pcs/ case

## 一次性床单

### Disposable bed sheet

#### 产品信息

名称：一次性床单  
 材料：PP/SMS  
 克重：20-45克  
 颜色：白色，蓝色，粉色  
 包装方式：片装或卷装  
 尺寸：:80\*180cm,80\*200m,120\*220cm

#### Product Information

Name: bed sheet  
 Material: SPP ,SMS,PP+PE Lamination fabric.  
 Weight: 20-45gsm  
 Color: white,blue,pink  
 Packing Style: individual package or roll fabric package.  
 Size:80\*180cm,80\*200m,120\*220cm  
 Feature: Comfortable, breathable,  
 dust-proof,fluid-resistant.



Product Name	Specification	Packaging
Bed sheet	SPP 20-45gsm	10pcs/ bag, 100pcs/ case
Bed sheet	SMS 20-45gsm	10pcs/ bag, 100pcs/ case
Bed sheet	PP+PE 30gsm	10pcs/ bag, 100pcs/ case