

Test Report

No.: TJHL1912007233SD

Date: JAN.03, 2020

Page: 1 of 2

ZHANGJIAGANG XIEHE APPARATUS & INSTRUMENTS CO., LTD
FUQIAN ECONOMIC DEVELOPMENT ZONE, YANGSHE TOWN, ZHANGJIAGANG CITY, JIANGSU
PROVINCE, CHINA, ZIP:215600

The following sample(s) was/were submitted and identified by the client as:

Sample Description : SCOOP STRETCHER
Manufacture : ZHANGJIAGANG XIEHE APPARATUS & INSTRUMENTS CO., LTD
Country of Origin : CHINA
Sample Receiving Date : DEC.27, 2019
Testing Period : DEC.27, 2019 TO JAN.03, 2020
Test Performed : SELECTED TEST(S) AS REQUESTED BY APPLICANT
Test Requested : PLEASE REFER TO THE RESULT SUMMARY.
Test Result(s) : FOR FURTHER DETAILS, PLEASE REFER TO THE
FOLLOWING PAGE(S)
Test Summary :

Test Requested	Conclusion
EN 1865-1:2010+A1:2015 CLASUE 4.6.8, CLASUE 4.6.11 & CLASUE 4.6.12	Pass

Signed for and on behalf of
SGS-CSTC Standards Technical Services (Tianjin) Co., Ltd.

Michael



Michael Wen
Authorized Signatory



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

Building C-9, No.80, Haiyun Avenue, TEDA, Tianjin, China 300457 t:(86-22) 59995061
中国·天津·经济技术开发区海云街80号C-9厂房 邮编: 300457 t:(86-22) 59995061

www.sgs.com.cn
sgs.china@sgs.com

Test Report

No.: TJHL1912007233SD

Date: JAN.03, 2020

Page: 2 of 2

Test Conducted:

EN 1865-1:2010+A1:2015 CLASUE 4.6.8, CLASUE 4.6.11 & CLASUE 4.6.12

Nos. of Specimen : 1 pc (s) .

Test Result : Pass. Details shown as following table

Clause	Test Requirement	Result	Note
4.6.8 Deformation of the frame	The frame of the pick up stretcher shall not bend or break when tested in accordance with 5.5.1. 250kg	Pass	--
4.6.11 Deformation of the lying area	There shall be no remaining deformation of the lying area when tested in accordance with 5.1.3. 150kg	Pass	--
4.6.12 Resistance to torsion	There shall be no remaining deformation when tested in accordance with 5.5.2.	Pass	--

Sample Photo:

Tested sample view 1



Tested sample view 2



SGS authenticate the photo on original report only

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

Building C-9, No.80, Haiyun Avenue, TEDA, Tianjin, China 300457 t:(86-22) 59995061
中国·天津·经济技术开发区海云街80号C-9厂房 邮编: 300457 t:(86-22) 59995061

www.sgsgroup.com.cn
sgs.china@sgs.com