


























Price list of stretcher

Item	Model	Photo for reference	Specification	Unit Price US\$	MOQ
Aluminum Alloy folding stretcher (2 folder)	YXZ-D-B1		Unfolded size: 185*50*20cm Folded size:92*50*9cm Packing size(2pcs/carton):96*54*19cm N.W: 6kg G.W:13kg Load bearing: 159kg	65	20
Aluminum Alloy folding stretcher (2 folder)	YXZ-D-B2		Unfolded size: 185*50*18cm Folded size:92*50*9cm Packing size(2pcs/carton):96*54*19cm N.W: 6kg G.W:13kg Load bearing: 159kg	65	20
Aluminum Alloy folding stretcher (2 folder)	YXZ-D-A1		Unfolded size: 185*50*3cm Folded size:92*50*7cm Packing size(2pcs/carton):96*54*16cm N.W: 4kg G.W:9kg Load bearing: 159kg	40	20
Powder coated steel folding stretcher (2 folder)	YXZ-D-A1		Unfolded size: 185*50*3cm Folded size:92*50*10cm Packing size(2pcs/carton):96*54*16cm N.W: 6.5kg G.W:14kg Load bearing: 159kg	30	20
Aluminum Alloy folding stretcher (2 folder)	YXZ-D-B3		Unfolded size: 185*50*24cm Folded size:92*50*11cm Packing size(1pcs/carton):96*54*12cm N.W: 6kg G.W: 9kg Load bearing: 159kg	70	20
Aluminum Alloy camp folding stretcher	YXZ-D-B4		Unfolded size: 194*66*41cm Folded size:90*12*10cm Packing size(1pcs/carton):92*26*15cm N.W: 8.5kg G.W: 9.5kg Load bearing: 159kg	85	20
Aluminum Alloy Military folding stretcher (4 folder)	YXZ-D-B5		Unfolded size: 185*48*21cm Folded size:46*48*14cm Packing size(1pcs/carton):50*48*15cm N.W: 7kg G.W: 8.2kg Load bearing: 159kg	65	20
Aluminum Alloy Military folding stretcher (4 folder)	YXZ-D-B6		Unfolded size: 185*48*21cm Folded size:46*48*14cm Packing size(1pcs/carton):50*48*15cm N.W: 4.5kg G.W: 6.2kg Load bearing: 159kg	50	20
Aluminum Alloy folding stretcher (2 folder)	YXZ-D-A5 YXZ-D-A6		Unfolded size: 185*50*3cm Folded size:92.5*50*7cm Packing size(2pcs/carton):96*54*16cm N.W: 6.5kg G.W: 14kg Load bearing: 159kg	68	20
Aluminum Alloy Military folding stretcher (2 folder)	YXZ-D-C1		Unfolded size: 220*55*15cm Folded size:180*12*15cm Packing size(1pcs/carton):183*14*16cm N.W: 7.5kg G.W: 9.5kg Load bearing: 159kg	68	20

<p>Aluminum Alloy Military folding stretcher (2 folder)</p>	YXZ-D-C2		<p>Unfolded size: 220*55*15cm Folded size:95*20*9cm Packing size(1pcs/carton):98*22*11cm N.W: 7.5kg G.W: 8.2kg Load bearing: 159kg</p>	60	20
<p>Aluminum Alloy Military folding stretcher (4 folder)</p>	YXZ-D-C3		<p>Unfolded size: 220*55*15cm Folded size:50*26*18cm Packing size(1pcs/carton):56*28*20cm N.W: 9kg G.W: 9.8kg Load bearing: 159kg</p>	130	20
<p>Aluminum Alloy Military folding stretcher (4 folder)</p>	YXZ-D-D1		<p>Unfolded size: 214*54.5*12cm Folded size:58*21*14cm Packing size(1pcs/carton):67*30*17cm N.W: 5kg G.W: 7kg Load bearing: 159kg</p>	65	20
<p>Aluminum Alloy Military folding stretcher</p>	YXZ-D-C4		<p>Unfolded size: 215*54.5*12cm Folded size:182*9*17cm Packing size(1pcs/carton):185*11*19cm N.W: 6.5kg G.W: 8kg Load bearing: 159kg</p>	45	20
<p>Aluminum Alloy folding stretcher (2 folder)</p>	YXZ-C-C4		<p>Unfolded size: 208*54.5*13cm Folded size:104*17*9cm Packing size(2pcs/carton):108*18.5*18cm N.W: 4.7kg G.W: 11kg Load bearing: 159kg</p>	45	20
<p>Aluminum Alloy folding stretcher (2 folder)</p>	YXZ-D-C5		<p>Unfolded size: 208*54.5*15cm Folded size:104*17*9cm Packing size(2pcs/carton):108*18.5*18cm N.W: 5kg G.W: 12kg Load bearing: 159kg</p>	58	20
<p>Aluminum Alloy folding stretcher (4 folder)</p>	YXZ-D-C6		<p>Unfolded size: 208*54.5*15cm Folded size:104*17*9cm Packing size(1pc/carton):108*18.5*18cm N.W: 5kg G.W: 12kg Load bearing: 159kg</p>	100	20
<p>Aluminum Alloy folding stretcher (2 folder)</p>	YXZ-D-C8		<p>Unfolded size: 210*50*12cm Folded into wheel chair:137cm Packing size(1pc/carton):100*52*14cm N.W: 8kg G.W: 9.5kg Load bearing: 159kg</p>	98	20
<p>Aluminum Alloy stretcher</p>	YXZ-D-C7		<p>Unfolded size:187*47*18cm Packing size(1pc/carton):190*50*5cm N.W: 5kg G.W: 6kg Load bearing: 159kg</p>	65	20

Aluminum Alloy stair stretcher	YXZ-D-C8 YXZ-D-C9		Unfolded size:90*51*91cm Folded size:90*17*59cm Packing size(1pc/carton):96*17*64cm N.W: 8kg G.W: 9.2kg Load bearing: 159kg	C8:100 C9:110	20
Aluminum Alloy stair stretcher	YXZ-D-C10		Unfolded size:90*51*91cm Folded size:90*17*59cm Packing size(1pc/carton):96*17*64cm N.W: 10.5kg G.W: 11.7kg Load bearing: 159kg	125	20
Aluminum Alloy scoop stretcher	YXZ-D-E1 YXZ-D-E2 YXZ-D-E4		Unfolded size:190*44*6cm Folded size:120*44*9cm Packing size(1pc/carton):170*47*9cm N.W: 8.5kg G.W: 10.2kg Load bearing: 159kg	E1:105 E2:110 E4:200	20
Aluminum Alloy scoop stretcher	YXZ-D-E3		Unfolded size:214*42*7cm Folded size:167*42*7cm Packing size(1pc/carton):174*45*10cm N.W: 8kg G.W: 11kg Load bearing: 159kg	255	20
Basket stretcher	YXZ-D-5A		Unfolded size:216*61*19cm Packing size(1pc/carton):223*66*24cm N.W: 18kg G.W: 20kg Load bearing: 159kg	260	20
Basket stretcher (Detachable)	YXZ-D-5B		Unfolded size:216*61*19cm Folded size:130*61*26cm Packing size(1pc/carton):135*64*28cm N.W: 20.5kg G.W:22.5kg Load bearing: 159kg	275	20
Spine board	YXZ-D-A1		Unfolded size:184*46*5cm Packing size(1pc/carton):186*46*7cm N.W: 7kg G.W:8.5kg Load bearing: 159kg	50	20
Spine board	YXZ-D-A2		Unfolded size:186*45*5.5cm Packing size(1pc/carton):186*47*8cm N.W: 7.5kg G.W:9kg Load bearing: 159kg	45	20
Spine board	YXZ-D-A3		Unfolded size:184*45*5cm Packing size(1pc/carton):186*47*7cm N.W: 6.5kg G.W:8kg Load bearing: 159kg	45	20
Spine board	YXZ-D-A4		Unfolded size:183*45*5cm Packing size(1pc/carton):185*47*7cm N.W: 6.5kg G.W:7.5kg Load bearing: 159kg	45	20

Head immobilizer	YXZ-D-W		<ol style="list-style-type: none"> Made of high-density plastic materials, large ear holes for monitoring the patient's ear canal X-Ray, CT scanning, MRI noninterference It is waterproof, easy to clean and anti-bacterial Size: 40*25*16cm Packing: 1pc/carton 43.5*27*14cm N.W/G.W: 0.9kg/1.1kg 	40	20
soft foam neck collar	YXZ-D-01		<ol style="list-style-type: none"> with lock system to fix position Designed with hole on throat position, it is convenient for the treatment of throat and trachea The height can be adjustable in 4 different levels shape fits for human neck with special material, allow X-ray, CAT and MRI scan size adjustable to fit all kinds of patients keep clean for reuse, prevent bacterial growth on it Unfolding size:42*26*18cm Packing(40pcs/big carton): 60*50*48cm/ carton N.W.:8kg G.W.:13kg 	7	80
Soft PVC stretcher	YXZ-D-D2		<p>Product size:178*54cm N.W:1kg Load bearing:159kg Packing(10pcs/carton):89*25*40 G.W: 11kg</p>	30	20
Soft PVC stretcher	YXZ-D-D3		<p>Product size:212*75cm N.W:1kg Load bearing:159kg Packing(10pcs/carton):89*25*40 G.W: 11kg</p>	38	20
Aluminum alloy ambulance stretcher	YXZ-D-G1		<p>Unfolded size:190*55*50cm Adjustable back angle:70 degree Packing size(1pc/carton):195*62*34cm N.W: 25.5kg G.W: 31.5kg Load bearing: 159kg</p>	218	10
Aluminum alloy ambulance stretcher (including the scoop stretcher)	YXZ-D-G2		<p>Unfolded size:190*55*50cm Adjustable back angle:70 degree Packing size(1pc/carton):195*62*34cm N.W: 29kg G.W: 35kg Load bearing: 159kg</p>	335	10
Aluminum alloy ambulance stretcher	YXZ-D-H1		<p>Unfolded size:190*55*26~92cm Adjustable back angle:70 degree Packing size(1pc/carton):195*62*28cm N.W: 36.8kg G.W: 43kg Load bearing: 159kg</p>	380	10

Aluminum alloy ambulance stretcher	YXZ-D-H2		Unfolded size:190*55*26~92cm Adjustable back angle:70 degree Packing size(1pc/carton):195*62*28cm N.W: 43kg G.W: 49kg Load bearing: 159kg	410	10
Aluminum alloy ambulance stretcher	YXZ-D-H3		Unfolded size:190*55*25~92cm Adjustable back angle:65 degree Packing size(1pc/carton):195*62*25.5cm N.W: 40kg G.W: 46kg Load bearing: 159kg	335	10
Aluminum alloy ambulance stretcher	YXZ-D-H4		Unfolded size:190*55*25~92cm Adjustable back angle:65 degree Packing size(1pc/carton):195*62*28cm N.W: 39kg G.W: 42kg Load bearing: 159kg	350	10
Aluminum alloy ambulance stretcher	YXZ-D-H5		Unfolded size:190*55*36~98cm Adjustable back angle:75 degree Packing size(1pc/carton):198*62*38cm N.W: 39kg G.W: 45kg Load bearing: 159kg	435	10
Aluminum alloy ambulance stretcher	YXZ-D-3E		Unfolded size:196*55*25~86cm Can be changed into a wheelchair Packing size(1pc/carton):198*64*26cm N.W: 34kg G.W: 40kg Load bearing: 159kg	500	10
Aluminum alloy ambulance stretcher	YXZ-D-3G		Unfolded size:196*55*72~88cm Packing size(1pc/carton):195*62*38cm N.W: 34.5kg G.W: 40.5kg Load bearing: 159kg	350	10
Aluminum alloy ambulance stretcher	YXZ-D-I		Unfolded size:196*55*34~88cm Adjustable back angle:70 degree Packing size(1pc/carton):198*62*38cm N.W: 38.2kg G.W: 44kg Load bearing: 159kg	450	10
Aluminum alloy ambulance stretcher	YXZ-D-I2		Unfolded size:196*55*34~88cm Adjustable back angle:75 degree Packing size(1pc/carton):198*62*38cm N.W: 42kg G.W: 48kg Load bearing: 159kg	600	10

Aluminum alloy stretcher base	YXZ-D-03		Unfolded size:196*59*45cm Packing size(1pc/carton):232*62*15cm N.W: 35kg G.W: 41kg Load bearing: 159kg	300	10
Stainless steel stretcher base	YXZ-D-04		Unfolded size:203*67*40cm Packing size(1pc/carton):238*68*18cm N.W: 27g G.W: 43kg Load bearing: 159kg	485	10
Stainless steel stretcher trolley	YXZ-D-J1		Unfolded size:195*60*80cm Packing size(1pc/carton):198*73*19cm N.W: 34g G.W: 40kg Load bearing: 159kg	205	10
Aluminum alloy stretcher trolley (Height adjustable)	YXZ-D-J		Unfolded size:195*55*88cm Folded size:195*55*54cm Packing size(1pc/carton):195*62*34cm N.W: 34g G.W: 40kg Load bearing: 159kg	220	10
Aluminum alloy stretcher trolley (Height adjustable) (Back adjustable)	YXZ-D-K		Unfolded size:195*55*88cm Folded size:195*55*54cm Adjustable back angle:70 degree Packing size(1pc/carton):195*62*34cm N.W: 34g G.W: 41kg Load bearing: 159kg	240	10
Vacuum mattress stretcher	YXZ-D-05		Product size:212*75*7cm Folded size:76*60*37cm Packing size(1pc/carton):95*28*58cm Load bearing: 159kg N.W:5kg G.W:7kg	268	10
Stretcher base	YXZ-D-07			75	10