



HAVING HEALTH,
HAVING EVERYTHING !



PRODUCT CATALOG

Anhui Gloria Import and Export Co.,Ltd

2022



Anhui Gloria Import and Export Co.,Ltd

E-mail: admin@gloriacn.com

Tel/Fax No: 86-556-5202676

Online Website: www.gloriacn.com/www.otssuperfine.com

Address: No 2620 ,B2 Building, Wuyue Plaza, Yingjiang District, Anqing City, Anhui, China. 246000

1st Factory Address: Tongcheng Qingcao Industrial Park, Anhui Province, China, 231490

2nd Factory Address: D1, Jiabao Industrial Zone, Bridge Development Area, Anqing, Anhui, China. 246000

- + STERILE BARRIE SYSTEM
- + STERILE CONTROL SYSTEM
- + STERILE PACKING MATERIAL
- + STERILE DENTAL SYSTEM
- + MEDICAL FACE MASK
- + MEDICAL CONSUMABLE SUPPLIES



COMPANY PROFILE

COMPANY PROFILE

- STERILE BARRIE SYSTEM
- STERILE CONTROL SYSTEM
- STERILE PACKING MATERIAL
- STERILE DENTAL SYSTEM
- MEDICAL CONSUMABLE SUPPLIES
- MEDICAL FACE MASK

FACTORY ADVANTAGE



ABOUT US

Anhui Gloria Import and Export Co.,Ltd is one of largest manufacturer and distributor of sterilization pouches and chemical indicators supplies in China, has three production bases, eight production workshops, covering a total area of 20,000 square meters, including purification workshop of 10,000 square meters providing high quality and price competitive products that servemore than 50 countries world-wide.We now providing reliable, efficient and flexible solutions for worldwide customers with more than 5,000 products.

Our innovative products can be found in many hospitals, surgery centers, healthcare institutions, medical devicemanufacturers,healthcare dealers and retail outlets over the world.

We consistently provide outstanding value and industry-leading quality sterilization products,we help customers deliver quality health care by making the products they want, with the high quality and the value they demand. Also, we offer our customers with OEM services and custom solutions to help solve sterilization issues they need. Gloria offers its quality stands for real China quality & updated tech of Made-in-China.We are making great efforts to achieve this win-win situation and sincerely welcome you to join us.

Leading manufacturer Medical packaging

Industry and trade integration

CLEAN WORKSHOP

Advanced production and processing equipment.

SCIENTIFIC MANAGEMENT SYSTEM

Product development, design, planning production and sales in one



COMPLETE QUALITY CONTROL SYSTEM

ISO9001, ISO13485, ISO, FDA, CE and QS systems



ENTERPRISE CERTIFICATE

All products have passed the certification of ISO9001, ISO11607, CE, QS and FDA system.



▶ TABLE OF CONTENTS

01 Sterilization Packaging Pouches

- Self-sealing Sterilization Pouches
- Heat-sealing Sterilization Pouches
- Tyvek Pouch & Reel
- Heat-sealing Sterilization Reel Gusseted
- Sterilization Paper Bag Gusseted Paper

02 Medical Consumable Supplies

- Medical Crepe Paper
- Non-woven Fabric
- Disposable Coverall
- Surgical Gown + Isolation Gown
- Mob Cap
- Shoe Cover + Sleeve Cover + Bed Cover

04 Sterile Dental System

- Dental Bibs
- Saliva Ejector
- Barrier Film
- Micro Applicator
- Syringe Tips
- Interdental Brush

03 Medical Gauze Sponge

- Non-sterile Gauze Swab
- Sterile Gauze Swab
- Cotton Ball
- Dental Cotton Roll
- Surgical Tape

05 Sterile Indicator Control System

- Steam/EO/Plasma Indicator Tape
- Steam/EO/Plasma Indicator Strip
- EO Indicator Label
- Bowie Dick Test Pack

06 Medical Face Mask

- 3Ply Medical TYPE II/TYPE IIR Mask
- KF94/KF99 Fish Face Mask
- FFP2 Face Mask

07 Disposable Sterile Dressing Pack

- Disposable Sterile Dressing Pack
- Blood Taking Set
- Delivery Procedural Set

»» Sterilization Packaging Pouches



Self-sealing Sterilization Pouches

- ◆ Medical grade paper 60g and CPP/PET 52um film
- ◆ Steam and EO Sterilization
- ◆ Blue to Black under Steam Sterilization, Pink to Yellow under EO Sterilization

Item No	Size	Packing
GCSF-057130	57mm*130mm	200pcs/box,30box/ctn
GCSF-070260	70mm*260mm	200pcs/box,15box/ctn
GCSF-090165	90mm*165mm	200pcs/box,20box/ctn
GCSF-090260	90mm*260mm	200pcs/box,10box/ctn
GCSF-135280	135mm*280mm	200pcs/box,10box/ctn
GCSF-190360	190mm*360mm	200pcs/box,5box/ctn
GCSF-200320	200mm*320mm	200pcs/box,5box/ctn
GCSF-200360	200mm*360mm	200pcs/box,5box/ctn
GCSF-230395	230mm*395mm	200pcs/box,5box/ctn
GCSF-305430	305mm*430mm	200pcs/box,5box/ctn
GCSF-400535	400mm*535mm	200pcs/box,5box/ctn
GCSF-500630	500mm*630mm	200pcs/box,5box/ctn



Heat-sealing Sterilization Pouches

- ◆ Medical grade paper 60g and CPP/PET 52um film
- ◆ Steam and EO Sterilization
- ◆ Blue to Black under Steam Sterilization, Pink to Yellow under EO Sterilization

Item No	Size	Packing
GCHR-050200	50mm*200m	12roll/ctn
GCHR-075200	75mm*200m	8roll/ctn
GCHR-100200	100mm*200m	6roll/ctn
GCHR-120200	120mm*200m	4roll/ctn
GCHR-125200	125mm*200m	4roll/ctn
GCHR-150200	150mm*200m	4roll/ctn
GCHR-200200	200mm*200m	2roll/ctn
GCHR-250200	250mm*200m	2roll/ctn
GCHR-300200	300mm*200m	2roll/ctn
GCHR-350200	350mm*200m	2roll/ctn
GCHR-400200	400mm*200m	2roll/ctn
GCHR-500200	500mm*200m	2roll/ctn



Sterilization >>> Packaging Pouches

Tyvek Pouch & Reel

- ◆ Dupont Tyvek paper 2FS/1059B/1073B + Laminated PE easy peel film
- ◆ Excellent microbial barrier properties
- ◆ Easy peel and aseptic presentation
- ◆ Superior tear strength and puncture resistance
- ◆ One side is transparent film, easy for identification on packed
- ◆ Clear and accurate color change from purple to red
- ◆ Suitable for plasma(H2O2) , ETO gas, Electron - beam sterilization
- ◆ Purple to red after Plasma sterilization



Item No.	Size	Packing
GCTR-07570	75mm*100m	8roll/ctn
GCTR-10070	100mm*70m	6roll/ctn
GCTR-15070	150mm*70m	4roll/ctn
GCTR-20070	200mm*70m	2roll/ctn
GCTR-25070	250mm*70m	2roll/ctn
GCTR-30070	300mm*70m	2roll/ctn
GCTR-35070	350mm*70m	2roll/ctn
GCTR-40070	400mm*70m	2roll/ctn
GCTR-075100	75mm*100m	8roll/ctn
GCTR-100100	100mm*100m	6roll/ctn
GCTR-150100	150mm*100m	4roll/ctn
GCTR-200100	200mm*100m	2roll/ctn
GCTR-250100	250mm*100m	2roll/ctn
GCTR-300100	300mm*100m	2roll/ctn
GCTR-350100	350mm*100m	2roll/ctn
GCTR-400100	400mm*100m	2roll/ctn



Heat-sealing Sterilization Reel Gusseted

- ◆ Medical grade paper 60g and CPP/PET 52um film
- ◆ Steam and EO Sterilization
- ◆ Blue to Black under Steam Sterilization, Pink to Yellow under EO Sterilization



Item No	Size	Packing
GCGR-075100	75mm*100m	8roll/ctn
GCGR-100100	100mm*100m	6roll/ctn
GCGR-150100	150mm*100m	4roll/ctn
GCGR-200100	200mm*100m	2roll/ctn
GCGR-250100	250mm*100m	2roll/ctn
GCGR-300100	300mm*100m	2roll/ctn
GCGR-350100	350mm*100m	2roll/ctn
GCGR-400100	400mm*100m	2roll/ctn

Sterilization Paper Bag Gusseted Paper

- ◆ Turn to Black from Blue under steam sterilization
- ◆ Sterilization with Steam, EO, FORM and Gamma
- ◆ Medical grade paper with technic 60gsm or 70gsm, an effective barrier against bacteria
- ◆ Three-dimensional can increase the volume of the pouch
- ◆ Clean and accurate indicator color change



Item No	Size	Pack(pcs/ctn)
GPB-A	90*50*150mm	1000
GPB-B	90*50*250mm	1000
GPB-C	110*30*190mm	1000
GPB-D	140*75*250mm	1000
GPB-F	140*50*330mm	1000
GPB-G	140*50*360mm	1000
GPB-H	180*95*380mm	500
GPB-J	100*50*150mm	1000
GPB-L	250*100*380mm	500
GPB-P	190*65*330mm	500
GPB-Q	110*45*190mm	1000

Medical >>> Consumable Supplies

Medical Crepe Paper

- ◆ 45g, 60g are available upon request
- ◆ Medical grade paper with no smell, no powder and non-toxic
- ◆ Excellent water resistance and bursting strength
- ◆ High bacterial filtration efficiency (BFE)
- ◆ Prevent microbial barrier



Item No	Size	Packing
GCCP-3030	30cm*30cm	2000pcs/ctn
GCCP-4040	40cm*40cm	1000pcs/ctn
GCCP-4545	45cm*45cm	1000pcs/ctn
GCCP-5050	50cm*50cm	1000pcs/ctn
GCCP-6060	60cm*60cm	1000pcs/ctn
GCCP-7575	75cm*75cm	500pcs/ctn
GCCP-9090	90cm*90cm	250pcs/ctn
GCCP-100100	100cm*100cm	250pcs/ctn
GCCP-120120	120cm*100cm	250pcs/ctn

Non-woven Fabric

- ◆ SMS (Spunbond+ meltblown+spunbond)
- ◆ SMMS (Spunbond+meltblown+meltblown+spunbond)
- ◆ SMMMS (Spunbond+meltblown+meltblown+meltblown+spunbond)
- ◆ Made from environmental friendly cellulose material, 100% Polypropylene with features of non toxic, no fiber exfoliation and high bacterial resistance.



Item No	Size	Packing
GCNF-5050	50cm*50cm	500pcs/ctn
GCNF-6060	60cm*60cm	500pcs/ctn
GCNF-7575	75cm*75cm	250pcs/ctn
GCNF-9090	90cm*90cm	250pcs/ctn
GCNF-100100	100cm*100cm	200pcs/ctn
GCNF-120120	120cm*120cm	125pcs/ctn
GCNF-120140	120cm*140cm	125pcs/ctn
GCNF-120160	120cm*160cm	125pcs/ctn
GCNF-120180	120cm*180cm	125pcs/ctn

Disposable Coverall >>>



GC-DM01

Package: 1pcs/ polybag,50pcs/ctn

Fabric: SPP/Microporous

Type: With Hood Without Boot;Elastic Hood, Cuffs,Ankles;

Half Elastic on Waist;Zipper on Front;

Color: White, Blue,Green, Yellow, Dark blue, other colors also available.

DISPOSABLE COVERALL



Fabric	<p>PP/SBPP anti dust, water resistance, breathable, different color available</p> <p>PP+PE anti dust, waterproof, different color available</p> <p>SMS good texture, fine tensile force, good water resistance, anti-bacterial</p> <p>SF/Microporous waterproof, breathable, soft touch, light weight, anti-dust, anti static</p>
Long coat size	S(Length 105x125cm), M(Length 110x130cm), L(Length 115x135cm) XL(Length 120x140cm), 2XL(Length 125x145cm) Can produce per customized size.
Short jacket size	S(Length 75x125cm), M(Length 80x130cm), L(Length 85x135cm) XL(Length 90x140cm), 2XL(Length 95x145cm) Can produce per customized size.
Color	White, blue, green, pink, black etc.
Packing	1pc/bag, 50pcs/carton; 5 or 10pcs/bag, 50 or 100pcs one carton
Application	Medical, Laboratory, dust free workshop etc. any place need protection
Style&Product code	<ul style="list-style-type: none"> with collar, without pocket, open wrist with collar, without pocket, elastic cuff knitted collar, with pockets, knitted cuff with collar, with pockets, knitted cuff with collar, with pockets, elastic cuff
We support customized specification	

Disposable Isolation Gown >>

Disposable Isolation Gown

GC-DM02

Packing: 10pcs/polybag, 100pcs/ctn

Material: PP non-woven/PP+PE






Size: 115*137cm, 120*140cm

Style: Elastic cuff, ties on back waist, back neck ties or hook&loop

Color: blue, white, yellow, other colors can customized



Disposable Gown

Fabric	<p>PP/SBPP anti dust, water resistance, breathable, different color available</p> <p>PP+PE anti dust, waterproof, different color available</p> <p>SMS good texture, fine tensile force, good water resistance, anti-bacterial</p> <p>SF/Microporous waterproof, breathable, soft touch, light weight, anti-dust, anti static</p>
Long coat size	<p>S(Length 105x125cm), M(Length 110x130cm), L(Length 115x137cm)</p> <p>XL(Length 120x140cm), 2XL(Length 125x145cm) Can produce per customized size.</p>
Packing	1pc/bag, 50pcs/carton; 5 or 10pcs/bag, 50 or 100pcs one carton
Application	Medical, Laboratory, dust free workshop etc. any place need protection
Style&Product code	<ul style="list-style-type: none">  with collar, without pocket, open wrist  with collar, without pocket, elastic cuff  knitted collar, with pockets, knitted cuff  with collar, with pockets, knitted cuff  with collar, with pockets, elastic cuff
We support customized specification	

Disposable Gown >>

Disposable Surgical Gown

GC-DM03

Package: 10pcs/poly bag, 100pcs/ctn

Material: SMS/ SMMS

Size: 115*137cm, 120*140cm

Style: Ties on Back Neck, Ties on Back Waist, Elastic or Knitted Cuff.

Cuff: Elastic-cuff or Knitted Cuff

Color: White, Blue, Green, Yellow, Dark blue, Dark green, other colors also available.



Disposable plastic gown/CPE apron

GC-DM04

Material: SPE, PE

Package: 1pcs/bag, 200pcs/carton

Color: White, Blue, Green, other colors also available.

Style: thumb hole on cuff, elastic cuff, no cuff



Disposable >>> Mob Cap

- ◆ Light weight, latex free, dust proof
- ◆ Automatic machine made, high capability
- ◆ Clip-shaped cap before unfolded
- ◆ Soft non-irritating elastic band
- ◆ Prevent hair falls



Disposable Mob Cap

GC-DM05

Package: 100pcs/poly bag, 1000pcs or 2000pcs/ctn

Fabric: PP non woven, SMS non woven, PE plastic

Style: Single elastic, Double elastic

Color: White, Blue, Green, Pink, Yellow, Black etc.



Disposable Bouffant Cap

GC-DM06

Package: 100pcs/poly bag, 1000pcs or 2000pcs/ctn

Fabric: PP, PP+PE, SMS, SF, PE Style: elastic

Color: White, Blue, Green, Pink, Yellow, Black etc.



Sleeve Cover >>>

Non-woven Sleeve Cover

GC-DM07

Style: Automatic machine made

Fabric: PP, SMS, PP+PE, Mircoporous

Color: White, Blue, Green, any color available

Package: 100pcs/bag, 10bags or 20bags/carton



PP/CPE Sleeve Cover

GC-DM08

Style: Machine made or hand made

Fabric: PP, CPE

Color: White, Blue, Green, any color available

Package: 100pcs/bag, 10bags or 20bags/carton.



Shoe Cover >>>

Disposable shoe cover

GC-DM09

Style: Without Non-Skid/Anti-Skid

Fabric: SPP/CPE/SPP+CPE

Color: Blue, Pink any color available

Package: 10pcs/roll, 10rolls/bag, 1000pcs or 2000pcs/ctn

Size: 16*41cm



Bed Cover

Bed Cover

GC-DM10

Fabric: Non-woven/SMS

Color: White, Green, Blue

Package: 10pcs/bag, 10bags/ctn

Size: 100*200cm, 110*210cm, 120*220cm

Weight: 30g/m², 35g/m², 40g/m², 50g/m², 60g/m²



Medical Gauze Sponge

Non-sterile Gauze Swab

GC-NGS



- ◆ Product Name: Non-sterile gauze swab (folded and unfolded, with and without X-ray)
- ◆ Material: 100% cotton
- ◆ Cotton Yarn: 40'S,32'S,21'S
- ◆ Different Gauze Mesh: 20*9,19*15,26*18,24*20,30*20
- ◆ Color: White, Green, Blue
- ◆ Packing: 100pcs/200pcs/Pack, 100packs/20Packs/CTN
- ◆ Gauze swab size: 2"x2", 3"x3", 4"x4", 4"x8"
- ◆ Number of Plies: 4ply, 8ply, 12ply, 16ply, 24ply

Sterile Gauze Swab

GC-SGS



- ◆ Product Name: Sterile gauze swab (folded and unfolded, with and without X-ray)
- ◆ Material: 100% cotton
- ◆ Cotton Yarn: 40'S,32'S,21'S
- ◆ Different Gauze Mesh: 20*9,19*15,26*18,24*20,30*20
- ◆ Color: White, Green, Blue
- ◆ Packing: 1pc/2pcs/pouch, 100pouches/50pouches/Box
- ◆ Blister Packing : 5pcs/10pcs/bag, 20bags/10bags/Box
- ◆ Gauze swab size: 2"x2", 3"x3", 4"x4", 4"x8"
- ◆ Sterile Method: By EO and Gamma
- ◆ Number of Plies: 4ply, 8ply, 12ply, 16ply, 24ply

Cotton Ball

GC-CB

- ◆ Material: 100% cotton
- ◆ Color: White, pink, blue, yellow, etc
- ◆ Size : Grams around 0.1g-1.2g/pc or customerized
- ◆ Type: Steril or Non-sterile
- ◆ Sterile: Blister pack
- ◆ Non-sterile: Sealed poly bag
- ◆ Packing: 0.3g/pc 500g/bag 20bags/ctn 0.5g/pc 500g/bag 20bags/ctn
1.0g/pc 500g/bag 20bags/ctn 1.0g/pc 100g/bag 10bags/ctn



Dental Cotton Roll

GC-CR

- ◆ Material: 100% cotton
- ◆ Diameter: 8x38mm; 10x38mm; 12x38mm; is available
- ◆ Feature: Surface flatness; Lint Free 150bags/ctn
- ◆ PACKING: 8x38mm 50pcs/bundle, 20bundle/bag, 20bags/ctn
10x38mm 50pcs/bundle, 20bundle/bag, 20bags/ctn
12x38mm 50pcs/bundle, 20bundle/bag, 20bags/ctn



Surgical Tape

GC-CR

- ◆ Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
- ◆ Length: 3m, 4m, 4.5m, 5m, 9m, 10m
- ◆ Glue: hot melt or acrylic acid glue



Specification	Packing Qty/Ctn	Ctn Size
1.25cm*5yds	24rolls/box, 30boxes/ctn	24*23.5*28.5
2.5cm*5yds	12rolls/box, 30boxes/ctn	24*23.5*28.5
5cm*5yds	6rolls/box, 30boxes/ctn	24*23.5*28.5
7.5cm*5yds	6rolls/box, 30boxes/ctn	24*23.5*41
10cm*5yds	6rolls/box, 30boxes/ctn	38.5*23.5*33.5

Sterile >> Dental System

Barrier Film

GC-BF

- ◆ Material: PE
- ◆ Perforation Line
- ◆ 4"*6"-1200pcs, length 180m,thickness: 0.05mm
- ◆ With easy peel finger lift
- ◆ Available in 5 colors: Blue/Clear/Pink/ Yellow/Green
- ◆ Suitable for barrier film shrink wrap package



Saliva Ejector

GC-SE

- ◆ Material: transparent PVC tube, copper interior wiring with chrome plated
- ◆ Sterilization: EO Gas
- ◆ Perforated wall to prevent tongue and mucous membranes of the patient to be sucked into the opening
- ◆ The curly saliva ejector has a further portion held in place by a clamping action between chain and front tooth to prevent pressure against the mouth bottom
- ◆ Comfortable ,Soft & Flexible



Dental Bibs

GC-DB

- ◆ Material: 100% Wooden Pulp and PE film which can prevent leakage of water stains
- ◆ Size: 33cm*45cm
- ◆ Color: pink, blue, green, purple, yellow, white, orange
- ◆ Features: water-proof, high water absorbent and anti-leakage protects patient from fluids
- ◆ packing: 125pcs/bag,500pcs/ctn



Micro Applicator

GC-MA

- ◆ Product Name:Micro Fiber Applicator Brush
- ◆ Material:100% Mao Bamboo
- ◆ Size:100*3mm
- ◆ Color:White, Yellow,Pink,Blue,Orange,Purple etc.
- ◆ Packing:100pcs/box
- ◆ Feature:Disposable,Eco-Friendly,Stocked
- ◆ The flexible disposable applicator can be used for bonding in the oral cavity, fluoride coating of teeth, saving medicine and cost.
- ◆ The brush head is non-absorbent, lint-free, and maintains an overhanging state to prevent liquid splashing, dripping, or waste.



Syringe Tips

GC-ST02

- ◆ Material:ABS,Plastic Air Water Syringe Tip
- ◆ Color: White,Blue,Yellow,Green,Purple,Orange
- ◆ Disinfecting Method: EO and Steam
- ◆ Packing:250pcs/box,40box/ctn or 250pcs/bag,40bag/ctn
- ◆ Application: Dental,Clinic,Hospital



Interdental Brush

GC-IB

- ◆ Insert the brush into interdental space slowly and scrub back and forth.
- ◆ After using this brush,you can clean it in the water and dry it for next time use.
- ◆ Gently slide between teeth
- ◆ Move around each tooth,above and below the gumline
- ◆ Use daily to help reduce plaque buildup

Dia. NO.	1(SSS)	2(SS)	3(S)	4(M)	5(L)
Shape					
Tapered	0.7MM	0.8MM	1.0MM	1.2MM	1.5MM





Sterile Indicator >>> Control System

Chemical Indicator Tape

Gloria Chemical indicator tapes are a safe solution for closing all wrapping materials. They consist of adhesive, backing and chemical indicator stripes. The resilient, pressure-sensitive adhesive is designed for use with a variety of wraps like disposable non-woven, paper, reusable 100% cotton and cotton/poly blends. After sterilization, the tape is designed to remove easily and cleanly from wraps. The backing is a beige colored crepe paper and the chemical indicator stripes will show a visual color change when exposed to a sterilization process and indicates if the pack has been processed.

The chemical indicator tapes are suitable for steam and ethylene oxide, Plasma and Formaldehyde sterilization processes and also provide a clean peel from the wrapping material.



► Specification Sheet

Chemical Indicator Tapes

Description	Reference No.	Dimension	Rolls/Carton	Color Change
Steam Indicator Tape	GCST-1250	12mm x 50m 0.5"	76rolls/ctn	From Beige to dark brown/black under steam sterilization
	GCST-1550	15mm x 50m 0.6"	60rolls/ctn	
	GCST-1950	19mm x 50m 0.75"	48rolls/ctn	
	GCST-2550	25mm x 50m 1"	36rolls/ctn	
EO Indicator Tape	GCST-1250	12mm x 50m 0.5"	76rolls/ctn	From Red to Green under EO sterilization
	GCST-1550	15mm x 50m 0.6"	60rolls/ctn	
	GCST-1950	19mm x 50m 0.75"	48rolls/ctn	
	GCST-2550	25mm x 50m 1"	36rolls/ctn	
Plasma Indicator Tape	GCST-1250	12mm x 50m 0.5"	76rolls/ctn	From Blue to pink under plasma sterilization
	GCST-1550	15mm x 50m 0.6"	60rolls/ctn	
	GCST-1950	19mm x 50m 0.75"	48rolls/ctn	
	GCST-2550	25mm x 50m 1"	36rolls/ctn	

Chemical Indicator Strips

Gloria chemical indicator strips are used to determine whether the sterilization process conditions were met inside each pack. They consist of special card with a chemical indicator ink. The water based and non-toxic chemical indicator in compliance with ISO 11140-1 will show a clear and accurate color change after successful sterilization.

The chemical indicator strips are suitable for steam and ethylene oxide, Plasma and Formaldehyde sterilization processes.



► Specification Sheet

Chemical Indicator Strips

Description	Reference No.	Dimension	Packing	Color Change
Steam Indicator Tape	GCS-Class 4	15mm x 100mm	200pcs/box 100boxes/ctn	From Blue to Black
	GCS-Class 5	20mm x 108mm	200pcs/box 100boxes/ctn	From Blue to Black
	GCS-Class 6	20mm x 140mm	200pcs/box 100boxes/ctn	From Yellow to Black
EO Indicator Strip	PMDCE	15mm x 100mm	200pcs/box 100boxes/ctn	From Red to Blue
Plasma Indicator Strip	PMDSP	15mm x 100mm	200pcs/box 100boxes/ctn	From Green to Yellow

Sterile Indicator >>>

Control System

Chemical Indicator Tape

Gloria Bowie Dick Test Pack are used in all steam sterilizers for detecting air leak problems, vacuum pump problems, High temperature, wet steam and inadequate steam penetration. It is advised to place packs on the bottom shelf above drain. Bowie Dick Test Packs should be used on daily basis and after major maintenance, relocating sterilizer, installation and after malfunctions. Bowie-Dick Test Packs verifies that the pre-vacuum sterilizer effectively removes air and provides

Efficient steam penetration. If no appreciable air is present in the sterilizer chamber after vacuums are drawn, the steam will penetrate to the middle of the test pack and turn the indicator sheet a uniform black colour.

Bowie-Dick Test Packs are manufactured according to EN ISO 11140-1 standards and Classified as Class 2.



► Specification Sheet

Chemical Indicator Tapes

Item No	Sterilization Method	Pcs/Carton	Color Change
GC-BD	115mm x 125mm	100	From Blue to Black

EO Indicator Label

Scope of application: it is used for ethylene oxide sterilization monitoring, which is used to mark whether the package has passed the process of ethylene oxide treatment. Meanwhile, it is used for monitoring and judging the effect of sterilization indication.

Indicates the label from the backing off, affixed to the sterilized goods package or directly affixed to the sterilized items, then placed in the ethylene oxide sterilization chamber after sterilization, the label changed from red to blue, proved qualified sterilization,

The instruction label only provides information about whether the substance has passed ethylene oxide sterilization, and does not reflect the process and effect of disinfection and sterilization of sterilized articles.

Storage temperature of 15-25 degrees, relative humidity is not more than 50%RH, avoid light, and do not coexist with contaminated or toxic chemicals.



► Specification Sheet

EO Indicator Label

Item No	Sterilization Method	Size	Package	Color Change
GC-EL	EO Gas	30mm Diameter	5,000 pcs/bag	30mm Diameter

Medical >>> Face Mask



European Standard

Classification	TYPE I ^a	TYPE II	TYPE II R
BFE%	≥95	≥98	≥98
Pa/cm2	<40	<40	<60
KPa	/	/	≥16.0
cfu/g	≤30	≤30	≤30

Medical TYPE II /TYPE II R Face Mask


GC-3P02

Packing: 50pcs/box 2000pcs/ctn

Size: 17.5*9.5cm Style: Ear loop

Fabric: PP non woven+Melt blown filter paper+PP non woven

Color: Blue



3 PLY



Protection



Breathable



Disposable

AMERICAN STANDARD

AAMI PB70 Barrier performance	Test Method	Test definition AATCC	Requirement
Level 1	AATCC 42 Impact penetration	AATCC 42 <i>Measures the resistance of fabrics to the liquid penetration of water by impact.</i>	Water impact <4.5g
Level 2	AATCC 42 Impact penetration AATCC 127 Hydrostatic Pressure	AATCC 42 <i>Measures the resistance of fabrics to the liquid penetration of water by impact.</i> AATCC 127 <i>Measures the resistance of fabrics to the liquid penetration of water by impact under constant and increasing hydrostatic pressure.</i>	Spary Impact < 1.0g Hydrostatic Pressure > 20cm
Level 3	AATCC 42 Impact penetration AATCC 127 Hydrostatic Pressure	AATCC 42 <i>Measures the resistance of fabrics to the liquid penetration of water by impact.</i> AATCC 127 <i>Measures the resistance of fabrics to the liquid penetration of water by impact under constant and increasing hydrostatic pressure.</i>	Spary Impact < 1.0g Hydrostatic Pressure > 50cm
Level 4	ASTM F1671 Viral penetration	ASTM F1671 <i>Measures the resistance of fabrics to bloodborne pathogens using viral penetration at 2 psi and ambient pressure</i>	Totally Impervious

KF94 / KF99 Face Mask >>>

KF94 Fish Face Mask with Earloop GC-KF01

Packing: 10pcs/bag, 2500pcs/ctn
 Size: 20.5*8cm Style: Ear loop
 Fabric: PP non woven+Melt blown filter paper
 Color: Any color available



PLEASE CONTACT US FOR MORE STYLES

시험 항목 및 결과

시험항목	시험기준		단위	결과		비고	
성상	흰색의 3단 가로접이식 본체에 코편이 있고, 양 측면에 흰색의 끈이 있는 부직포 마스크		-	흰색의 3단 가로접이식 본체에 코편이 있고, 양 측면에 흰색의 끈이 있는 부직포 마스크		자사기준 및 시험방법	
형상	본체	가로	205 ± 10	mm	본체 가로	206	자사기준 및 시험방법
		세로	162 ± 10		본체 세로	162	
	머리끈	좌	155 ± 10		머리끈 좌	155	
		우	155 ± 10		머리끈 우	153	
순도시험 (색소)	위에서 관찰할 때 색을 나타내지 않는다.		-	위에서 관찰할 때 색을 나타내지 않는다.		의약외품에 관한 기준 및 시험방법	
순도시험 (산 또는 알칼리)	홍색을 나타내지 않는다. / 적색을 나타내지 않는다.		-	홍색을 나타내지 않는다. / 적색을 나타내지 않는다.		의약외품에 관한 기준 및 시험방법	
순도시험 (형광증백제)	형광을 나타내지 않아야한다.		-	형광을 나타내지 않는다.		의약외품에 관한 기준 및 시험방법	
순도시험 (포름알데히드)	검액이 나타내는 색은 비교액이 나타내는 색보다 진하지 않다.		-	검액이 나타내는 색은 비교액이 나타내는 색보다 진하지 않다.		의약외품에 관한 기준 및 시험방법	

KF94 Fish Face Mask with Adjust Earloop

GC-KF03

Packing: 10pcs/bag, 2000pcs/ctn
 Size: 20.5*8cm Style: Adjustable earloop
 Fabric: PP non woven+Melt blown filter paper
 Color: Any color available



KF99 Fish Face Mask with Earloop

GC-KF08

Packing: 10pcs/bag, 2500pcs/ctn
 Size: 20.5*8cm Style: Ear loop
 Fabric: PP non woven+Melt blown filter paper
 Color: Mainly black and white



KF94 Fish Face Mask with Headloop

GC-KF02

Packing: 10pcs/bag, 2000pcs/ctn
 Size: 20.5*8cm Style: Head loop
 Fabric: PP non woven+Melt blown filter paper
 Color: Any color available



FFP2 >>>

Face Mask



Sodium Chloride Testing (sheet):

Penetration of filter material; Sodium Chloride Testing				
Condition	No. of Sample	Sodium Chloride Testing 95 L/min max 0%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
(A.R.)	19	0,15	FFP1 ≤ 20 % FFP2 ≤ 6% FFP3 ≤ 1%	Filtering half masks fulfill the requirements of the standard CE EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1, FFP2 and FFP3 classes.
(A.R.)	20	0,17		
(A.R.)	21	0,15		
(S.W.)	22	0,33		
(S.W.)	23	0,27		
(S.W.)	24	0,31		
(M.S. T.C)	25	0,52		
(MS. T.C)	26	0,59		
(M.S. T.C)	27	0,57		

FFP2 Adult Face Mask (Without Valve)

GC-FP01

Packing: 1 pcs/bag, 20 pcs/box, 1000 pcs/ctn
 Certificate: CE (EN149:2001+A1:2009)
 CE 0370 Size: 15.5*10.5cm
 Fabric: PP non woven+Melt blown filter paper
 Style: Ear loop Color: Any color available



FFP2 Children Face Mask (Without Valve)

GC-FPM05

Packing: 1 pcs/bag, 10 pcs/box, 1200 pcs/ctn
 Certificate: CE (EN149:2001+A1:2009)
 CE 1463 Size: 13.5*10cm
 Fabric: PP non woven+Melt blown filter paper
 Style: Ear loop Color: Any color available



Disposable Sterile Dressing Pack >>

Sterile Dressing Change Kit for Single Use

GC-DP01

Style: Surgical Supplies



Advantages:

1. Reasonable configuration, stable quality and high cost performance;
2. Design of different application scenarios, comprehensive products and more selectivity;
3. high degree of customization, can meet the needs of different markets and clinical.



Accessories:

- 6×Sterile cotton ball
- 1×Plastic tray
- 5×Iodophor cotton ball 0.4g
- 2×Tweezers
- 1×44*32cm surgical pad
- 5×7.5*7.5cm 8 layers of gauze

Disposable Sterile Dressing Pack

GC-DP02

Style: Surgical Supplies

Accessories:

- 6×Iodophor cotton ball 0.4g
- 1×44*32cm surgical pad
- 4×Sterile cotton ball
- 2×Tweezers
- 5×7.5*7.5cm 8 layers of gauze
- 1×Plastic tray



Medical Consumables Medical Grade Dressing Kit

GC-DP03

Style: Surgical Supplies

Accessories:

- 1×44*32cm surgical pad
- 5×Sterile cotton ball
- 1×Tweezers
- 5×7.5*7.5cm 8 layers of gauze



Disposable Sterile Dressing Pack >>

Dressing Change Tray with Chloraprep

GC-DP04




Accessories:
Gauze, Window, Towel, Vinyl gloves, Chloraprep, Mask with earloops, Cotton-tip applicator, Tape

Product name	Medical Consumables Medical Grade Dressing Kit
Material	Polyester / Cotton
size	Customzied
Sterile	EO sterile
Color	Customzied
Certificate	CE ISO
Attention	Disposable ,for single use only

Blood Taking Set >>

Blood Taking Set



Application to treat affected site after blood taking procedure in order to provide healing and prevent further harm to patient.

Accessories:
Plastic Tray with 3-Compartments, Plastic Forceps, Cotton Balls, Dressing Towel, Disposable Bag

Delivery Procedural Set >>

Delivery Procedural Set



A procedural set consist of general medical which prepared to suit delivery procedure of newborn. To facilitate the delivery procedure of newborn.

Accessories:
Plastic Tray, Sterile Field, Patient Drape, Under pad, Gloves, Cotton Balls, Cotton String, Cord Clamps, Amniotic Hook, Dressing Towel, Biohazard Bag, Green Crepe Paper, Umbilical Cord Scissor.