



JIANGSU DEFENG MEDICAL EQUIPMENT CO.,LTD
SUZHOU DEFENG TECHNOLOGY CO.,LTD

DONGFU RD, DONGLAI TOWN, ZHANGJIAGANG CITY, JIANGSU PROVINCE, CHINA 215627

Tel: (+86)-512-58107888 ext. 815

ATTN: Susan Xu

Skype: susanxumei
















Fax: (+86)-512-58107300






Mobile: +86-13915674322

Email: sales3@cndeftech.com

IMAGE	DETAILS	IMAGE	DETAILS
	<p>Name: Emergency stretcher Model: DDJ-1A Feature: single fold, two wheels Material: durable aluminum Cloth: 600D textile with PVC backing</p> <p>Dimension: 185×50×24cm Folded: 92×50×9cm Capacity: 159kg N.W.: 5kg Packing: 95×51×16cm, 2pcs/CTN G.W.: 12kg</p>		<p>Name: Emergency stretcher Model: DDJ-1B Feature: single fold, basic type Material: durable aluminum Cloth: 600D textile with PVC backing</p> <p>Dimension: 185×50×3cm Folded: 92×50×6cm Capacity: 159kg N.W.: 3kg Packing: 95×51×16cm, 2pcs/CTN G.W.: 8kg</p>
	<p>Name: Emergency stretcher Model: DDJ-1C Feature: single fold Material: durable aluminum Cloth: 600D textile with PVC backing</p> <p>Dimension: 185×50×5cm Folded: 93×52×6cm Capacity: 159kg N.W.: 3kg Packing: 97×55×18cm, 2pcs/CTN G.W.: 8kg</p>		<p>Name: Emergency stretcher Model: DDJ-1D Feature: single fold, with backrest Material: durable aluminum Cloth: 600D textile with PVC backing</p> <p>Dimension: 185×54×24cm Folded: 94×54×14cm Capacity: 159kg N.W.: 9kg Packing: 98×54×29cm, 2pcs/CTN G.W.: 18kg</p>
	<p>Name: Emergency stretcher Model: DDJ-1E Feature: chair and stretcher, 2 wheels Material: durable aluminum</p> <p>Cloth: 600D textile with PVC backing</p> <p>Dimension: 190×58×20cm Folded: 105×55×20cm Capacity: 159kg N.W.: 14kg Packing: 107×57×22cm, 1pc/CTN G.W.: 16kg</p>		<p>Name: Disaster stretcher Model: DDJ-1F Feature: easy for stack Material: durable aluminum</p> <p>Cloth: 600D textile with PVC backing</p> <p>Dimension: 187×47×18cm Folded: / Capacity: 159kg N.W.: 5kg Packing: 190×50×5cm, 1pc/CTN G.W.: 6kg</p>
	<p>Name: Scoop stretcher Model: DDJ-3A Feature: lacquered surface two 5" quick-release belts</p> <p>Material: durable aluminum</p> <p>Dimension: 204×44×7cm Folded: 169×44×7cm Capacity: 159kg N.W.: 8kg Packing: 172×48×8cm, 1pc/CTN G.W.: 9kg</p>		<p>Name: Scoop stretcher Model: DDJ-3B Feature: oxydix surface two 5" quick-release belts</p> <p>Material: durable aluminum</p> <p>Dimension: 204×44×7cm Folded: 169×44×7cm Capacity: 159kg N.W.: 8kg Packing: 172×48×8cm, 1pc/CTN G.W.: 9kg</p>
	<p>Name: Scoop stretcher Model: DDJ-3C Feature: lacquered, 3 pieces two 5" quick-release belts</p> <p>Material: durable aluminum</p> <p>Dimension: 204×44×7cm</p> <p>Folded: 169×44×7cm Capacity: 159kg N.W.: 8kg Packing: 172×48×8cm, 1pc/CTN G.W.: 9kg</p>		<p>Name: Scoop stretcher Model: DDJ-3D Feature: plastic two 5" quick-release belts</p> <p>Material: special PE</p> <p>Dimension: 204×44×7cm</p> <p>Folded: 169×44×7cm Capacity: 159kg N.W.: 8kg Packing: 172×48×8cm, 1pc/CTN G.W.: 9kg</p>

	<p>Name Scoop stretcher Model DDJ-3E Feature both scoop and spine board X-ray allowed Material HDPE</p> <p>Dimension 187×42×6cm Folded unfolable Capacity 205kg N.W. 8kg Packing 189×43×8cm, 1pc/CTN G.W. 9kg</p>		<p>Name Scoop stretcher Model DDJ-3F Feature oxydix surface three 5" quick-release belts Material durable aluminum</p> <p>Dimension 204×44×7cm Folded 169×44×7cm Capacity 159kg N.W. 8kg Packing 172×48×8cm, 1pc/CTN G.W. 9kg</p>
	<p>Name Pole stretcher Model DDJ-2A Feature single fold, popular type Material durable aluminum Cloth 600D textile with PVC backing</p> <p>Dimension 205×52×14cm Folded 102×18×11cm Capacity 159kg N.W. 6kg Packing 114×20×24cm, 2pcs/CTN G.W. 13kg</p>		<p>Name Pole stretcher Model DDJ-2B Feature single fold, 4 wheels Material durable aluminum Cloth 600D textile with PVC backing</p> <p>Dimension 205×52×14cm Folded 102×18×11cm Capacity 159kg N.W. 7kg Packing 114×20×24cm, 2pcs/CTN G.W. 15kg</p>
	<p>Name Pole stretcher Model DDJ-2C Feature unfold Material durable aluminum</p> <p>Cloth 600D textile with PVC backing Dimension 205×52×14cm Folded 205×11×14cm Capacity 159kg N.W. 5kg Packing 208×25×30cm, 6pcs/CTN G.W. 32kg</p>		<p>Name Pole stretcher Model DDJ-2D Feature unfold, 4 wheels Material durable aluminum</p> <p>Cloth 600D textile with PVC backing Dimension 205×52×14cm Folded 205×11×14cm Capacity 159kg N.W. 6kg Packing 208×25×30cm, 6pcs/CTN G.W. 38kg</p>
	<p>Name Pole stretcher Model DDJ-2E Feature unfold</p> <p>Material durable aluminum Cloth 600D textile with PVC backing Dimension 205×52×14cm Folded 205×11×14cm Capacity 159kg N.W. 5kg Packing 208×25×30cm, 6pcs/CTN G.W. 32kg</p>		<p>Name Camp cot Model DDJ-9 Feature single fold, carry bag</p> <p>Material heavyduty aluminum Cloth 600D textile Dimension 194×66×41CM Folded 90×12×20CM Capacity 159kgs N.W. 8.5kgs Packing 92×26×15CM 1PC/CTN G.W. 9.5kgs</p>
	<p>Name Pole stretcher Model DDJ-4A Feature double fold, 2 wheels</p> <p>Material durable aluminum Cloth 600D textile with PVC backing Dimension 186×50×22cm Folded 46×50×15cm Capacity 159kg N.W. 7kg Packing 95×51×16cm, 1pc/CTN G.W. 8kg</p>		<p>Name Pole stretcher Model DDJ-4B Feature double fold, basic type</p> <p>Material durable aluminum Cloth 600D textile with PVC backing Dimension 213×52×12cm Folded 58×20×15cm Capacity 159kg N.W. 6kg Packing 60×25×15cm, 1pc/CTN G.W. 7kg</p>

 	<p>Name Army stretcher Model DDJ-4C Feature double fold, army type Material durable aluminum</p> <p>Cloth 600D textile with PVC backing Dimension 213×52×16cm</p> <p>Folded 54×18×20cm Capacity 159kg N.W. 9kg Packing 60×25×15cm, 1pc/CTN G.W. 10kg</p>	 	<p>Name Pole stretcher Model DDJ-4E Feature double fold Material durable aluminum</p> <p>Cloth 600D textile with PVC backing Dimension 210×53×12cm</p> <p>Folded 61×30×12cm Capacity 159kg N.W. 6kg Packing 63×25×18cm, 1pc/CTN G.W. 7kg</p>
 	<p>Name Pole stretcher Model DDJ-4D Feature double fold, deluxe type Material durable aluminum Cloth 600D textile with PVC backing Dimension 215×50×15cm Folded 55×18×20cm</p> <p>Capacity 159kg N.W. 7kg Packing 57×25×20cm, 1pc/CTN G.W. 8kg</p>	 	<p>Name Tactical litter Model DDJ-4F Feature double fold, army type Material durable aluminum Cloth Mesh textile Dimension 215×50×15cm Folded 55×18×20cm</p> <p>Capacity 159kg N.W. 8kg Packing 60×25×15cm, 1pc/CTN G.W. 9kg</p>
	<p>Name Animal Stretcher Model DDJ-2J Feature suitable for animal rescue two 5" quick-release belts Dimension 120×52.2×7cm Folded / Capacity 130kg</p> <p>N.W. 1.5kg Packing 122×54×8cm, 1pc/CTN G.W. 3kg</p>	 	<p>Name transfer board Model DDJ-5 Feature single fold two 5" quick-release belts Material durable aluminum Dimension 185×50×4cm Folded 94×50×8cm</p> <p>Capacity 159kg N.W. 7kg Packing 95×51×26cm, 2pcs/CTN G.W. 15kg</p>
	<p>Name Spine board Model DDJ-6A Material HDPE Feature universal type well for X-ray, MRI, CT</p> <p>two 5" quick-release belts Dimension 183×45×6cm Capacity 159kg N.W. 8kg Packing 188×48×7cm, 1pc/CTN G.W. 9kg</p>		<p>Name Spine board Model DDJ-6B Material HDPE Feature individual leg design well for X-ray, MRI, CT</p> <p>two 5" quick-release belts Dimension 183×44×5cm Capacity 159kg N.W. 8kg Packing 188×48×7cm, 1pc/CTN G.W. 9kg</p>
	<p>Name Spine board Model DDJ-6C Material HDPE Feature small width well for X-ray, MRI, CT two 5" quick-release belts</p> <p>Dimension 183×41×5cm Capacity 159kg N.W. 8kg Packing 188×48×6cm, 1pc/CTN G.W. 9kg</p>		<p>Name Spine board Model DDJ-6D Material HDPE Feature pediatric type well for X-ray, MRI, CT two 5" quick-release belts</p> <p>Dimension 138×35×5cm Capacity 159kg N.W. 4kg Packing 140×37×6cm, 1pc/CTN G.W. 5kg</p>













	<p>Name Spine board Model DDJ-6E Material LLDPE Feature made from rotation molding well for X-ray, MRI, CT two 5" quick-release belts</p> <p>Dimension 183×46×6cm</p> <p>Capacity 159kg N.W. 6kg Packing 188×48×7cm, 1pc/CTN G.W. 7kg</p>		<p>Name Spine board Model DDJ-6F Material LLDPE Feature strong but lightweight well for X-ray, MRI, CT two 5" quick-release belts</p> <p>Dimension 183×45×4cm</p> <p>Capacity 159kg N.W. 6kg Packing 188×48×7cm, 1pc/CTN G.W. 7kg</p>
	<p>Name transfer board Model DDJ-6K Feature single fold two 5" quick-release belts</p> <p>Material engineering plastic</p> <p>Dimension 188×50×4.5cm Folded 94×50×9cm</p> <p>Capacity 159kg N.W. 9kg Packing 96×52×10cm, 1pc/CTN G.W. 11kg</p>		<p>Name Basket stretcher Model DDJ-7E Material carton steel Feature one section, black color all accessories included</p> <p>Dimension 212*62*18cm Folded /</p> <p>Capacity 272kgs N.W. 13kgs Packing 215*64*22cm 1PC/CTN G.W. 15kgs</p>
	<p>Name Basket stretcher Model DDJ-7A Material PE Feature two sections all accessories included</p> <p>Dimension 219*64*18cm Folded 130*67*29cm</p> <p>Capacity 272kgs N.W. 20kgs Packing 130*68*31cm 1PC/CTN G.W. 24kgs</p>		<p>Name Basket stretcher Model DDJ-7B Material PE Feature one section all accessories included</p> <p>Dimension 221*63*19cm Folded /</p> <p>Capacity 159kgs N.W. 18kgs Packing 225*65*24cm 1PC/CTN G.W. 20kgs</p>
	<p>Name Basket stretcher Model DDJ-7C Material Stainless steel Feature one section all accessories included</p> <p>Dimension 212*62*18cm Folded /</p> <p>Capacity 272kgs N.W. 13kgs Packing 215*64*22cm 1PC/CTN G.W. 15kgs</p>		<p>Name Basket stretcher Model DDJ-7D Material Stainless steel Feature two sections all accessories included</p> <p>Dimension 212*62*18cm Folded 130*62*27cm</p> <p>Capacity 272kgs N.W. 13kgs Packing 130*68*31cm 1PC/CTN G.W. 15kgs</p>
	<p>Name Soft stretcher Model DDJ-8A Feature 6 handles blue color</p> <p>Cloth 600D textile</p> <p>Dimension 163*65cm Folded yes</p> <p>Capacity 159kgs N.W. 2kgs Packing 89*25*40cm 10PC/CTN G.W. 22kgs</p>		<p>Name Soft stretcher Model DDJ-8B</p> <p>Feature</p> <p>Material</p> <p>Cloth</p> <p>Dimension 198*97cm Folded /</p> <p>Capacity 159kgs N.W. 8.5kgs Packing 60*30*30cm 5PC/CTN G.W. 9.5kgs</p>
	<p>Name Soft stretcher Model DDJ-8C Feature 10 handles yellow color</p> <p>Material Nylon cloth</p> <p>Dimension 212*75cm Folded 88*24cm</p> <p>Capacity 100kg N.W. 2kg Packing 66*44*40cm, 10pcs/CTN G.W. 28kg</p>		<p>Name Full-body vacuum mattress Model DDJ-10 Feature fit for X-ray with carry bag and pump</p> <p>Material TPU cloth</p> <p>Dimension 212×75×7CM Folded 70×60×37CM</p> <p>Capacity 159kgs N.W. 5kgs Packing 95×28×58CM G.W. 7kgs</p>

	<p>Name Military stretcher Model DDJ-11A(DDJ-2H) Feature unfoldable Material durable aluminum Cloth camouflage color canvas cloth Dimension 220*54*15cm Folded / Capacity 159kg</p> <p>N.W. 8kg Packing 191*17*13cm, 1pc/CTN G.W. 9kg</p>		<p>Name Military stretcher Model DDJ-11F(DDJ-4F) Feature double fold Material durable aluminum Cloth camouflage color canvas cloth Dimension 188*45*4cm Folded 46*45*15cm Capacity 159kg</p> <p>N.W. 6kg Packing 55*18*52cm, 1pc/CTN G.W. 7kg</p>
	<p>Name Stair stretcher Model DFT-1A Feature foldable, 2 wheels Material durable aluminum Cloth plastic Dimension 69*51*90CM Folded 18*51*90CM</p> <p>Capacity 159kgs N.W. 11kgs Packing 20*52*92CM 1PC/CTN G.W. 13kgs</p>		<p>Name Stair stretcher Model DFT-1B Feature 4 wheels Material durable aluminum Cloth plastic Dimension 69*51*90CM Folded 18*51*90CM</p> <p>Capacity 159kgs N.W. 11kgs Packing 20*52*92CM 1PC/CTN G.W. 13kgs</p>
	<p>Name Stair stretcher Model DFT-1C Feature foldable, 2 wheels Material durable aluminum Seat aluminum board Dimension 90*51*90CM Folded 90*17*59CM</p> <p>Capacity 159kgs N.W. 11kgs Packing 96*17*64CM 1PC/CTN G.W. 12kgs</p>		<p>Name Stair stretcher Model DFT-2C Feature foldable, 4 wheels Material durable aluminum Cloth plastic Dimension 69*52*90CM Folded 18*54*90CM</p> <p>Capacity 159kgs N.W. 9.5kgs Packing 20*55*92CM 1PC/CTN G.W. 11.5kgs</p>
	<p>Name Stair stretcher Model DFT-2D Feature foldable, 4 wheels Material durable aluminum Cloth plastic Dimension 80*50*160CM Folded 80*50*28CM</p> <p>Capacity 159kgs N.W. 10kgs Packing 82*52*30CM 1PC/CTN G.W. 12kgs</p>		<p>Name Stair stretcher Model DFT-3D Feature foldable, 4 wheels Material durable aluminum Cloth plastic Dimension 80*50*160CM Folded 80*50*28CM</p> <p>Capacity 159kgs N.W. 10kgs Packing 82*52*30CM 1PC/CTN G.W. 12kgs</p>
	<p>Name Stair stretcher Model DFT-2E Feature go downstairs Material stainless steel Cloth plastic Dimension 131*50*155CM Folded 15*51*30CM</p> <p>Capacity 159KGS N.W. 15KGS Packing 108*54*33CM 1PC/CTN G.W. 17KGS</p>		<p>Name Stair stretcher Model DFT-3E Feature electric, go upstairs & downstairs Material stainless steel Cloth plastic Dimension 131*50*155CM Folded 15*51*30CM</p> <p>Capacity 159KGS N.W. 20KGS Packing 108*54*33CM 1PC/CTN G.W. 22KGS</p>
	<p>Name Ambulance stretcher Model DDC-1 Feature low position 5cm mattress Material durable aluminum Dimension 190*55*54CM Folded /</p> <p>Capacity 159kgs N.W. 23kgs Packing 195*56*35CM 1PC/CTN G.W. 28kgs</p>		<p>Name Ambulance stretcher Model DDC-2 Feature low position 5cm mattress Material durable aluminum Dimension 190*55*50CM Folded /</p> <p>Capacity 159kgs N.W. 30kgs Packing 195*56*35CM 1PC/CTN G.W. 35kgs</p>

	<p>Name: Ambulance stretcher Model: DDC-3 Feature: high-low position 5cm mattress Material: durable aluminum Dimension: 193*57*(27-93)cm Folded: / Capacity: 159kgs N.W.: 25kgs Packing: 195*65*27CM 1PC/CTN G.W.: 28kgs</p>		<p>Name: Ambulance stretcher Model: DDC-3B Feature: high-low position 5cm mattress Material: durable aluminum Dimension: 193*57*(27-90)cm Folded: / Capacity: 159kgs N.W.: 38kgs Packing: 196*65*27CM 1PC/CTN G.W.: 44kgs</p>
	<p>Name: Ambulance stretcher Model: DDC-3C Feature: high-low position 5cm mattress Material: durable aluminum Dimension: 193*57*(27-90)cm Folded: / Capacity: 159kgs N.W.: 23kgs Packing: 196*65*27CM 1PC/CTN G.W.: 35kgs</p>		<p>Name: Ambulance stretcher Model: DDC-4 Feature: high-low position 5cm mattress Material: durable aluminum Dimension: 197*57*(22-90)cm Folded: / Capacity: 159kgs N.W.: 35kgs Packing: 205*64*27CM 1PC/CTN G.W.: 40kgs</p>
	<p>Name: Ambulance stretcher Model: DDC-5 Feature: height adjusted 5cm mattress Material: durable aluminum Dimension: 198*53*(34-98)cm Folded: / Capacity: 159kgs N.W.: 35kgs Packing: 200*55*35CM 1PC/CTN G.W.: 45kgs</p>		<p>Name: Ambulance stretcher Model: DDC-6 Feature: high-low position 5cm mattress Material: durable aluminum Dimension: 198*53*(34-98)cm Folded: / Capacity: 159kgs N.W.: 40kgs Packing: 200*55*35CM 1PC/CTN G.W.: 45kgs</p>
	<p>Name: Stretcher base Model: DDC-10 Feature: basic type can put scoop stretcher below Material: stainless steel Cloth: stainless steel Dimension: 203*67*40cm Folded: / Capacity: 220kg N.W.: 27kg Packing: 238*68*18cm, 1pc/CTN G.W.: 30kg</p>		
	<p>Name: First-aid stretcher Model: DQC-1 Feature: with I.V. pole top stretcher can be separated Material: durable aluminum Cloth: plastic PVC Dimension: 192*56*(50-88)cm Capacity: 159kgs N.W.: 30kgs Packing: 194*56*31CM 1PC/CTN G.W.: 35kgs</p>		<p>Name: First-aid trolley Model: DQC-2 Feature: backrest adjustable top stretcher can be separated Material: durable aluminum Cloth: plastic PVC Dimension: 192*56*(50-88)cm Capacity: 159kgs N.W.: 30kgs Packing: 194*56*31CM 1PC/CTN G.W.: 35kgs</p>
	<p>Name: Steel trolley Model: DQC-11 Material: durable steel Feature: two wheels with brake easy to operate Dimension: 122*50*55CM Folded: 8*60*25CM Capacity: 300kgs N.W.: 13kgs Packing: 10*64*27CM 1PC/CTN G.W.: 15kgs MOQ: 20PCS Lead time: 15 days</p>		<p>Name: Coffin trolley Model: DQC-12 Material: durable stainless steel Feature: flexible frame with 4 carrying handles Dimension: 145*80*61CM Folded: 56*21*73CM Capacity: 450kgs N.W.: 16kgs Packing: 57*58*21CM 1PC/CTN G.W.: 18kgs MOQ: 20PCS Lead time: 15 days</p>







	<p>Name Coffin trolley Model DQC-13 Material durable steel with epoxy coated Feature flexible frame with 4 carrying handles Dimension 145*80*61CM Folded 56*21*73CM Capacity 450kgs N.W. 16kgs Packing 57*58*21CM 1PC/CTN G.W. 18kgs MOQ 20PCS Lead time 15 days</p>		<p>Name Coffin trolley Model DQC-14 Material durable aluminum alloy Feature flexible frame with 4 carrying handles Dimension 145*80*61CM Folded 56*21*73CM Capacity 450kgs N.W. 16kgs Packing 57*58*21CM 1PC/CTN G.W. 18kgs MOQ 20PCS Lead time 15 days</p>
	<p>Name Mortuary stretcher Model DQC-15 Feature low-middle-high positions 125mm diameter castors Material durable aluminum Dimension 195*56*80CM Folded 195*56*28CM Capacity 250kgs N.W. 26kgs Packing 197*58*30CM 1PC/CTN G.W. 30kgs MOQ 20PCS Lead time 15 days</p>		<p>Name Basic embalming table Model DQC-16 Feature 150mm diameter castors require easy assembly Material Stainless steel 301 Dimension 205*75*85CM Folded unfordable Capacity 450kgs N.W. 37kgs Packing 195*56*35CM 1PC/CTN G.W. 40kgs MOQ 75 PCS Lead time 20 days</p>
	<p>Name Folding embalming table Model DQC-17 Feature 150mm diameter castors easy for storage & delivery Material Stainless steel 301 Dimension 205*75*85CM Folded 205*75*25CM Capacity 450kgs N.W. 37kgs Packing 195*56*35CM 1PC/CTN G.W. 40kgs MOQ 75 PCS Lead time 20 days</p>		<p>Name Casket Lowering Device Model DQC-201 Material Stainless Steel Style American Style Application Adult Dimension 62" x 25" to 94" x 38" Loading capacity 500KG N.W. 75kgs 3 Packing 152*18*10cm ;96*21*21cm ;96*21*21cm G.W. 80kgs MOQ 5pcs Lead time 3-10 days</p>
	<p>Name I.V. stand Model DDC-102 Material stainless steel hook and tube powder coated steel base Holder 2 and 4 hooks available Height 123-230cm adjustable Diameter 26cm/19cm of tube Castor 5cm diameter N.W. 4kgs Packing 130*20*40CM 10PCS/CTN G.W. 45kgs MOQ 100 PCS Lead time 20 days</p>		<p>Name I.V. stand Model DDC-103 Material stainless steel hook and tube chromium coated steel base Holder 2 and 4 hooks available Height 125-230cm adjustable Diameter 26cm/19cm of tube Castor 5cm diameter N.W. 4kgs Packing 130*20*40CM 10PCS/CTN G.W. 45kgs MOQ 100 PCS Lead time 20 days</p>
	<p>Name I.V. stand Model DDC-104 Material stainless steel tube plastic ABS hook and base Holder 2 and 4 hooks available Height 120-224cm adjustable Diameter 26cm/19cm of tube Castor 5cm diameter N.W. 4kgs Packing 130*20*40CM 10PCS/CTN G.W. 45kgs MOQ 100 PCS Lead time 20 days</p>		<p>Name I.V. stand Model DDC-105 Material aluminum hook and tube plastic ABS base Holder 2 and 4 hooks available Height 125-230cm adjustable Diameter 26cm/19cm of tube Castor 5cm diameter N.W. 4kgs Packing 130*20*40CM 10PCS/CTN G.W. 45kgs MOQ 100 PCS Lead time 20 days</p>
	<p>Name I.V. stand Model DDC-106 Material stainless steel hook and tube powder coated steel base Holder 2 and 4 hooks available Height 118-205cm adjustable Diameter 27cm/22cm of tube Castor 5cm diameter N.W. 8kgs Packing 130*20*40CM 10PCS/CTN G.W. 85kgs MOQ 100 PCS Lead time 20 days</p>		<p>Name I.V. stand Model DDC-107 Material aluminum hook and tube aluminum base Holder 2 and 4 hooks available Height 120-215cm adjustable Diameter 27cm/22cm of tube Castor 5cm diameter N.W. 4kgs Packing 130*20*40CM 10PCS/CTN G.W. 45kgs MOQ 100 PCS Lead time 20 days</p>

	<p>Name I.V. pole Model DDC-108 Material stainless steel hook and tube powder coated steel base Holder 2 and 4 hooks available Height 118-205cm adjustable Diameter 27cm/22cm of tube Castor 5cm diameter</p> <p>N.W. 8kgs Packing 130*20*40CM 10PCS/CTN G.W. 85kgs MOQ 100 PCS Lead time 20 days</p>		<p>Name I.V. pole Model DDC-109 Material aluminum hook and tube aluminum base Holder 2 and 4 hooks available Height 120-215cm adjustable Diameter 27cm/22cm of tube Castor 5cm diameter</p> <p>N.W. 4kgs Packing 130*20*40CM 10PCS/CTN G.W. 45kgs MOQ 100 PCS Lead time 20 days</p>
	<p>Name restraint strap Model DDC-201</p> <p>MOQ 100 PCS Lead time 20 days</p>		<p>Name restraint strap Model DDC-202</p> <p>MOQ 100 PCS Lead time 20 days</p>
	<p>Name restraint strap Model DDC-203</p> <p>MOQ 100 PCS Lead time 20 days</p>		<p>Name restraint strap Model DDC-204</p> <p>MOQ 100 PCS Lead time 20 days</p>
	<p>Name Spider strap Model DDC-205</p> <p>MOQ 100 PCS Lead time 20 days</p>		<p>Name Spider strap Model DDC-206</p> <p>MOQ 100 PCS Lead time 20 days</p>
	<p>Name first aid kit Model DF-010</p> <p>MOQ 50 PCS</p>		<p>Name first aid kit Model DF-014</p> <p>MOQ 50 PCS</p>
	<p>Name first aid kit Model DF-015</p> <p>MOQ 50 PCS</p>		<p>Name first aid kit Model DF-016</p>

	<p>Name Head immobilizer Model DTB-1 Feature adult type fit most backboards compatible with X-ray, CT, MRI Material polyethylene & PU foam Color red/orange/yellow/blue Dimension 43.5*27*14CM N.W. 0.9kg Packing 55*46*56CM 8 pcs/CTN G.W. 10kgs MOQ 80pcs Lead time 10 days</p>		<p>Name Head immobilizer Model DTB-2 Feature pediatric type fit most backboards compatible with X-ray, CT, MRI Material polyethylene & PU foam Color red/orange/yellow/blue Dimension 35*28*10CM N.W. 0.5kg Packing 72*58*36CM 12 pcs/CTN G.W. 8kgs MOQ 120pcs Lead time 10 days</p>
	<p>Name Cervical collar Model DJT-1A (adult) DJT-1B (child) Material Special PE Feature Universal type 4 different size adjustable fit X-ray, CT and MRI Adult and child model N.W. 0.2kg Packing 60*50*42cm 40pcs/CTN G.W. 12kg MOQ 480pcs Lead time 10 days</p>		<p>Name Cervical collar Model DJT-2 Feature Philadelphia type high stabilizing effect water resistant maneuverable chine light weight Material ployethylene foam & rigid plastic Size S/M/L Packing 52*46*40cm 480PCS/CTN G.W. 10kgs MOQ 800pcs Lead time 10 days</p>
	<p>Name Cervical collar Model DJT-3 Material EVA Feature length adjustable Packing 50*50*60cm, 100pcs/CTN G.W. 20kg MOQ 1000pcs Lead time 10 days</p>		<p>Name Cervical collar Model DJT-4S DJT-4M DJT-4L Material foam & polyethylene coating Feature washable, non toxic with velcro straps X-ray translucent Length 51cm/56cm/62cm Packing 60*50*42cm, 80pcs/CTN MOQ 1000pcs Lead time 10 days</p>
	<p>Name Liquid silicone resuscitator Model DMR-1A (adult) DMR-1B (child) DMR-1C (infant) Face mask material: liquid silicone Resuscitator bag: 1650ml/600ml/300ml Reservoir bag: 2000ml/1600ml/1600ml Oxygen tube: PVC/200cm PP box 28*16*15cm Packing 57*34*46cm, 12pcs/CTN G.W. 12kg MOQ 480pcs Lead time 10-15 days</p>		<p>Name Solid silicone resuscitator Model DMR-2A (adult) DMR-2B (child) DMR-2C (infant) Face mask material: solid silicone Resuscitator bag: 1650ml/600ml/300ml Reservoir bag: 2000ml/1600ml/1600ml Oxygen tube: PVC/200cm PP box 28*16*15cm Packing 57*34*46cm, 12pcs/CTN G.W. 12kg MOQ 480pcs Lead time 10-15 days</p>
	<p>Name SEBS resuscitator Model DMR-3A (adult) DMR-3B (child) DMR-3C (infant) Face mask material: SEBS Resuscitator bag: 1650ml/600ml/300ml Reservoir bag: 2000ml/1600ml/1600ml Oxygen tube: PVC/200cm PP box 28*16*15cm Packing 57*34*46cm, 12pcs/CTN G.W. 12kg MOQ 480pcs Lead time 10-15 days</p>		<p>Name PVC resuscitator Model DMR-4A (adult) DMR-4B (child) DMR-4C (infant) Face mask material: PVC Resuscitator bag: 1650ml/600ml/300ml Reservoir bag: 2000ml/1600ml/1600ml Oxygen tube: PVC/200cm PP box 28*16*15cm Packing 57*34*46cm, 12pcs/CTN G.W. 12kg MOQ 480pcs Lead time 10-15 days</p>
	<p>Name CPR face shield Model DDJ-100 Feature one-way valve included with instrctions inside individually packed in zip bag Material PP & PVC Dimension 20*20cm N.W. 9kg Packing 61*45*40cm, 500pcs/CTN G.W. 10kg MOQ 10000pcs Lead time 10 days</p>		<p>Name Pocket CPR mask Model DDJ-101 Material PVC for inner CPR maask & gloves PP for hard cover Feature clear, non toxic different color PP box PVC gloves included alcohol pad included N.W. 8kg Packing 45*41*32cm, 60pcs/CTN MOQ 2000pcs Lead time 15 days</p>

	<p>Name Mylar rescue blanket Model EB150418.1 Material Polyethylene Terephthalate Feature silver/silver color up to 90% effective reflection a variety of applications Dimension 78"52"12 ± m Packing Detail Routine manual folding, self-styled bag, containing a 1 color paper card Folded Size 8*10cm Packing 39*35*28cm, 250pcs/CTN G.W. 13 kg/CTN Lead time 25 days MOQ 5000pcs</p>		<p>Name Mylar rescue blanket Model EB150418.1 Material Polyethylene Terephthalate Feature gold/silver color up to 90% effective reflection a variety of applications Dimension 78"52"12 ± m Packing Detail Routine manual folding, self-styled bag, containing a 1 color paper card Folded Size 8*10cm Packing 39*35*28cm, 250pcs/CTN G.W. 13 kg/CTN Lead time 25 days Quantity 5000pcs</p>
	<p>Name Poly-tissue rescue blanket Model EB16001 Material poly-tissue laminate Feature disposable foam layer retains body heat bright yellow for maximum safety Dimension 56" X 90" (142cmX229cm) Folded 20X30cm N.W. 0.4kg Packing 24pcs/CTN G.W. 12kg MOQ 5000pcs Lead time 15 days</p>		<p>Name Poly-foam rescue blanket Model EB17001 Material Poly-foam Feature single use keep patient warm and dry bright yellow for maximum safety Dimension 58"X90"(147cmX229cm) Folded 20X30cm N.W. 0.4kg Packing 18pcs/CTN G.W. 10kg MOQ 4000pcs Lead time 15 days</p>
	<p>Name CPR board (cardiac board) Model DDJ-102 Feature blue/red color Material polypropylene plastic Color Blue/Orange Dimension 60*42*8cm N.W. 1.3kg Packing 61*60*45cm, 10pcs/CTN G.W. 15kg MOQ 1000pcs Lead time 10 days</p>		<p>Name manual aspirator Model DDJ-103 Material PP & PE & PVC Feature portable, lightweight Vaccum -600mmHG, 185*64*168mm Flow 20L/min, MAX N.W. 230g Packing 52*43*40cm, 24pcs/CTN G.W. 11kg MOQ 1000pcs Lead time 10 days</p>
	<p>Name Electric suction unit Model DDJ-104 Voltage AC 220V 50HZ, DC 12V Power 44VA Pressure ≥0.075MPa Pumping ≥15L/min Reservoir 1000ml Noise ≤65dB (A) N.W 6kg Packing 28*20*29cm, 1pc/CTN G.W. 7kg MOQ 500pcs Lead time 30 days</p>		<p>Name Nebulizer Model DDJ-105 Voltage 110V/220V/230V, 50/60HZ, 1.5A power <140VA Pressure Max 2.2bar Air flow max 12L/min Noise 50dB(A) Efficiency 0.25ml/min - 0.5ml/min N.W. 3kg Packing 70*36*23cm, 5pcs/CTN G.W. 16kg MOQ 100pcs Lead time 30days</p>
	<p>Name Berman airway(white) Model DDJ-106 Material soft LDPE Size 4#, 5#, 6#, 7#, 8#, 9#, 10#, 11# Feature disposable 100% latex free available in 8 sizes Packing 46*36*30cm, 400-6000pcs/CTN G.W. 10kg MOQ 3000pcs Lead time 10-15 days</p>		<p>Name Berman airway(colorful) Model DDJ-107 Material soft LDPE Size 4#, 5#, 6#, 7#, 8#, 9#, 10#, 11# Feature disposable 100% latex free available in 8 sizes Packing 46*36*30cm, 400-6000pcs/CTN G.W. 10kg MOQ 5000pcs Lead time 10-15 days</p>
	<p>Name Guedel airway Model DDJ-108 Material semirigid polyethylene Size 4#, 5#, 6#, 7#, 8#, 9#, 10#, 11#, 12# Feature non-toxic 100% latex free available in 9 sizes Packing 46*36*30cm, 450-4000pcs/CTN G.W. 10kg MOQ 3000pcs Lead time 10-15 days</p>		<p>Name Non-rebreathing mask Model DDJ-109A(Adult) DDJ-109B(Child) Material medical PVC Feature adjustable nose clip high concentration available with anti-crush 7ft tubing Packing G.W. MOQ 5000pcs Lead time 10-15 days</p>

	<p>Name Simple oxygen mask Model DDJ-110 Material medical PVC Feature adjustable nose clip high concentration available with anti-crush 7ft tubing</p> <p>Packing G.W.</p> <p>MOQ 5000pcs Lead time 30 days</p>		<p>Name Air splint Model DSP-4 Feature 6 sets, with pump Hand & Wrist 15" Half arm 25" Full arm 32" Foot & Ankle 15" Half leg 25" Full leg 32"</p> <p>Packing 55*45*28cm, 10pcs/CTN G.W. 14kg MOQ 100pcs Lead time 10 days</p>
	<p>Name Extrication kits (Body splint set) Model DSP-1 Feature with carry bag Material oxford fabric & bamboo Dimension 82*77CM N.W. 4kg Packing 86*27*49cm 4pcs/CTN G.W. 19.6kgs</p> <p>MOQ 40pcs Lead time 10 days</p>		<p>Name Leg traction splint set Model DSP-2A(adult) DSP-2B(child) Feature with carry bag Dimension 136*22*18cm(adult) 118*18*17cm(child) Folded 90*22*18cm(adult) 80*18*17cm(child) N.W. 2/1.5kg Packing 98(91)*23*20CM 1pc/CTN G.W. 10/9kg MOQ 4PCS Lead time 10 days</p>
	<p>Name Sam splint - original Model DSP-3A Dimension 4.25" X 36" Material IXPE Feature Adaptable, light weight non-toxic, X-ray translucent aluminum thickness 1.2mm Color orange/blue Weight 4 oz Packing 60*50*37cm, 126roll/CTN G.W. 22kg MOQ 3000pcs Lead time 10 days</p>		<p>Name Sam splint - blank Model DSP-3B Dimension 4.25" X 36" Material IXPE Feature Adaptable, light weight non-toxic, X-ray translucent aluminum thickness 1.2mm Color orange/grey Weight 4 oz Packing 60*50*37cm, 126roll/CTN G.W. 22kg MOQ 3000pcs Lead time 10 days</p>
	<p>Name Sam splint - extra wide Model DSP-3C Dimension 5.5" X 36" Material IXPE Feature Adaptable, light weight non-toxic, X-ray translucent aluminum thickness 1.2mm Color orange/blue Weight 5.9 oz Packing 60*50*37cm, 140pcs/CTN G.W. 22kg MOQ 3000pcs Lead time 10 days</p>		<p>Name Sam splint - Junior Model DSP-3D Dimension 4.25" X 18" Material IXPE Feature Adaptable, light weight non-toxic, X-ray translucent aluminum thickness 1.2mm Color orange/grey Weight 2.4 oz Packing 60*50*37cm, 280pcs/CTN G.W. 22kg MOQ 3000pcs Lead time 10 days</p>
	<p>Name Sam splint Model DSP-3E Dimension 4.25" X 9" Material IXPE Feature Adaptable, light weight non-toxic, X-ray translucent aluminum thickness 1.2mm Color orange/blue Weight 1.1 oz Packing 60*50*37cm, 500pcs/CTN G.W. 22kg MOQ 5000pcs Lead time 10 days</p>		<p>Name Sam splint - finger Model DSP-3F Dimension 4.25" X 1.96" Material IXPE Feature Adaptable, light weight non-toxic, X-ray translucent aluminum thickness 1.2mm Color orange/blue Weight 0.2 oz Packing 60*50*37cm, 3000pcs/CTN G.W. 20kg MOQ 5000pcs Lead time 10 days</p>
	<p>Name Stack finger splint Model DSP-5 Material high quality rigid plastic Size #1 L 1.95", CIRC 2.35" Size #2 L 2.05", CIRC 2.45" Size #3 L 2.25", CIRC 2.65" Size #4 L 2.4", CIRC 2.75" Size #5 L 2.4", CIRC 2.85" Size #5.5 L 2.4", CIRC 3.35" Size #6 L 2.45", CIRC 3.6" Size #7 L 2.45", CIRC 3.9" MOQ 3000pcs Lead time 10 days</p>		<p>Name Baseball Finger Splint Model DSP-6A Material malleable aluminum with foam padding Size S, 8.9cm M, 10cm L, 12.7cm MOQ 1000pcs Lead time 10 days</p>

	<p>Name Fold Over Finger Splint Model DSP-6B Material malleable aluminum with foam padding Size XS, 2.5cm S, 3.8cm M, 5.7cm L, 7.6cm MOQ 1000pcs Lead time 10 days</p>		<p>Name Four Prong Finger Splint Model DSP-6C Material malleable aluminum with foam padding Size S, 3.8cm M, 7.6cm L, 8.9cm MOQ 1000pcs Lead time 10 days</p>
	<p>Name Curved Finger Splint Model DSP-6D Material malleable aluminum with foam padding Size S, 3.8cm M, 7.6cm L, 15cm XL, 23cm MOQ 1000pcs Lead time 10 days</p>		<p>Name Frog Finger Splint Model DSP-6E Material malleable aluminum with foam padding Size S, 3.8cm M, 7.6cm L, 8.9cm MOQ 1000pcs Lead time 10 days</p>
	<p>Name Metal splint Model DSP-7 Material lightweight aluminum Size 9.5*60cm Feature X-ray transparent easy to clean cut and shape to size Packing 52*45*25cm, 400pcs/CTN G.W. 30kg MOQ 4000pcs Lead time 10 days</p>		<p>Name Hand hold pulse oxymeter Model DDJ-121 Display LCD mode Spo2 1-100%, accuracy in 80-99% Pulse 30-250 BPM Power two AA batteries Feature weak infusion 0.5% minimum 2013 new product N.W. 150g Packing 10*7*4CM 1 pc/CTN G.W. 180g MOQ 20pcs Lead time 10 days</p>