



SUPER UNION GROUP

YANGZHOU SUPER UNION IMPORT AND EXPORT CO., LTD.

Address: Room 2003 - 2004 Huacheng Building, NO.456 Wenchang West Road, Yangzhou city, Jiangsu province, China.
Post Code: 225009
Tel: +86 514 80783522
Fax: +86 514 80783519
Mobile: +86 136 0144 3135
Email: sales@ysumed.com; info@ysumed.com
Website: www.yzsumed.com

YANGZHOU SUPER UNION MEDICAL MATERIAL CO., LTD.

Products: Gauze sponge, gauze bandage, gauze roll, lap sponge, paraffin gauze, gauze ball and surgical dressings/first aid kit

JIANGSU HEALTHMED PRODUCTS CO., LTD.

Products: High elastic bandage, PBT bandage, Crepe bandage, POP bandage, Undercast padding, Selvage gauze bandage and surgical dressings

JIANGSU BEST MEDICAL PRODUCTS CO., LTD.

Products: Cotton roll, zigzag cotton, dental cotton roll, cotton ball, cotton swab and surgical dressings

JIANGSU UNION MEDICAL PRODUCTS CO., LTD.


Products: Zinc oxide plaster, silk tape, non woven tape, PE tape, band aid, non woven & transparent wound dressing, eye pad, alcohol swab and surgical dressings

JIANGSU AOSPITAL NON WOVEN PRODUCTS CO., LTD.

Products: Face mask, shoe cover, bouffant cap, doctor cap, clip cap, surgical gown, patient gown, bed sheet, bed cover and non woven products

JIANGSU SOUTHMED PRODUCTS CO., LTD.

Products: Syringe, needles, infusion sets, transfusion sets, surgical sutures, scalp vein sets and disposable products



COMPANY INTRODUCTION

Super Union Group is a professional supplier of medical product development, covering thousands of products in the medical field. We have our own factory that specialized in manufacturing gauze, cotton, non woven products, all kinds of plasters, bandages, tapes and the other medical products.

Based on our principles of honesty and joint venture with our customers, our company has been constantly expanding to take a leading position in the medical industry, our high efficient team had developed new products every year, thus maintaining the fast growth trend of the company, in order to enhance our management level, as well as to ensure that such high level products in the medical industry to meet international quality requirements. We have gotten ISO13485, CE, FDA and SA8000.

We also manufacture according to your specified requirements and provide OEM and ODM service.

Super Union have won an excellent reputation for quality, efficiency and competitive price. Presently, we have established business relationship with hundreds of clients from over seventy countries and regions in the world, and we hope to do business with you in the near future.



15
MEDICAL
DEVICE



Sterile Gauze Swab

- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 22,20,17,15,13,12,11 threads etc
- ◆ With or without x-ray detectable threads
- ◆ Folded edge or non - folded edge
- ◆ 8ply, 12ply, 16ply, 24ply, 32ply etc
- ◆ 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm, 10x40cm, 20x20cm, 30x30cm etc
- ◆ 1's, 2's, 3's, 5's, 10's packed into pouch (Sterile)
- ◆ Box: 100, 50, 25, 10, 4 pouches/box
- ◆ Pouch: paper+paper, paper+film
- ◆ Gamma, EO, Steam

08/40S 24X20MESH, FOLDED EDGE, 10PCS/POUCH, 10POUCHES/BOX, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17F4816-10S	4**8"-16ply	52*28*46cm	80pouches
SA17F4416-10S	4**4"-16ply	55*30*46cm	160pouches
SA17F3316-10S	3**3"-16ply	53*28*46cm	200pouches
SA17F2216-10S	2**2"-16ply	43*39*46cm	400pouches
SA17F4812-10S	4**8"-12ply	52*28*42cm	80pouches
SA17F4412-10S	4**4"-12ply	55*30*42cm	160pouches
SA17F3312-10S	3**3"-12ply	53*28*42cm	200pouches
SA17F2212-10S	2**2"-12ply	43*39*42cm	400pouches
SA17F4808-10S	4**8"-8ply	52*28*32cm	80pouches
SA17F4408-10S	4**4"-8ply	55*30*32cm	160pouches
SA17F3308-10S	3**3"-8ply	53*28*32cm	200pouches
SA17F2208-10S	2**2"-8ply	43*39*32cm	400pouches



03/40S 30X20MESH, 1PC/POUCH, 12POUCHES/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SA20-1840-12B-DL	18*40cm	56*31*24cm	576
SA20-3640-12B-DL	36*40cm	56*35*24cm	576
SA20-1840-12B-SL	18*40cm	39*28*35cm	576
SA20-3640-12B-SL	36*40cm	39*28*35cm	576
SA20-1840-12BK	18*40cm	50*40*36cm	2700
SA20-3640-12BK	36*40cm	50*40*36cm	1800
SA20-1840-100B	18*40cm	68*46*55cm	5000
SA20-3640-100B	36*40cm	55*35*55cm	2000

09/40S 19X15MESH, FOLDED EDGE, 10PCS/POUCH, 10POUCHES/BOX, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13F4816-10S	4**8"-16ply	52*28*40cm	80pouches
SA13F4416-10S	4**4"-16ply	55*30*40cm	160pouches
SA13F3316-10S	3**3"-16ply	53*28*40cm	200pouches
SA13F2216-10S	2**2"-16ply	43*39*40cm	400pouches
SA13F4812-10S	4**8"-12ply	52*28*32cm	80pouches
SA13F4412-10S	4**4"-12ply	55*30*32cm	160pouches
SA13F3312-10S	3**3"-12ply	53*28*32cm	200pouches
SA13F2212-10S	2**2"-12ply	43*39*32cm	400pouches
SA13F4808-10S	4**8"-8ply	52*26*28cm	80pouches
SA13F4408-10S	4**4"-8ply	51*28*32cm	160pouches
SA13F3308-10S	3**3"-8ply	55*26*28cm	200pouches
SA13F2208-10S	2**2"-8ply	43*26*39cm	400pouches

07/40S 30X20MESH, FOLDED EDGE, 10PCS/POUCH, 10POUCHES/BOX, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SA20F4816-10S	4**8"-16ply	52*26*50cm	80pouches
SA20F4416-10S	4**4"-16ply	55*30*50cm	160pouches
SA20F3316-10S	3**3"-16ply	53*28*50cm	200pouches
SA20F2216-10S	2**2"-16ply	43*39*50cm	400pouches
SA20F4812-10S	4**8"-12ply	52*28*46cm	80pouches
SA20F4412-10S	4**4"-12ply	55*30*46cm	160pouches
SA20F3312-10S	3**3"-12ply	53*28*46cm	200pouches
SA20F2212-10S	2**2"-12ply	43*39*46cm	400pouches
SA20F4808-10S	4**8"-8ply	52*28*36cm	80pouches
SA20F4408-10S	4**4"-8ply	55*30*36cm	160pouches
SA20F3308-10S	3**3"-8ply	53*28*36cm	200pouches
SA20F2208-10S	2**2"-8ply	43*39*36cm	400pouches

16/40S 18X11MESH, FOLDED EDGE, 2PCS/POUCH, 50POUCHES/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SA11F4812-2S	4**8"-12ply	55*27*34cm	400pouches
SA11F4412-2S	4**4"-12ply	68*33*34cm	1000pouches
SA11F3312-2S	3**3"-12ply	55*27*34cm	1000pouches
SA11F2212-2S	2**2"-12ply	43*39*34cm	2000pouches
SA11F4808-2S	4**8"-8ply	55*27*28cm	400pouches
SA11F4408-2S	4**4"-8ply	68*28*33cm	1000pouches
SA11F3308-2S	3**3"-8ply	55*27*28cm	1000pouches
SA11F2208-2S	2**2"-8ply	43*39*28cm	2000pouches

05/40S 30X20MESH, 25PCS/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SA201010-25BK	10*10,25PCS,25p/box	62*25*45cm	1000
SA201010-25-12BK	10*10,25PCS,12p/box	65*27*42cm	960
SA201010-25-100B	10*10,25PCS,100p/box	68*28*37cm	1000
SA202020-25-100B	20*20,25PCS,100p/box	75*47*53cm	600
SA203030-25-100B	30*30,25PCS,100p/box	71*47*71cm	400
SA204040-25-100B	40*40,25PCS,100p/box	91*47*47cm	200
SA201010-25-100	10*10,25PCS,100p/ctn	28*15*22cm	100
SA201010-25-4B	10*10,25PCS,4p/box	56*31*53cm	480pouches

11/40S 24X20MESH, FOLDED EDGE, 5PCS/POUCH, 20POUCHES/BOX, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17F4816-5S	4**8"-16ply	52*28*46cm	160pouches
SA17F4416-5S	4**4"-16ply	55*30*46cm	320pouches
SA17F3316-5S	3**3"-16ply	53*28*46cm	400pouches
SA17F2216-5S	2**2"-16ply	43*39*46cm	800pouches
SA17F4812-5S	4**8"-12ply	52*28*42cm	160pouches
SA17F4412-5S	4**4"-12ply	55*30*42cm	320pouches
SA17F3312-5S	3**3"-12ply	53*28*42cm	400pouches
SA17F2212-5S	2**2"-12ply	43*39*42cm	800pouches
SA17F4808-5S	4**8"-8ply	52*28*32cm	160pouches
SA17F4408-5S	4**4"-8ply	55*30*32cm	320pouches
SA17F3308-5S	3**3"-8ply	53*28*32cm	400pouches
SA17F2208-5S	2**2"-8ply	43*39*32cm	800pouches

10/40S 30X20MESH, FOLDED EDGE, 5PCS/POUCH, 20POUCHES/BOX, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SA20F4816-5S	4**8"-16ply	52*26*50cm	160pouches
SA20F4416-5S	4**4"-16ply	55*30*50cm	320pouches
SA20F3316-5S	3**3"-16ply	53*28*50cm	400pouches
SA20F2216-5S	2**2"-16ply	43*39*50cm	800pouches
SA20F4812-5S	4**8"-12ply	52*28*46cm	160pouches
SA20F4412-5S	4**4"-12ply	55*30*46cm	320pouches
SA20F3312-5S	3**3"-12ply	53*28*46cm	400pouches
SA20F2212-5S	2**2"-12ply	43*39*46cm	800pouches
SA20F4808-5S	4**8"-8ply	52*28*36cm	160pouches
SA20F4408-5S	4**4"-8ply	55*30*36cm	320pouches
SA20F3308-5S	3**3"-8ply	53*28*36cm	400pouches
SA20F2208-5S	2**2"-8ply	43*39*36cm	800pouches

17/40S 30X20MESH, FOLDED EDGE, 1PC/POUCH, 100POUCHES/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SA20F4812-2S	4**8"-12ply	55*27*28cm	400pouches
SA20F4412-2S	4**4"-12ply	68*33*28cm	1000pouches
SA20F3312-2S	3**3"-12ply	55*27*28cm	1000pouches
SA20F2212-2S	2**2"-12ply	43*39*28cm	2000pouches
SA20F4808-2S	4**8"-8ply	55*27*46cm	800pouches
SA20F4408-2S	4**4"-8ply	68*33*46cm	2000pouches
SA20F3308-2S	3**3"-8ply	55*27*46cm	2000pouches
SA20F2208-2S	2**2"-8ply	43*39*46cm	4000pouches

**12/40S 19X15MESH,FOLDED EDGE,5PCS/POUCH,
20POUCHES/BOX,BLISTER POUCH**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13F4816-5S	4**8"-16ply	52*28*40cm	160pouches
SA13F4416-5S	4**4"-16ply	55*30*40cm	320pouches
SA13F3316-5S	3**3"-16ply	53*28*40cm	400pouches
SA13F2216-5S	2**2"-16ply	43*39*40cm	800pouches
SA13F4812-5S	4**8"-12ply	52*28*32cm	160pouches
SA13F4412-5S	4**4"-12ply	55*30*32cm	320pouches
SA13F3312-5S	3**3"-12ply	53*28*32cm	400pouches
SA13F2212-5S	2**2"-12ply	43*39*32cm	800pouches
SA13F4808-5S	4**8"-8ply	52*26*28cm	160pouches
SA13F4408-5S	4**4"-8ply	51*30*28cm	320pouches
SA13F3308-5S	3**3"-8ply	55*26*28cm	400pouches
SA13F2208-5S	2**2"-8ply	43*26*39cm	800pouches

**14/40S 24X20MESH,FOLDED EDGE,2PCS/POUCH,
50POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17F4812-2S	4**8"-12ply	55*27*46cm	400pouches
SA17F4412-2S	4**4"-12ply	68*32*46cm	1000pouches
SA17F3312-2S	3**3"-12ply	55*26*46cm	1000pouches
SA17F2212-2S	2**2"-12ply	43*29*46cm	2000pouches
SA17F4808-2S	4**8"-8ply	55*27*38cm	400pouches

**13/40S 30X20MESH,FOLDED EDGE,2PCS/POUCH,
50POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA20F4812-2S	4**8"-12ply	55*27*50cm	400pouches
SA20F4412-2S	4**4"-12ply	68*32*50cm	1000pouches
SA20F3312-2S	3**3"-12ply	55*26*50cm	1000pouches
SA20F2212-2S	2**2"-12ply	43*39*50cm	2000pouches
SA20F4808-2S	4**8"-8ply	55*27*42cm	400pouches
SA20F4408-2S	4**4"-8ply	68*32*42cm	1000pouches
SA20F3308-2S	3**3"-8ply	55*26*42cm	1000pouches
SA20F2208-2S	2**2"-8ply	43*39*42cm	2000pouches

**19/40S 19X15MESH,FOLDED EDGE,1PC/POUCH,
100POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13F4812-2S	4**8"-12ply	55*27*40cm	800pouches
SA13F4412-2S	4**4"-12ply	68*33*40cm	2000pouches
SA13F3312-2S	3**3"-12ply	55*27*40cm	2000pouches
SA13F2212-2S	2**2"-12ply	43*39*40cm	4000pouches
SA13F4808-2S	4**8"-8ply	55*27*34cm	800pouches
SA13F4408-2S	4**4"-8ply	68*33*34cm	2000pouches
SA13F3308-2S	3**3"-8ply	55*27*34cm	2000pouches
SA13F2208-2S	2**2"-8ply	43*39*34cm	4000pouches

**27/40S 24X20MESH,NON-FOLDED EDGE,2PCS/POUCH,
50POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17NF4808-2S	4**8"-8ply	55*27*30cm	400pouches
SA17NF4408-2S	4**4"-8ply	68*30*33cm	1000pouches
SA17NF3308-2S	3**3"-8ply	55*30*53cm	2000pouches
SA17NF2208-2S	2**2"-8ply	43*39*30cm	2000pouches

**15/40S 19X15MESH,FOLDED EDGE,2PCS/POUCH,
50POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13F4812-2S	4**8"-12ply	55*27*38cm	400pouches
SA13F4412-2S	4**4"-12ply	68*33*38cm	1000pouches
SA13F3312-2S	3**3"-12ply	55*27*38cm	1000pouches
SA13F2212-2S	2**2"-12ply	40*39*38cm	2000pouches
SA13F4808-2S	4**8"-8ply	55*27*30cm	400pouches
SA13F4408-2S	4**4"-8ply	68*33*30cm	1000pouches
SA13F3308-2S	3**3"-8ply	55*27*30cm	1000pouches
SA13F2208-2S	2**2"-8ply	43*39*30cm	2000pouches

**24/40S 24X20MESH,NON-FOLDED EDGE,5PCS/POUCH,
BLISTER POUCH,20POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17NF4412-5S	4**4"-12ply	55*38*45cm	480pouches

**22/40S 24X20MESH,NON-FOLDED EDGE,10PCS/POUCH,
BLISTER POUCH**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17NF4808-10S	4**8"-8ply	41*26*44cm	100pouches
SA17NF4408-10S	4**4"-8ply	30*28*32cm	100pouches
SA17NF3308-10S	3**3"-8ply	35*26*30cm	100pouches
SA17NF2208-10S	2**2"-8ply	50*35*34cm	100pouches
SA17NF131308-10S	13**13"-24ply	44*30*36cm	50pouches

**18/40S 24X20MESH,FOLDED EDGE,1PC/POUCH,
100POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17F4812-2S	4**8"-12ply	55*27*50cm	800pouches
SA17F4412-2S	4**4"-12ply	68*33*50cm	2000pouches
SA17F3312-2S	3**3"-12ply	55*27*50cm	2000pouches
SA17F2212-2S	2**2"-12ply	43*39*50cm	4000pouches
SA17F4808-2S	4**8"-8ply	55*27*40cm	800pouches
SA17F4408-2S	4**4"-8ply	68*33*40cm	2000pouches
SA17F3308-2S	3**3"-8ply	55*27*40cm	2000pouches
SA17F2208-2S	2**2"-8ply	43*39*40cm	4000pouches

**20/40S 18X11MESH,FOLDED EDGE,1PC/POUCH,
100POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA11F4812-2S	4**8"-12ply	55*27*34cm	800pouches
SA11F4412-2S	4**4"-12ply	68*33*34cm	2000pouches
SA11F3312-2S	3**3"-12ply	55*27*34cm	2000pouches
SA11F2212-2S	2**2"-12ply	43*39*34cm	4000pouches
SA11F4808-2S	4**8"-8ply	55*27*34cm	800pouches
SA11F4408-2S	4**4"-8ply	68*33*34cm	2000pouches
SA11F3308-2S	3**3"-8ply	55*27*34cm	2000pouches
SA11F2208-2S	2**2"-8ply	43*39*34cm	4000pouches

**23/40S 24X20MESH,NON-FOLDED EDGE,5PCS/POUCH,
BLISTER POUCH**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17NF4808-5S	4**8"-8ply	41*26*23cm	100pouches
SA17NF4408-5S	4**4"-8ply	28*17*26cm	100pouches
SA17NF3308-5S	3**3"-8ply	26*17*20cm	100pouches
SA17NF2208-5S	2**2"-8ply	35*11*20cm	100pouches

**14/40S 24X20MESH,FOLDED EDGE,2PCS/POUCH,
50POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17F4408-2S	4**4"-8ply	68*30*38cm	1000pouches
SA17F3308-2S	3**3"-8ply	55*26*38cm	1000pouches
SA17F2208-2S	2**2"-8ply	43*39*38cm	2000pouches

**28/40S 24X20MESH,NON-FOLDED EDGE,1PC/POUCH,
100POUCHES/BOX**


Code no	Model	Carton size	Q'ty(pks/ctn)
SA17NF4808-1S	4**8"-8ply	55*27*34cm	400pouches
SA17NF4408-1S	4**4"-8ply	68*33*34cm	1000pouches
SA17NF3308-1S	3**3"-8ply	55*34*53cm	2000pouches
SA17NF2208-1S	2**2"-8ply	43*39*34cm	2000pouches

**29/40S 19X15MESH,NON-FOLDED EDGE,2PCS/POUCH,
50POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13NF4808-2S	4**8"-8ply	55*26*27cm	400pouches
SA13NF4408-2S	4**4"-8ply	68*26*33cm	1000pouches
SA13NF3308-2S	3**3"-8ply	55*26*53cm	2000pouches
SA13NF2208-2S	2**2"-8ply	43*26*39cm	2000pouches

Tampon Gauze

- ◆ Sterile or non sterile
- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's,40's
- ◆ Mesh of 22,20,17,15,13,12,11 threads ect
- ◆ 5cmx5m,7.5cmx5m,5cmx10m,7.5cmx10m,3.5cmx7m,7cmx7m ect
- ◆ Gamma,EO,Steam



**25/40S 19X15MESH,NON-FOLDED EDGE,10PCS/POUCH,
BLISTER POUCH**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13NF4408-10S	4**4"-8ply	58*28*30cm	200pouches
SA13NF3308-10S	3**3"-8ply	50*28*47cm	400pouches
SA13NF2208-10S	2**2"-8ply	50*24*35cm	400pouches

**26/40S 19X15MESH,NON-FOLDED EDGE,5PCS/POUCH,
BLISTER POUCH**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA17NF4808-5S	4**8"-8ply	58*58*30cm	400pouches
SA17NF4408-5S	4**4"-8ply	50*28*47cm	400pouches
SA17NF3308-5S	3**3"-8ply	50*24*35cm	800pouches

**30/40S 19X15MESH,NON-FOLDED EDGE,1PC/POUCH,
100POUCHES/BOX**

Code no	Model	Carton size	Q'ty(pks/ctn)
SA13NF4808-1S	4**8"-8ply	55*27*30cm	800pouches
SA13NF4408-1S	4**4"-8ply	68*30*33cm	2000pouches
SA13NF3308-1S	3**3"-8ply	55*30*53cm	4000pouches
SA13NF2208-1S	2**2"-8ply	43*30*39cm	4000pouches

01/40S,24X20MESH,ZIG-ZAG,1PC/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SL1710005M	10cm*5m-4ply	59*39*29cm	160
SL1707005M	7cm*5m-4ply	59*39*29cm	180
SL1705005M	5cm*5m-4ply	59*39*29cm	180
SL1705010M	5cm*10m-4ply	59*39*29cm	140
SL1707010M	7cm*10m-4ply	59*39*29cm	120

**02/40S,24X20MESH,WITH INDIFORM,ZIG-ZAG,
1PC/POUCH**

Code no	Model	Carton size	Q'ty(pks/ctn)
SLI1710005M	10cm*5m-4ply	59*39*47cm	140
SLI1707005M	7cm*5m-4ply	59*39*47cm	160
SLI1705005M	5cm*5m-4ply	59*39*47cm	160
SLI1702505M	2.5cm*5m-4ply	59*39*47cm	160
SLI1701005M	1cm*5m-4ply	59*39*47cm	200

03/40S,28X26MESH,1PC/ROLL,1PC/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SI2214007	14cm*7m	52*50*52cm	400pouches
SL2207007	7cm*7m	60*48*52cm	600pouches
SL2203507	3.5cm*7m	65*62*43cm	1000pouches

03/40S,24X20MESH,10PCS/BAG

Code no	Model	Carton size	Q'ty(pks/ctn)
SOTGA	MouseA	43*34*38cm	400pcs/ctn
SOTGB	MouseB	40*24*30cm	400pcs/ctn

Non Sterile Gauze Swab

- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 22,20,17,15,13,12,11 threads etc
- ◆ With or without x-ray detectable threads
- ◆ Folded edge or non-folded edge
- ◆ 8ply, 12ply, 16ply, 24ply, 32ply etc
- ◆ 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm, 10x40cm, 20x20cm, 30x30cm etc
- ◆ 100pcs/pack, 200pcs/pack(Non-sterile)



08/40S 19X10MESH, UNFOLDED EDGE, 200PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
NSGNF	2**2*-12ply	52*27*42cm	100
NSGNF	3**3*-12ply	52*32*42cm	40
NSGNF	4**4*-12ply	52*42*42cm	40
NSGNF	4**8*-12ply	52*42*28cm	20
NSGNF	4**8*-12ply+X-RAY	52*42*42cm	20

07/40S 19X10MESH, FOLDED EDGE, 200PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
NSGF2212	2**2*-12ply	52*27*46cm	100
NSGF3312	3**3*-12ply	40*32*46cm	40
NSGF4412	4**4*-12ply	52*42*46cm	40
NSGF4812	4**8*-12ply	52*42*32cm	20
NSGF4812X	4**8*-12ply+X-RAY	52*42*46cm	20

01/40S 30X20MESH, FOLDED EDGE, 100PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
A20F4816	4**8*-16ply	52*22*46cm	10
A20F4416	4**4*-16ply	52*22*46cm	20
A20F3316	3**3*-16ply	46*32*40cm	40
A20F2216	2**2*-16ply	52*22*46cm	80
A20F4812	4**8*-12ply	52*22*42cm	10
A20F4412	4**4*-12ply	52*22*42cm	20
A20F3312	3**3*-12ply	40*32*42cm	40
A20F2212	2**2*-12ply	52*22*42cm	80
A20F4808	4**8*-8ply	52*32*42cm	20
A20F4408	4**4*-8ply	52*32*52cm	50
A20F3308	3**3*-8ply	40*32*40cm	50
A20F2208	2**2*-8ply	52*27*32cm	100

02/40S 24X20MESH, FOLDED EDGE, 100PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
A17F4816	4**8*-16ply	52*22*42cm	10
A17F4416	4**4*-16ply	52*22*52cm	25
A17F3316	3**3*-16ply	40*22*40cm	25
A17F2216	2**2*-16ply	52*22*27cm	50
A17F4812	4**8*-12ply	52*38*42cm	20
A17F4412	4**4*-12ply	52*38*52cm	50
A17F3312	3**3*-12ply	40*38*40cm	50
A17F2212	2**2*-12ply	52*27*38cm	100
A17F4808	4**8*-8ply	52*28*42cm	20
A17F4408	4**4*-8ply	52*28*52cm	50
A17F3308	3**3*-8ply	40*28*40cm	50
A17F2208	2**2*-8ply	52*27*28cm	100

03/40S 19X15MESH, FOLDED EDGE, 100PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
A13F4816	4**8*-16ply	52*36*42cm	20
A13F4416	4**4*-16ply	52*36*52cm	50
A13F3316	3**3*-16ply	40*36*40cm	50
A13F2216	2**2*-16ply	52*27*36cm	100
A13F4812	4**8*-12ply	52*28*42cm	20
A13F4412	4**4*-12ply	52*28*52cm	50
A13F3312	3**3*-12ply	40*28*40cm	50
A13F2212	2**2*-12ply	52*28*27cm	100
A13F4808	4**8*-8ply	52*22*42cm	20
A13F4408	4**4*-8ply	52*22*52cm	50
A13F3308	3**3*-8ply	40*22*40cm	50
A13F2208	2**2*-8ply	52*22*27cm	100

04/40S 24X20MESH, NON-FOLDED EDGE, 100PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
A17NF4816	4**8*-16ply	52*38*42cm	20
A17NF4416	4**4*-16ply	52*38*52cm	50
A17NF3316	3**3*-16ply	40*38*40cm	50
A17NF2216	2**2*-16ply	52*27*38cm	100
A17NF4812	4**8*-12ply	52*32*42cm	20
A17NF4412	4**4*-12ply	52*32*52cm	50
A17NF3312	3**3*-12ply	40*32*40cm	50
A17NF2212	2**2*-12ply	52*27*32cm	100
A17NF4808	4**8*-8ply	52*22*42cm	20
A17NF4408	4**4*-8ply	52*22*52cm	50
A17NF3308	3**3*-8ply	40*22*40cm	50
A17NF2208	2**2*-8ply	52*22*27cm	100

05/40S 19X15MESH, NON-FOLDED EDGE, 100PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
A13NF4816	4**8*-16ply	52*32*42cm	20
A13NF4416	4**4*-16ply	52*32*52cm	50
A13NF3316	3**3*-16ply	40*32*40cm	50
A13NF2216	2**2*-16ply	52*27*32cm	100
A13NF4812	4**8*-12ply	52*26*42cm	20
A13NF4412	4**4*-12ply	52*26*52cm	50
A13NF3312	3**3*-12ply	40*26*40cm	50
A13NF2212	2**2*-12ply	52*26*27cm	100
A13NF4808	4**8*-8ply	52*42*42cm	50
A13NF4408	4**4*-8ply	52*42*42cm	100
A13NF3308	3**3*-8ply	42*32*40cm	100
A13NF2208	2**2*-8ply	42*22*27cm	100

06/40S 18X11MESH, NON-FOLDED EDGE, 100PCS/PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
A11NF4816	4**8*-16ply	52*28*42cm	20
A11NF4416	4**4*-16ply	52*28*52cm	50
A11NF3316	3**3*-16ply	40*28*40cm	50
A11NF2216	2**2*-16ply	52*27*28cm	100
A11NF4812	4**8*-12ply	52*22*42cm	20
A11NF4412	4**4*-12ply	52*22*52cm	50
A11NF3312	3**3*-12ply	40*22*40cm	50
A11NF2212	2**2*-12ply	52*22*27cm	100
A11NF4808	4**8*-8ply	52*37*42cm	50
A11NF4408	4**4*-8ply	52*37*42cm	100
A11NF3308	3**3*-8ply	40*32*37cm	100
A11NF2208	2**2*-8ply	37*22*27cm	100

Sterile Lap Sponge

- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 22,20,17,15,13,12,11 threads etc
- ◆ With or without x-ray detectable/x-ray film, with or without blue loop
- ◆ Pre-washed or non-washed
- ◆ 4ply, 6ply, 8ply, 12ply, 16ply, 24ply etc
- ◆ 20x20cm,22.5x22.5cm,30x30cm,40x40cm,45x45cm,50x50cm etc
- ◆ 5's, 2's, 1's in sterile pouch or double package
- ◆ Gamma,EO,Steam



01/40S 24X20MESH,WITH LOOP AND X-RAY DETECTABLE, NON-WASHED, 5PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC17454512-5S	45cm*45cm-12ply	50*32*45cm	30pouches
SC17404012-5S	40cm*40cm-12ply	57*27*40cm	20pouches
SC17303012-5S	30cm*30cm-12ply	50*32*40cm	60pouches
SC17454508-5S	45cm*45cm-8ply	50*32*30cm	30pouches
SC17404008-5S	40cm*40cm-8ply	57*27*40cm	30pouches
SC17303008-5S	30cm*30cm-8ply	50*32*40cm	90pouches
SC17454504-5S	45cm*45cm-4ply	50*32*45cm	90pouches
SC17404004-5S	40cm*40cm-4ply	57*27*40cm	60pouches
SC17303004-5S	30cm*30cm-4ply	50*32*40cm	180pouches

02/40S 24X20MESH,WITH LOOP AND X-RAY DETECTABLE, PRE-WASHED, 5PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC17454512PW-5S	45cm*45cm-12ply	57*30*32cm	30pouches
SC17404012PW-5S	40cm*40cm-12ply	57*30*28cm	30pouches
SC17303012PW-5S	30cm*30cm-12ply	52*29*32cm	50pouches
SC17454508PW-5S	45cm*45cm-8ply	57*30*32cm	40pouches
SC17404008PW-5S	40cm*40cm-8ply	57*30*28cm	40pouches
SC17303008PW-5S	30cm*30cm-8ply	52*29*32cm	60pouches
SC17454504PW-5S	45cm*45cm-4ply	57*30*32cm	50pouches
SC17404004PW-5S	40cm*40cm-4ply	57*30*28cm	50pouches
SC17303004PW-5S	30cm*30cm-4ply	52*29*32cm	100pouches

04/40S 30X20MESH,WITH LOOP AND X-RAY DETECTABLE FILM, NON-WASHED, 3PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC20101004-3S	10cm*10cm-4ply	51*26*27cm	400pouches
SC20103004-3S	10cm*30cm-4ply	62*28*30cm	400pouches
SC20106004-3S	10cm*60cm-4ply	62*28*48cm	400pouches
SC20202004-3S	20cm*20cm-4ply	62*28*41cm	400pouches
SC20203004-3S	20cm*30cm-4ply	62*28*26cm	200pouches
SC20303004-3S	30cm*30cm-4ply	62*28*45cm	200pouches
SC20304504-3S	30cm*45cm-4ply	66*34*42cm	200pouches
SC20404004-3S	40cm*40cm-4ply	60*46*38cm	200pouches
SC20457004-3S	45cm*70cm-4ply	66*34*38cm	100pouches
SC20505004-3S	50cm*50cm-4ply	62*38*30cm	100pouches

07/40S 19X11MESH,WITH LOOP AND X-RAY DETECTABLE, NON-WASHED FILM, 5PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC12454512-5S	45cm*45cm-12ply	50*32*45cm	60pouches
SC12303012-5S	30cm*30cm-12ply	50*32*46cm	120pouches

05/40S 19X15MESH,WITH LOOP AND X-RAY DETECTABLE, NON-WASHED, 5PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC13454512-5S	45cm*45cm-12ply	50*32*56cm	60pouches
SC13404012-5S	40cm*40cm-12ply	57*27*52cm	50pouches
SC13303012-5S	30cm*30cm-12ply	50*32*54cm	120pouches
SC13454508-5S	45cm*45cm-8ply	50*32*56cm	90pouches
SC13404008-5S	40cm*40cm-8ply	57*27*52cm	75pouches
SC13303008-5S	30cm*30cm-8ply	50*32*54cm	180pouches
SC13454504-5S	45cm*45cm-4ply	50*32*56cm	120pouches
SC13404004-5S	40cm*40cm-4ply	57*27*52cm	100pouches
SC13303004-5S	30cm*30cm-4ply	50*32*54cm	240pouches

06/40S 19X15MESH,WITH LOOP AND X-RAY DETECTABLE, PRE-WASHED FILM, 5PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC13454512PW-5S	45cm*45cm-12ply	48*30*32cm	30pouches
SC13404012PW-5S	40cm*40cm-12ply	48*30*28cm	30pouches
SC13303012PW-5S	30cm*30cm-12ply	44*29*32cm	50pouches
SC13454508PW-5S	45cm*45cm-8ply	48*30*32cm	40pouches
SC13404008PW-5S	40cm*40cm-8ply	48*30*28cm	40pouches
SC13303008PW-5S	30cm*30cm-8ply	44*29*32cm	60pouches
SC13454504PW-5S	45cm*45cm-4ply	48*30*32cm	50pouches
SC13404004PW-5S	40cm*40cm-4ply	48*30*28cm	50pouches
SC13303004PW-5S	30cm*30cm-4ply	44*29*32cm	100pouches

03/40S 30X20MESH,WITH LOOP AND X-RAY DETECTABLE FILM, PRE-WASHED, 5PCS/BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SC20454504PW-5S	45cm*45cm-4ply	51*30*36cm	50pouches



Non Sterile Lap Sponge

- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 22,20,17,15,13,12,11 threads etc
- ◆ With or without x-ray detectable/x-ray chip, with or without blue loop
- ◆ Pre-washed or non-washed
- ◆ 4ply, 6ply, 8ply, 12ply, 16ply, 24ply etc
- ◆ 20x20cm, 22.5x22.5cm, 30x30cm, 40x40cm, 45x45cm, 50x50cm etc
- ◆ 10s, 25s, 50s, 100s per pack (Non sterile)



01/40S 30X20MESH, WITH LOOP AND X-RAY DETECTABLE LINE, 50PCS/PE-BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
C20457004	45cm*70cm-4ply	50*32*38cm	300
C20505004	50cm*50cm-4ply	52*34*52cm	400
C20454504	45cm*45cm-4ply	46*46*37cm	400
C20404004	40cm*40cm-4ply	62*42*37cm	600
C20304504	30cm*45cm-4ply	47*47*37cm	600
C20304004	30cm*40cm-4ply	47*42*37cm	600
C20303004	30cm*30cm-4ply	47*32*37cm	600
C20252504	25cm*25cm-4ply	51*38*32cm	1200
C20203004	20cm*30cm-4ply	52*32*37cm	1000
C20202004	20cm*20cm-4ply	52*42*37cm	2000
C20104504	10cm*45cm-4ply	47*32*42cm	1800
C20106004	10cm*60cm-4ply	62*32*42cm	1800

02/40S 19X15MESH, WITH LOOP AND X-RAY DETECTABLE LINE, PRE-WASHED, 50PCS/PE-BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
C13454512PW	45cm*45cm-12ply	57*30*42cm	200
C13404012PW	40cm*40cm-12ply	48*30*38cm	200
C13303012PW	30cm*30cm-12ply	52*36*40cm	500
C13303012PW-5P	30cm*30cm-12ply	57*25*46cm	100 pk
C13454508PW	45cm*45cm-8ply	57*42*42cm	400
C13454508PW-5P	45cm*45cm-8ply	60*28*50cm	80 pk
C13404008PW	40cm*40cm-8ply	48*42*36cm	400
C13303008PW	30cm*30cm-8ply	57*36*45cm	600
C13454504PW	45cm*45cm-4ply	57*42*42cm	800
C13454504PW-5P	45cm*45cm-4ply	54*39*52cm	200 pk
C13404004PW	40cm*40cm-4ply	48*42*38cm	800
C13303004PW	30cm*30cm-4ply	57*40*45cm	1200
C13303004PW-5P	30cm*30cm-4ply	57*38*40cm	200 pk

03/40S 24X20MESH, WITH LOOP AND X-RAY DETECTABLE FILM, PRE-WASHED, 50PCS/PE-BAG OR 25PCS/PE-BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
C17454512PW	45cm*45cm-12ply	59*37*44cm	200
C17454512PW-5P	45cm*45cm-12ply	59*37*45cm	200
C17404012PW	40cm*40cm-12ply	38*37*40cm	200
C17303012PW	30cm*30cm-12ply	52*40*41cm	500
C17303012PW-5P	30cm*30cm-12ply	57*28*37cm	100
C17454508PW	45cm*45cm-8ply	59*44*46cm	400
C17404008PW	40cm*40cm-8ply	48*44*42cm	400
C17303008PW	30cm*30cm-8ply	57*40*45cm	600
C17454504PW	45cm*45cm-4ply	59*44*46cm	800
C17454504PW-5P	45cm*45cm-4ply	58*25*28cm	40 pk
C17404004PW	40cm*40cm-4ply	48*44*42cm	800
C17303004PW	30cm*30cm-4ply	57*45*45cm	1200
C17104504PW	10cm*45cm-4ply	40*22*42cm	500
C17104504PW-5P	10cm*45cm-4ply	40*18*39cm	100 pk
C17209004PW	20cm*90cm-4ply	42*42*62cm	500

04/40S 24X20MESH, WITH LOOP AND X-RAY DETECTABLE, NON-WASHED, 50PCS/PE-BAG OR 25PCS/PE-BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
C17292932	29cm*29cm-32ply	60*31*47cm	200
C1732532524	32.5cm*32.5cm-24ply	66*34*36cm	200
C17292924	29cm*29cm-24ply	60*34*37cm	250
C17232324	23cm*23cm-24ply	60*38*49cm	500
C17202024	20cm*20cm-24ply	51*40*42cm	500
C17292916	29cm*29cm-16ply	60*31*47cm	400
C17454512	45cm*45cm-12ply	49*32*47cm	200
C17404012	40cm*40cm-12ply	49*42*42cm	300
C17303012	30cm*30cm-12ply	62*36*32cm	400
C17303012-5P	30cm*30cm-12ply	60*32*33cm	80
C17454508	45cm*45cm-8ply	62*38*47cm	400
C17404008	40cm*40cm-8ply	55*33*42cm	400
C17303008	30cm*30cm-8ply	42*32*46cm	800
C1722522508	22.5cm*22.5cm-8ply	52*24*46cm	800
C17404006	40cm*40cm-6ply	48*42*42cm	400
C17454504	45cm*45cm-4ply	62*38*47cm	800
C17404004	40cm*40cm-4ply	56*42*46cm	800
C17303004	30cm*30cm-4ply	62*32*27cm	1000
C17104504	10cm*45cm-4ply	47*42*40cm	2000
C17154504	15cm*45cm-4ply	62*38*32cm	800
C17253504	25cm*35cm-4ply	54*39*52cm	1600
C17304504	30cm*45cm-4ply	62*32*48cm	800

05/40S 19X15MESH, WITH LOOP AND X-RAY DETECTABLE, NON-WASHED, 50PCS/PE-BAG OR 25PCS/PE-BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
C13292932	29cm*29cm-32ply	53*31*48cm	250
C13202032	20cm*20cm-32ply	52*22*32cm	250
C13292924	29cm*29cm-24ply	53*31*37cm	250
C13232324	23cm*23cm-24ply	57*27*48cm	500
C13202024	20cm*20cm-24ply	52*26*42cm	500
C13454516	45cm*45cm-16ply	46*45*45cm	200
C13303016	30cm*30cm-16ply	60*32*47cm	400
C13292916	29cm*29cm-16ply	58*30*47cm	400
C13232316	23cm*23cm-16ply	57*25*36cm	500
C1322522516	22.5cm*22.5cm-16ply	57*35*46cm	1000
C13202016	20cm*20cm-16ply	52*34*42cm	1000
C13454512	45cm*45cm-12ply	62*47*40cm	400
C13404012	40cm*40cm-12ply	52*42*40cm	400
C13303012	30cm*30cm-12ply	62*32*32cm	400
C13303012-5p	30cm*30cm-12ply	60*32*35cm	80 pk
C1322522512	22.5cm*22.5cm-12ply	57*38*47cm	800
C13454508	45cm*45cm-8ply	62*27*46cm	400
C13454508-5P	45cm*45cm-8ply	59*26*50cm	80 pk
C13404008	40cm*40cm-8ply	52*30*42cm	400
C13303008	30cm*30cm-8ply	62*32*36cm	800
C1322522508	22.5cm*22.5cm-8ply	57*38*42cm	1000
C13454504	45cm*45cm-4ply	62*46*34cm	800
C13454504-5P	45cm*45cm-4ply	61*37*50cm	200 pk
C13404004	40cm*40cm-4ply	52*30*42cm	800
C13303004	30cm*30cm-4ply	62*32*36cm	1600
C13303004-5P	30cm*30cm-4ply	55*32*32cm	200 pk

Non Sterile Gauze Bandage

- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 22,20,17,15,13,12,11 threads etc
- ◆ Width: 5cm, 7.5cm, 14cm, 15cm, 20cm
- ◆ Length: 10m, 10yards, 7m, 5m, 5yards, 4m, 4yards, 3m, 3yards
- ◆ 10rolls/pack, 12rolls/pack(Non-sterile)
- ◆ 1roll packed into pouch/box(Sterile)
- ◆ Gamma,EO,Steam



01/21S 30X20MESH, 1PC/WHITE PAPER PACKAGE, 12ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D21201010M	10CM*10M	51*31*52cm	25
D21201510M	15CM*10M	60*32*50cm	20

02/40S 19X11MESH, 1PC/WHITE PAPER PACKAGE, 1ROLL/BOX, 12BOXES/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
D1205010YBS	2**10yards	39*36*32cm	600
D1275010YBS	3**10yards	39*36*44cm	600
D1210010YBS	4**10yards	39*36*57cm	600

03/40S 30X20MESH, 1PC/WHITE PAPER PACKAGE, 12ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D2005005M	2**5M	54*33*52cm	200
D2007005M	3**5M	54*33*40cm	100
D2010005M	4**5M	54*33*52cm	100
D2015005M	6**5M	77*28*33cm	50
D2005004Y	2**4yards	52*29*46cm	200
D2007504Y	3**4yards	46*29*40cm	100
D2010004Y	4**4yards	52*29*46cm	100
D2015004Y	6**4yards	68*24*32cm	50

04/40S 30X20MESH, 1PC/WHITE PAPER PACKAGE, 10ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D2007010M-10	7CM*10M	37*37*44cm	60
D2010010M-10	10CM*10M	37*37*62cm	60
D2015010M-10	15CM*10M	37*37*47cm	30
D2020010M-10	20CM*10M	37*37*42cm	20
D2012010M	120CM*10M	122*37*35cm	100
D2005005M	5CM*5M	42*39*62cm	288
D2007005M	7CM*5M	42*39*62cm	196
D2010005M	10CM*5M	42*39*62cm	144

04/40S 30X20MESH, 1PC/WHITE PAPER PACKAGE, 10ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D2015005M	15CM*5M	42*39*62cm	96
D2020005M	20CM*5M	42*39*62cm	72
D2012005M	120CM*5M	122*27*25cm	100

05/40S 24X20MESH, 1PC/WHITE PAPER PACKAGE, 12ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D1705010M	2**10M	52*36*43cm	100
D1707510M	3**10M	40*36*43cm	50
D1710010M	4**10M	52*36*43cm	50
D1715010M	6**10M	47*36*43cm	30
D1720010M	8**10M	42*36*43cm	20
D1705010Y	2**10yards	52*37*44cm	100
D1707510Y	3**10yards	40*37*44cm	50
D1710010Y	4**10yards	52*37*44cm	50
D1715010Y	6**10yards	47*37*44cm	30
D1720010Y	8**10yards	42*37*44cm	20
D1705006Y	2**6yards	52*27*32cm	100
D1707506Y	3**6yards	40*27*32cm	50
D1710006Y	4**6yards	52*27*32cm	50
D1715006Y	6**6yards	47*27*32cm	30
D1720006Y	8**6yards	42*27*32cm	20
D1705005M	2**5M	52*27*32cm	100
D1707505M	3**5M	40*27*32cm	50
D1710005M	4**5M	52*27*32cm	50
D1715005M	6**5M	47*27*32cm	30
D1720005M	8**5M	42*27*32cm	20
D1705005Y	2**5yards	52*25*30cm	100
D1707505Y	3**5yards	40*25*30cm	50
D1710005Y	4**5yards	52*25*30cm	50
D1715005Y	6**5yards	47*25*30cm	30
D1720005Y	8**5yards	42*25*30cm	20
D1708004M-10	8CM*4M	46*24*42cm	100
D1705010M-10	5CM*10M	52*36*36cm	100

06/40S 24X20MESH, 1PC/WHITE PAPER PACKAGE, 1ROLL/BLUE, 12BOXES/SHRINK

Code no	Model	Carton size	Q'ty(pks/ctn)
D1710005M-60	10CM*5M	65*25*58cm	24
D1714007M-60	14CM*7M	49*30*58cm	9
D1702510YBS	1**10yards	63*32*35cm	1200
D1705010YBS	2**10yards	50*32*38cm	720
D1707510YBS	3**10yards	35*33*33cm	720
D1710010YBS	4**10yards	38*32*35cm	360
D1705006YBS	2**6yards	50*32*35cm	720
D1707506YBS	3**6yards	50*32*38cm	720
D1710006YBS	6**6yards	38*33*33cm	360

07/40S 19X15MESH, 1PC/WHITE PAPER PACKAGE, 12ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D1305006Y	2**6yards	52*24*28cm	100
D1307506Y	3**6yards	40*24*28cm	50
D1310006Y	4**6yards	52*24*28cm	50
D1315006Y	6**6yards	47*24*28cm	30
D1305005M	2**5M	52*25*30cm	100
D1307505M	3**5M	40*25*30cm	50
D1310005M	4**5M	52*25*30cm	50
D1315005M	6**5M	47*25*30cm	30
D1305005Y	2**5yards	52*22*26cm	100
D1307505Y	3**5yards	40*22*26cm	50
D1310005Y	4**5yards	52*22*26cm	50
D1315005Y	6**5yards	47*22*26cm	30
D1305004M	2**4M	52*21*25cm	100
D1307504M	3**4M	40*21*25cm	50
D1310004M	4**4M	52*21*25cm	50
D1315004M	6**4M	47*21*25cm	30

07/40S 19X15MESH, 1PC/WHITE PAPER PACKAGE, 12ROLLS/BLUE PAPER PACKAGE

Code no	Model	Carton size	Q'ty(pks/ctn)
D1305004Y	2**4yards	52*20*24cm	100
D1307504Y	3**4yards	40*20*24cm	50
D1310004Y	4**4yards	52*20*24cm	50
D1315004Y	6**4yards	47*20*24cm	30
D1305003M	2**3M	36*23*27cm	100
D1307503M	3**3M	40*23*36cm	100
D1310003M	4**3M	52*23*26cm	100
D1315003M	6**3M	36*23*32cm	40
D1305003Y	2**3yards	36*23*27cm	100
D1307503Y	3**3yards	40*23*36cm	100
D1310003Y	4**3yards	52*23*26cm	100
D1315003Y	6**3yards	36*23*32cm	40
D1308004M	8CM*4M	40*25*42cm	100

08/40S 30X30MESH, 1PC/WHITE PAPER PACKAGE, 12ROLLS/PE-BAG

Code no	Model	Carton size	Q'ty(pks/ctn)
D2402505M-SP	1**5M	32*20*28cm	600
D2405005M-SP	2**5M	34*28*30cm	500
D2407505M-SP	3**5M	34*20*30cm	300
D2410005M-SP	4**5M	34*28*30cm	300

09/40S 43X27MESH, 1PC/CANDY PACKAGE, 12ROLLS/PE-BAG

Code no	Model	Carton size	Q'ty(pks/ctn)
D2805004M-C	2**4M	37*37*44cm	960
D2807504M-C	3**4M	42*37*45cm	600
D2810004M-C	4**4M	40*36*37cm	480
D2815004M-C	6**4M	37*37*40cm	360

Sterile Gauze Bandage

01/32S 28X26MESH, 1PC/PAPER BAG, 50ROLLS/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SD322414007M-1S	14cm*7m	63*40*40cm	400

02/40S 28X26MESH, 1PC/PAPER BAG, 50ROLLS/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SD2414007M-1S	14cm*7m	66.5*35*37.5cm	400

03/40S 24X20MESH, 1PC/PAPER BAG, 50ROLLS/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SD1714007M-1S	14cm*7m	35*20*32cm	100
SD1710005M-1S	10cm*5m	45*15*21cm	100

04/40S 19X15MESH, 1PC/PE-BAG

Code no	Model	Carton size	Q'ty(pks/ctn)
SD1390005M-8P-S	90cm*5m-8ply	52*28*42cm	200
SD1380005M-4P-XS	80cm*5m-4ply+X ray	55*29*37cm	200





Sterile Paraffin Gauze

- ◆ 100% cotton
- ◆ Cotton yarn of 21's, 32's
- ◆ Mesh of 22,20,17 etc
- ◆ 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm, 10x30cm, 10x40cm, 10cmx5m, 7m etc
- ◆ Package: in 1's, 10's, 12's packed into pouch.
10's, 12's, 36's/Tin
- ◆ Box: 10, 50 pouches/box
- ◆ Gamma sterilization

01/PARAFFIN GAUZE, 1PC/POUCH, 10 POUCHES/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SP44-10T	10*10cm	59*25*31cm	100tin
SP44-12T	10*10cm	59*25*31cm	100tin
SP44-36T	10*10cm	59*25*31cm	100tin
SP44-500T	10*500cm	59*25*31cm	100tin
SP44-700T	10*700cm	59*25*31cm	100tin
SP44-800T	10*800cm	59*25*31cm	100tin
SP22-10B	5*5cm	45*21*41cm	2000pouches
SP33-10B	7.5*7.5cm	60*33*33cm	2000pouches
SP44-10B	10*10cm	40*29*33cm	1000pouches
SP48-10B	10*20cm	40*29*33cm	1000pouches
SP412-10B	10*30cm	53*29*33cm	1000pouches
SP416-10B	10*40cm	53*29*33cm	1000pouches
SP102-1B	10cm*2m	53*27*32cm	150rolls
SP152-1B	15cm*2m	53*27*32cm	100rolls
SP202-1B	20cm*2m	53*27*32cm	60rolls

02/PARAFFIN GAUZE, WITH CHLORHEXIDINE ACETATE 0.5% OR NEOMYCIN SULPHATE 0.5% 1PC/POUCH, 10 POUCHES/BOX

Code no	Model	Carton size	Q'ty(pks/ctn)
SPCA44-10T	10*10cm	59*25*31cm	100tin
SPCA44-36T	10*10cm	59*25*31cm	100tin
SPCA44-500T	10*500cm	59*25*31cm	100tin
SPCA44-700T	10*700cm	59*25*31cm	100tin
SPCA22-10B	5*5cm	45*21*41cm	2000pouches
SPCA33-10B	7.5*7.5cm	60*33*33cm	2000pouches
SPCA44-10B	10*10cm	40*29*33cm	1000pouches
SPCA48-10B	10*20cm	40*29*33cm	1000pouches
SPCA412-10B	10*30cm	53*29*33cm	1000pouches



Gauze Roll

- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 22, 20, 17, 15, 13, 11 threads etc
- ◆ With or without x-ray
- ◆ 1ply, 2ply, 4ply, 8ply.
- ◆ Zigzag gauze roll, pillow gauze roll, rounded gauze roll
- ◆ 36"x100m, 36"x100yards, 36"x50m, 36"x5m, 36"x100m etc
- ◆ Packing: 1roll/blue kraft paper or polybag
- ◆ 10rolls, 12rolls, 20rolls/ctn



01/GAUZE ROLL

Code no	Model	Carton size	Q'ty(pks/ctn)
R2036100Y-4P	30*20mesh, 40s/40s	66*44*44cm	12rolls
R2036100M-4P	30*20mesh, 40s/40s	65*44*46cm	12rolls
R2036100Y-2P	30*20mesh, 40s/40s	58*44*47cm	12rolls
R2036100M-2P	30*20mesh, 40s/40s	58*44*49cm	12rolls
R173650M-4P	24*20mesh, 40s/40s	50*42*46cm	12rolls
R133650M-4P	19*15mesh, 40s/40s	68*36*46cm	20rolls
R123650M-4P	19*10mesh, 40s/40s	56*33*46cm	20rolls
R113650M-4P	19*8mesh, 40s/40s	54*32*46cm	20rolls
R83650M-4P	12*8mesh, 40s/40s	42*24*46cm	20rolls
R1736100Y-2P	24*20mesh, 40s/40s	57*42*47cm	12rolls
R1336100Y-2P	19*15mesh, 40s/40s	77*37*47cm	20rolls
R1236100Y-2P	19*10mesh, 40s/40s	67*32*47cm	20rolls
R1136100Y-2P	19*8mesh, 40s/40s	62*30*47cm	20rolls
R836100Y-2P	12*8mesh, 40s/40s	58*28*47cm	20rolls
R1736100M-2P	24*20mesh, 40s/40s	57*42*47cm	12rolls
R1336100M-2P	19*15mesh, 40s/40s	77*36*47cm	20rolls
R1236100M-2P	19*10mesh, 40s/40s	67*33*47cm	20rolls
R1136100M-2P	19*8mesh, 40s/40s	62*32*47cm	20rolls
R836100M-2P	12*8mesh, 40s/40s	58*24*47cm	20rolls
R1736100Y-4P	24*20mesh, 40s/40s	57*39*46cm	12rolls
R1336100Y-4P	19*15mesh, 40s/40s	70*29*47cm	20rolls
R1236100Y-4P	19*10mesh, 40s/40s	67*28*46cm	20rolls
R1136100Y-4P	19*8mesh, 40s/40s	62*26*46cm	20rolls
R836100Y-4P	12*8mesh, 40s/40s	58*25*46cm	20rolls
R1736100M-4P	24*20mesh, 40s/40s	57*42*46cm	12rolls
R1336100M-4P	19*15mesh, 40s/40s	77*36*46cm	20rolls
R1236100M-4P	19*10mesh, 40s/40s	67*33*46cm	20rolls
R1136100M-4P	19*8mesh, 40s/40s	62*32*46cm	20rolls
R13365M-4PLY	19x15mesh, 40s/40s	36"x5m-4ply	400rolls

01/GAUZE ROLL

Code no	Model	Carton size
R20361000	30*20mesh, 40s/40s	diameter:38cm
R17361000	24*20mesh, 40s/40s	diameter:36cm
R13361000	19*15mesh, 40s/40s	diameter:32cm
R12361000	19*10mesh, 40s/40s	diameter:30cm
R11361000	19*8mesh, 40s/40s	diameter:28cm
R20362000	30*20mesh, 40s/40s	diameter:53cm
R17362000	24*20mesh, 40s/40s	diameter:50cm
R13362000	19*15mesh, 40s/40s	diameter:45cm
R12362000	19*10mesh, 40s/40s	diameter:40cm
R11362000	19*8mesh, 40s/40s	diameter:36cm
R17363000	24*20mesh, 40s/40s	diameter:57cm
R17366000	24*20mesh, 40s/40s	diameter:112cm

02/PILLOW GAUZE ROLL

Code no	Model	Carton size	Q'ty(pks/ctn)
RRR1736100Y-10R	24*20mesh, 40s/40s	74*38*46cm	10rolls
RRR1536100Y-10R	20*16mesh, 40s/40s	74*33*46cm	10rolls
RRR1336100Y-10R	20*12mesh, 40s/40s	74*29*46cm	10rolls
RRR1336100Y-30R	20*12mesh, 40s/40s	90*46*48cm	30rolls
RRR1336100Y-40R	24*12mesh, 40s/40s	110*48*50cm	40rolls

03/ZIG-ZAG GAUZE ROLL

Code no	Model	Carton size	Q'ty(pks/ctn)
RZZ1765100M	24*20mesh, 40s/40s	70*38*44cm	20pcs
RZZ1790100M	24*20mesh, 40s/40s	62*35*42cm	20pcs
RZZ17120100M	24*20mesh, 40s/40s	42*35*42cm	10pcs
RZZ1365100M	19*15mesh, 40s/40s	70*38*35cm	20pcs

Sterile Non-Woven Sponge

- ◆ Made of spunlace non woven material, 70% viscose + 30% polyester
- ◆ Weight: 30, 35, 40, 50gsm/sq
- ◆ With or without x-ray detectable
- ◆ 4ply, 6ply, 8ply, 12ply
- ◆ 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm etc
- ◆ 1's, 2's, 5's, 10's packed into pouch (Sterile)
- ◆ Box: 100, 50, 25, 10, 4pouches/box
- ◆ Pouch: paper+paper, paper+film
- ◆ Gamma, EO, Steam



01/55G/M2, 1PC/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB55440401-50B	4**8"-4ply	43*30*40cm	18
SB55330401-50B	3**3"-4ply	46*37*40cm	36
SB55220401-50B	2**2"-4ply	40*29*35cm	36
SB55440401-25B	4**8"-4ply	40*29*45cm	36
SB55330401-25B	3**3"-4ply	40*34*49cm	72
SB55220401-25B	2**2"-4ply	40*36*30cm	72
SB55440401-10B	4**8"-4ply	57*24*45cm	72
SB55330401-10B	3**3"-4ply	35*31*37cm	72
SB55220401-10B	2**2"-4ply	36*24*29cm	72

02/40G/M2, 5PCS/POUCH, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB40480405-20B	4**8"-4ply	42*36*53cm	240pouches
SB40440405-20B	4**4"-4ply	55*36*44cm	480pouches
SB40330405-20B	3**3"-4ply	50*36*42cm	600pouches
SB40220405-20B	2**2"-4ply	43*36*50cm	1000pouches
SB40480805-20B	4**8"-8ply	42*39*53cm	240pouches
SB40440805-20B	4**4"-8ply	55*39*44cm	480pouches
SB40330805-20B	3**3"-8ply	50*39*42cm	600pouches
SB40220805-20B	2**2"-8ply	43*39*50cm	1000pouches

03/40G/M2, 2PCS/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB40480402-50B	4**8"-4ply	55*27*40cm	400pouches
SB40440402-50B	4**4"-4ply	68*33*40cm	1000pouches
SB40330402-50B	3**3"-4ply	55*27*40cm	1000pouches
SB40220402-50B	2**2"-4ply	50*35*40cm	2000pouches
SB40480402-25B	4**8"-4ply	55*27*40cm	400pouches
SB40440402-25B	4**4"-4ply	68*33*40cm	1000pouches
SB40330402-25B	3**3"-4ply	55*27*40cm	1000pouches
SB40220402-25B	2**2"-4ply	55*35*40cm	2000pouches
SB40480402-12B	4**8"-4ply	53*28*53cm	480pouches
SB40440402-12B	4**4"-4ply	53*28*33cm	960pouches
SB40330402-12B	3**3"-4ply	45*28*33cm	960pouches
SB40220402-12B	2**2"-4ply	53*35*41cm	1920pouches

04/40G/M2, 1PC/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB40480401-100B	4**8"-4ply	55*27*40cm	800pouches
SB40440401-100B	4**4"-4ply	68*33*40cm	2000pouches
SB40330401-100B	3**3"-4ply	55*27*40cm	2000pouches
SB40220401-100B	2**2"-4ply	50*35*40cm	4000pouches
SB40480401-50B	4**8"-4ply	55*27*40cm	800pouches
SB40440401-50B	4**4"-4ply	68*33*40cm	2000pouches
SB40330401-50B	3**3"-4ply	55*27*40cm	2000pouches
SB40220401-50B	2**2"-4ply	50*35*40cm	4000pouches
SB40480401-25B	4**8"-4ply	53*28*53cm	1000pouches
SB40440401-25B	4**4"-4ply	53*28*33cm	1000pouches
SB40330401-25B	3**3"-4ply	53*27*45cm	2000pouches
SB40220401-25B	2**2"-4ply	53*35*41cm	4000pouches
SB40480401-12B	4**8"-4ply	53*28*53cm	960pouches
SB40440401-12B	4**4"-4ply	53*28*33cm	960pouches
SB40330401-12B	3**3"-4ply	53*27*45cm	1920pouches
SB40220401-12B	2**2"-4ply	53*35*41cm	3840pouches
SB40480401-10B	4**8"-4ply	53*28*53cm	800pouches
SB40440401-10B	4**4"-4ply	53*28*33cm	800pouches
SB40330401-10B	3**3"-4ply	53*27*45cm	1600pouches
SB40220401-10B	2**2"-4ply	53*35*41cm	3600pouches

05/40G/M2, 10PCS/POUCH, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB40480410-10B	4**8"-4ply	42*36*53cm	120pouches
SB40440410-10B	4**4"-4ply	55*36*44cm	240pouches
SB40330410-10B	3**3"-4ply	50*36*42cm	300pouches
SB40220410-10B	2**2"-4ply	43*36*50cm	500pouches
SB40480610-10B	4**8"-6ply	68*24*28cm	50pouches
SB40440610-10B	4**4"-6ply	68*24*33cm	100pouches
SB40330610-10B	3**3"-6ply	53*24*51cm	200pouches
SB40220610-10B	2**2"-6ply	48*24*43cm	250pouches

06/30G/M2, 5PCS/POUCH, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB40480405-20B	4**8"-4ply	52*26*28cm	160pouches
SB40440405-20B	4**4"-4ply	55*30*30cm	320pouches
SB40330405-20B	3**3"-4ply	52*28*30cm	400pouches
SB40220405-20B	2**2"-4ply	52*35*30cm	800pouches

07/30G/M2, 1PC/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB30441201-50B	4**4"-12ply	71.5*41.5*43cm	1800pouches
SB30331201-50B	3**3"-12ply	56*38*51cm	2000pouches
SB30221201-50B	2**2"-12ply	43*38*51cm	2500pouches
SB30481201P-25B	4**8"-12ply	60.5*42*53cm	900pouches
SB30441201P-25B	4**4"-12ply	60.5*41.5*32.5cm	900pouches
SB30331201P-40B	3**3"-12ply	50*45.5*27cm	960pouches
SB30221201P-50B	2**2"-12ply	57*27*21.5cm	900pouches
SB30441201-25B	4**4"-12ply	55*29*43cm	900pouches
SB30331201-25B	3**3"-12ply	57*45*35cm	1800pouches
SB30221201-25B	2**2"-12ply	51*29*36cm	1800pouches
SB30441201-10B	4**4"-12ply	55*26*43cm	720pouches
SB30331201-10B	3**3"-12ply	45*26*35cm	720pouches
SB30221201-10B	2**2"-12ply	35*26*28cm	720pouches

08/30G/M2, 2PCS/POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB30480402-50B	4**8"-4ply	55*27*34cm	400pouches
SB30440402-50B	4**4"-4ply	68*33*40cm	1000pouches
SB30330402-50B	3**3"-4ply	55*27*34cm	1000pouches
SB30220402-50B	2**2"-4ply	50*35*34cm	2000pouches
SB30480402-25B	4**8"-4ply	55*27*34cm	400pouches
SB30440402-25B	4**4"-4ply	68*33*34cm	1000pouches
SB30330402-25B	3**3"-4ply	55*27*34cm	1000pouches
SB30220402-25B	2**2"-4ply	55*35*34cm	2000pouches
SB30480402-12B	4**8"-4ply	53*28*53cm	480pouches
SB30440402-12B	4**4"-4ply	53*28*33cm	960pouches
SB30330402-12B	3**3"-4ply	45*28*33cm	960pouches
SB30220402-12B	2**2"-4ply	53*35*41cm	1920pouches

09/30G/M2, 10PCS/POUCH, BLISTER POUCH

Code no	Model	Carton size	Q'ty(pks/ctn)
SB30480410-10B	4**8"-4ply	52*26*28cm	80pouches
SB30440410-11B	4**4"-4ply	55*30*30cm	160pouches
SB30330410-12B	3**3"-4ply	52*28*30cm	200pouches
SB30220410-13B	2**2"-4ply	52*35*30cm	400pouches
SB30480610-10B	4**8"-6ply	68*21*28cm	50pouches
SB30440610-10B	4**4"-6ply	68*21*33cm	100pouches
SB30330610-10B	3**3"-6ply	53*21*51cm	200pouches
SB30220610-10B	2**2"-6ply	48*21*43cm	250pouches

Non Sterile Non-Woven Sponge

- ◆ Made of spunlace non woven material, 70% viscose +30% polyester
- ◆ Weight: 30, 35, 40,50gsm/sq
- ◆ With or without x-ray detectable
- ◆ 4ply, 6ply, 8ply, 12ply
- ◆ 5x5cm,7.5x7.5cm,10x10cm,10x20cm etc
- ◆ 60pcs, 100pcs, 200pcs/pack(Non-sterile)



01/40G/M2, 200PCS OR 100PCS/PAPER BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
B404812-60	4**8"-12ply	52*48*42cm	20
B404412-60	4**4"-12ply	52*48*52cm	50
B403312-60	3**3"-12ply	40*48*40cm	50
B402212-60	2**2"-12ply	48*27*27cm	50
B404808-100	4**8"-8ply	52*28*42cm	10
B404408-100	4**4"-8ply	52*28*52cm	25
B403308-100	3**3"-8ply	40*28*40cm	25
B402208-100	2**2"-8ply	52*28*27cm	50
B404806-100	4**8"-6ply	52*40*42cm	20
B404406-100	4**4"-6ply	52*40*52cm	50
B403306-100	3**3"-6ply	40*40*40cm	50
B402206-100	2**2"-6ply	40*27*27cm	50
B404804-100	4**8"-4ply	52*28*42cm	20
B404404-100	4**4"-4ply	52*28*52cm	50
B403304-100	3**3"-4ply	40*28*40cm	50
B402204-100	2**2"-4ply	28*27*27cm	50
B404804-200	4**8"-4ply	52*28*42cm	10
B404404-200	4**4"-4ply	52*28*52cm	25
B403304-200	3**3"-4ply	40*28*40cm	25
B402204-200	2**2"-4ply	28*27*27cm	25

02/30G/M2, 200PCS OR 100PCS/PAPER BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
B304812-100	4**8"-12ply	52*28*42cm	10
B304412-100	4**4"-12ply	52*28*52cm	25
B303312-100	3**3"-12ply	40*28*40cm	25
B302212-100	2**2"-12ply	28*27*27cm	25
B304808-100	4**8"-8ply	52*42*42cm	20
B304408-100	4**4"-8ply	52*42*52cm	50
B303308-100	3**3"-8ply	42*40*40cm	50
B302208-100	2**2"-8ply	42*27*27cm	50
B304806-100	4**8"-6ply	52*32*42cm	20
B304406-100	4**4"-6ply	52*32*52cm	50
B303306-100	3**3"-6ply	40*32*40cm	50
B302206-100	2**2"-6ply	32*27*27cm	50
B304804-100	4**8"-4ply	52*42*42cm	40
B304404-100	4**4"-4ply	52*42*52cm	100
B303304-100	3**3"-4ply	40*42*40cm	100
B302204-100	2**2"-4ply	42*27*27cm	100
B304804-200	4**8"-4ply	52*42*42cm	20
B304404-200	4**4"-4ply	52*42*52cm	50
B303304-200	3**3"-4ply	40*42*40cm	50
B302204-200	2**2"-4ply	42*27*27cm	50

Gamgee Dressing

- ◆ Material:100% cotton(Sterile and Non sterile)
- ◆ Size: 7*10cm, 10*10cm,10*20cm,20*25cm,35*40cm or customized.
- ◆ Weight of cotton:200gsm/300gsm/350gsm/400gsm or customized
- ◆ Type: non selvage/single selvage/double selvage
- ◆ Sterilization method:Gamma ray/EO gas/Steam

PACKING REFERENCE FOR SOME SIZES:

Code no	Model	Packing	Carton size
SUGD1010S	10*10cm sterile	1pc/pack, 10packs/bag, 60bags/ctn	42*28*36cm
SUGD1020S	10*20cm sterile	1pc/pack, 10packs/bag, 24bags/ctn	48*24*32cm
SUGD2025S	20*25cm sterile	1pc/pack, 10packs/bag, 20bags/ctn	48*30*38cm
SUGD3540S	35*40cm sterile	1pc/pack, 10packs/bag, 6bags/ctn	66*22*37cm
SUGD0710N	7*10cm non sterile	100pcs/bag, 20bags/ctn	37*40*35cm
SUGD1323N	13*23cm non sterile	50pcs/bag, 16bags/ctn	54*46*35cm
SUGD1020N	10*20cm non sterile	50pcs/bag, 20bags/ctn	52*40*52cm
SUGD2020N	20*20cm non sterile	25pcs/bag, 20bags/ctn	52*40*35cm
SUGD3030N	30*30cm non sterile	25pcs/bag, 8bags/ctn	62*30*35cm



Gauze Ball

- ◆ Sterile and Non sterile
- ◆ 100% cotton, high absorbency and softness
- ◆ Cotton yarn of 21's, 32's, 40's
- ◆ Mesh of 20, 17 threads etc
- ◆ With or without x-ray detectable, elastic ring
- ◆ 8x8cm, 9x9cm, 15x15cm, 18x18cm, 20x20cm,
- ◆ 25x30cm, 30x40cm, 35x40cm ect
- ◆ Package: 100pcs/polybag(Non-sterile)
- ◆ 5pcs, 10pcs packed into blister pouch(Sterile)
- ◆ Gamma, EO, Steam

02/40S, 24X20 MESH, WITH OR WITHOUT X-RAY LINE, WITH OR WITHOUT RUBBER RING, 100PCS/PE-BAG

Code no	Model	Carton size	Q'ty(pcs/ctn)
E1712	8*8cm	58*30*38cm	30000
E1716	9*9cm	58*30*38cm	20000
E1720	15*15cm	58*30*38cm	10000
E1725	18*18cm	58*30*38cm	8000
E1730	20*20cm	58*30*38cm	6000
E1740	25*30cm	58*30*38cm	5000
E1750	30*40cm	58*30*38cm	4000



Krinkle gauze bandage

- ◆ Material: 100% cotton
- ◆ Size: 3.4"x3.6yards-6ply,4.6"x4.1yards-6ply
- ◆ Feature: Sterile, Soft Pouch Ideal for multiple wound care applications
- ◆ Packing: Blister pack or Vacuum pack

Code no	Size	Packing	Carton size
SUKGB4641	4.6"x4.1yards-6ply	1roll/blister, 100rolls/ctn	50*35*26cm
SUKGB4541	4.5"x4.1yards-6ply	1roll/blister, 100rolls/ctn	50*35*26cm



Cotton Roll

- ◆ Material: 100% pure cotton
- ◆ Packing: 1roll/blue kraft paper or polybag
- ◆ It is suit for medical and daily use.
- ◆ Type: ordinary, pre-cut

Code no	Specification	Packing	Carton size
SUCTR25G	25g/roll	500 rolls/ctn	56*36*56cm
SUCTR40G	40g/roll	400 rolls/ctn	56*37*56cm
SUCTR50G	50g/roll	300 rolls/ctn	61*37*61cm
SUCTR80G	80g/roll	200 rolls/ctn	61*31*61cm
SUCTR100G	100g/roll	200 rolls/ctn	61*31*61cm
SUCTR125G	125g/roll	100 rolls/ctn	61*36*36cm
SUCTR200G	200g/roll	50 rolls/ctn	41*41*41cm
SUCTR250G	250g/roll	50 rolls/ctn	41*41*41cm
SUCTR400G	400g/roll	40 rolls/ctn	55*31*36cm
SUCTR454G	454g/roll	40 rolls/ctn	61*37*46cm
SUCTR500G	500g/roll	20 rolls/ctn	61*38*48cm
SUCTR1000G	1000g/roll	20 rolls/ctn	66*34*52cm



Cutting Cotton Roll

- ◆ Material: 100% cotton
- ◆ Size: 100G(φ13x4.5cm), 250G(φ17x7cm)
- ◆ Packing: 1 roll pressed and packed into color bag.
- ◆ Feature: Cutting, pressed, popular in South America.

Code no	Specification	Packing	Carton size
SUCCR100G	100G	150 rolls/ctn	67*41*47cm
SUCCR250G	250G	60 rolls/ctn	70*37*53cm



Cotton Ball

- ◆ 100% pure cotton
- ◆ Sterile and non-sterile
- ◆ Color: white, red, blue, pink, green etc
- ◆ Weight: 0.5g, 1.0g, 1.5g, 2.0g, 3g etc

Code no	Specification	Packing
SUCTB001	0.5g	100pcs/bag 200bags/ctn
SUCTB002	1g	100pcs/bag 100bags/ctn
SUCTB003	2g	100pcs/bag 50bags/ctn
SUCTB004	3.5g	100pcs/bag 20bags/ctn
SUCTB005	5g	100pcs/bag 10bags/ctn
SUCTB006	0.5g	5pcs/blister, 20blister/bag 20bags/ctn
SUCTB007	1g	5pcs/blister, 20blister/bag 10bags/ctn
SUCTB008	2g	5pcs/blister, 10blister/bag 10bags/ctn
SUCTB009	3.5g	5pcs/blister, 10blister/bag 10bags/ctn
SUCTB010	5g	5pcs/blister, 10blister/bag 10bags/ctn



Zigzag Cotton

- ◆ Material: 100% pure cotton
- ◆ Color: natural white
- ◆ Packing: 1roll/polybag
- ◆ It is suit for medical and daily use

Code no	Specification	Packing	Carton size
SUZZC25G	25g/roll	500 rolls/ctn	66*48*53cm
SUZZC50G	50g/roll	200 rolls/ctn	59*46*48cm
SUZZC100G	100g/roll	120 rolls/ctn	59*46*48cm
SUZZC200G	200g/roll	80 rolls/ctn	59*46*66cm
SUZZC250G	250g/roll	30 rolls/ctn	50*30*47cm



Dental Cotton Roll

- ◆ 100% pure cotton
- ◆ 50pcs/pack, 20packs/bag

Code no	Size	Packing	Carton size
SUDCR001	8mmx3.8cm	20bags/ctn	50*32*40cm
SUDCR002	10mmx3.8cm	20bags/ctn	60*38*40cm
SUDCR003	12mmx3.8cm	10bags/ctn	43*37*40cm
SUDCR004	14mmx3.8cm	10bags/ctn	50*32*40cm



Cotton Swab

- ◆ Code No.:SUCTS001
- ◆ Material: 100% pure cotton
- ◆ Apply: the ears, nose, skin, clean and makeup, beauty
- ◆ Type: Sterile and non-sterile
- ◆ Wooden stick or plastic stick
- ◆ Diameter: 0.5mm, 1mm, 2mm, 2.5mm etc
- ◆ Stick length: 7.5cm, 10cm or 15cm etc
- ◆ Package: 100pcs/polybag(Non-sterile)
- ◆ 3pcs, 5pcs, 10pcs packed into pouch(Sterile)
- ◆ Sterile way: EO GAS



Povidone Iodine Swabstick

- ◆ Code No.:SUPLS002
- ◆ Material: 100% combed cotton+plastic stick
- ◆ Main ingredients:saturated with 10% povidone-iodine, 1% available iodine
- ◆ Type: Sterile
- ◆ Size: 10cm
- ◆ Diameter: 10mm
- ◆ Package: 1pc/pouch, 50bags/box, 1000bags/ctn, ctn size:44x31x35cm
3pcs/pouch, 25bags/box, 500bags/ctn, ctn size:44x31x35cm



Cotton Pad

- ◆ Code No.:SUCTP003
- ◆ Material: 100% cotton
- ◆ Size: diameter 58mm or other sizes
- ◆ Color: white
- ◆ Type: Round, plain, press/unpress
- ◆ Grams: 0.3g-0.65g
- ◆ Package: string bag or zipbag



Waxing(Muslin) Roll/Strips

- ◆ Code No.:SUWXR004
- ◆ Material: non woven, cotton
- ◆ Weight: 70gsm - 90gsm
- ◆ Application: Body hair removal, beauty salon and Spa
- ◆ Size for different styles:
Muslin roll: 3"x100yards; 3.25"x40yards; 3.5"x40yards; 3"x50yards etc
Muslin strips: 7x20cm; 7.5x22.5cm; 7x23cm etc
- ◆ Type: bleached, unbleached; straight, zig zag cut



Makeup Brush

- ◆ Code No.:SUMKB005
- ◆ Material: Bristle, brass, wooden handle
- ◆ Size: 7.5", 6" etc
- ◆ Packing:
1pc/opp bag, 100pcs/middle bag, 1000pcs/ctn
1pc/opp bag, 100pcs/middle bag, 500pcs/ctn



Mascara Brush

- ◆ Code No.:SUMCB006
- ◆ Material: injection sticks + nylon brush head
- ◆ Specifications: brush head 2.6cm, Total length 9.8cm
- ◆ Color: black



Wax Applicator

- ◆ Material: wood
- ◆ Specifications

Code no	Item	Packing	Carton size
SUWAT001	150x18x1.6mm	500pcs/box, 10boxes/ctn	50.5x32.5x20.5cm
SUWAT002	114x10x2mm	1000pcs/box, 10boxes/ctn	62x45x15cm
SUWAT003	3.8x100mm	100pcs/bag, 10bags/box, 20boxes/ctn	59x43.5x18.5cm



Crepe Bandage

- ◆ Material: 97% cotton+3% spandex
- ◆ Color: White, skin, with aluminium clip or elastic clip
- ◆ Weight: 70g, 75g, 80g, 85g, 90g, 95g, 100g etc
- ◆ With or without red/blue line
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc
- ◆ Length: 10m, 10yards, 5m, 5yards, 4m, 4yards etc
- ◆ Packing: 1roll/individually packed

Code no	Size	Packing	Carton size
SUCPB0545	5cmx4.5m	960rolls/ctn	54*32*44cm
SUCPB7545	7.5cmx4.5m	480rolls/ctn	54*32*44cm
SUCPB1045	10cmx4.5m	360rolls/ctn	54*32*44cm
SUCPB1545	15cmx4.5m	240rolls/ctn	54*32*44cm
SUCPB2045	20cmx4.5m	120rolls/ctn	54*32*44cm



100% Cotton Crepe Bandage

- ◆ Material: 100% cotton
- ◆ Color: White, with aluminium clip or elastic clip
- ◆ Weight: 80g, 85g, 90g, 95g, 100g, 105g etc
- ◆ With or without red/blue line
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc
- ◆ Length: 10m, 10yards, 5m, 5yards, 4m, 4yards etc
- ◆ Packing: 1roll/individually packed

Code no	Size	Packing	Carton size
SUCCB0545	5cmx4.5m	960rolls/ctn	54*32*44cm
SUCCB7545	7.5cmx4.5m	480rolls/ctn	54*32*44cm
SUCCB1045	10cmx4.5m	360rolls/ctn	54*32*44cm
SUCCB1545	15cmx4.5m	240rolls/ctn	54*32*44cm
SUCCB2045	20cmx4.5m	120rolls/ctn	54*32*44cm



High Elastic Bandage

- ◆ Material: Polyester/cotton; rubber/spandex
- ◆ Color: light skin/dark skin/natural white etc
- ◆ Weight: 80g, 85g, 90g, 100g, 105g, 110g, 120g etc
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc
- ◆ Length: 5m, 5yards, 4m etc
- ◆ With latex or latex free
- ◆ Packing: 1roll/individually packed

Code no	Size	Packing	Carton size
SUHEB0545	5cmx4.5m	960rolls/ctn	54*43*44cm
SUHEB7545	7.5cmx4.5m	480rolls/ctn	54*32*44cm
SUHEB1045	10cmx4.5m	480rolls/ctn	54*42*44cm
SUHEB1545	15cmx4.5m	240rolls/ctn	54*32*44cm
SUHEB2045	20cmx4.5m	120rolls/ctn	54*42*44cm



PBT(Conforming)Bandage

- ◆ Material: viscose, cotton, polyamide
- ◆ Color: white
- ◆ Weight: 30g, 40g, 45g, 50g, 55g etc
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc
- ◆ Length: 5m, 5yards, 4m, 4yards etc
- ◆ Packing: 1roll/individually packed

Code no	Size	Packing	Carton size
SUPTB0545	5cmx4.5m	720rolls/ctn	43*35*36cm
SUPTB7545	7.5cmx4.5m	480rolls/ctn	43*35*36cm
SUPTB1045	10cmx4.5m	360rolls/ctn	43*35*36cm
SUPTB1545	15cmx4.5m	240rolls/ctn	43*35*36cm
SUPTB2045	20cmx4.5m	120rolls/ctn	43*35*36cm



POP Bandage

- ◆ Material: viscose, cotton, polyamide
- ◆ Weight: 30g, 40g, 45g, 50g, 55g etc
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc
- ◆ Length: 5m, 5yards, 4m, 4yards etc
- ◆ Packing: 1roll/individually packed

Code no	Size	Packing	Carton size
SUPPB0527	5cmx2.7m	240rolls/ctn	57*33*26cm
SUPPB7527	7.5cmx2.7m	240rolls/ctn	57*33*26cm
SUPPB1027	10cmx2.7m	120rolls/ctn	57*33*26cm
SUPPB12527	12.5cmx2.7m	120rolls/ctn	57*33*26cm
SUPPB1527	15cmx2.7m	120rolls/ctn	57*33*26cm
SUPPB2027	20cmx2.7m	60rolls/ctn	57*33*26cm



Under Cast Padding For POP

- ◆ Material: 100% polyester or 100% cotton
- ◆ Color: white
- ◆ Weight: 60-140gsm etc
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm, 20cm etc
- ◆ Length: 2.7m, 3m, 3.6m, 4m, 4.5m, 5m etc
- ◆ Packing: 1roll/individually packed

Code no	Size	Packing	Carton size
SUUCP0527	5cmx2.7m	720rolls/ctn	66*33*48cm
SUUCP7527	7.5cmx2.7m	480rolls/ctn	66*33*48cm
SUUCP1027	10cmx2.7m	360rolls/ctn	66*33*48cm
SUUCP1527	15cmx2.7m	240rolls/ctn	66*33*48cm
SUUCP2027	20cmx2.7m	120rolls/ctn	66*33*48cm



Selvage Gauze Bandage

- ◆ Material: 100% cotton
- ◆ Mesh: 30x20, 24x20 etc
- ◆ Width: 5cm, 7.5cm, 10cm, 12cm, 15cm etc
- ◆ Length: 10m, 10yards, 5m, 5yards, 4m etc
- ◆ Packing: 1roll/polybag

Code no	Size	Packing	Carton size
SUSGB055	5cmx5m	960rolls/ctn	36*30*43cm
SUSGB065	6cmx5m	880rolls/ctn	36*30*46cm
SUSGB755	7.5cmx5m	1080rolls/ctn	50*33*41cm
SUSGB085	8cmx5m	720rolls/ctn	36*30*52cm
SUSGB105	10cmx5m	480rolls/ctn	36*30*43cm
SUSGB125	12cmx5m	480rolls/ctn	36*30*50cm
SUSGB155	15cmx5m	360rolls/ctn	36*32*45cm



Heavy Elastic Adhesive Bandage(EAB)

- ◆ Material: 100% cotton elastic fabric
- ◆ Color: White with yellow middle line etc
- ◆ Width: 5cm, 7.5cm, 10cm, 15cm etc
- ◆ Length: 4.5m etc
- ◆ Glue: Hot melt adhesive, latex free

Code no	Size	Packing	Carton size
SUHEB0545	5cmx4.5m	1roll/polybag, 216rolls/ctn	50*38*38cm
SUHEB7545	7.5cmx4.5m	1roll/polybag, 144rolls/ctn	50*38*38cm
SUHEB1045	10cmx4.5m	1roll/polybag, 108rolls/ctn	50*38*38cm
SUHEB1545	15cmx4.5m	1roll/polybag, 72rolls/ctn	50*38*38cm



Adhesive Elastic Bandage

- ◆ Material: non woven/cotton
- ◆ Color: blue, red, green, yellow etc
- ◆ Width: 2.5cm, 5m, 7.5cm, 10cm etc
- ◆ Length: 5m, 5yards, 4m, 4yards, 3m etc
- ◆ Packing: 1roll/candybag or blister

Code no	Size	Packing	Carton size
SUAEB0545	5cmx4.5m	1roll/polybag, 216rolls/ctn	50*38*38cm
SUAEB7545	7.5cmx4.5m	1roll/polybag, 144rolls/ctn	50*38*38cm
SUAEB1045	10cmx4.5m	1roll/polybag, 108rolls/ctn	50*38*38cm
SUAEB1545	15cmx4.5m	1roll/polybag, 72rolls/ctn	50*38*38cm



Orthopedic Casting Tape

- ◆ Material: fiberglass/polyester
- ◆ Color: red, blue, yellow, pink, green, purple etc
- ◆ Width: 5cmx4yards, 7.5cmx4yards, 10cmx4yards, 12.5cmx4yards, 15cmx4yards

Code no	Size	Packing	Carton size
SUOCT054	5cmx4yards	10pcs/box, 16boxes/ctn	55.5*49*44cm
SUOCT754	7.5cmx4yards	10pcs/box, 12boxes/ctn	55.5*49*44cm
SUOCT104	10cmx4yards	10pcs/box, 10boxes/ctn	55.5*49*44cm
SUOCT154	15cmx4yards	10pcs/box, 8boxes/ctn	55.5*49*44cm
SUOCT204	20cmx4yards	10pcs/box, 8boxes/ctn	55.5*49*44cm



Net Bandage

- ◆ Material: Polymide+rubber, nylon+latex
- ◆ Length: normal 25m after stretched
- ◆ Width: 0.6cm, 1.7cm, 2.2cm, 3.8cm, 4.4cm, 5.2cm etc
- ◆ Package: 1pc/box

Code no	Size	Packing	Carton size
SUNTB001	0.5, 0.7cm x 25m	1pc/box, 180boxes/ctn	68*38*28cm
SUNTB002	1.0, 1.7cm x 25m	1pc/box, 120boxes/ctn	68*38*28cm
SUNTB003	2.0, 2.0cm x 25m	1pc/box, 120boxes/ctn	68*38*28cm
SUNTB004	3.0, 2.3cm x 25m	1pc/box, 84boxes/ctn	68*38*28cm
SUNTB005	4.0, 3.0cm x 25m	1pc/box, 84boxes/ctn	68*38*28cm
SUNTB006	5.0, 4.2cm x 25m	1pc/box, 56boxes/ctn	68*38*28cm
SUNTB007	6.0, 5.8cm x 25m	1pc/box, 32boxes/ctn	68*38*28cm

Tubular Bandage

Code no	Size	Carton size	Packing	G.W/kg	N.W/kg
SUTB055	5cmx5m	33*38*30cm	72rolls/ctn	8.5	6.5
SUTB755	7.5cmx5m	33*38*30cm	48rolls/ctn	8.5	6.5
SUTB105	10cmx5m	33*38*30cm	36rolls/ctn	8.5	6.5
SUTB155	15cmx5m	33*38*30cm	24rolls/ctn	8.5	6.5
SUTB205	20cmx5m	42*30*30cm	18rolls/ctn	8.5	6.5
SUTB255	25cmx5m	28*47*30cm	15rolls/ctn	8.8	6.8
SUTB0510	5cmx10m	54*28*29cm	40rolls/ctn	9.2	7.2
SUTB7510	7.5cmx10m	41*41*29cm	30rolls/ctn	10.1	8.1
SUTB1010	10cmx10m	54*28*29cm	20rolls/ctn	9.2	7.2
SUTB1510	15cmx10m	54*33*29cm	16rolls/ctn	10.6	8.6
SUTB2010	20cmx10m	54*46*29cm	16rolls/ctn	13.5	11.5
SUTB2510	25cmx10m	54*41*29cm	12rolls/ctn	12.8	10.8
SUTB0525	5cmx25m	46*28*46cm	20rolls/ctn	11	9
SUTB7525	7.5cmx25m	46*33*46cm	16rolls/ctn	12.8	10.8
SUTB1025	10cmx25m	46*33*46cm	12rolls/ctn	12.8	10.8
SUTB1525	15cmx25m	46*33*46cm	8rolls/ctn	12.8	10.8
SUTB2025	20cmx25m	46*23*46cm	4rolls/ctn	9.2	7.2
SUTB2525	25cmx25m	46*28*46cm	4rolls/ctn	11	9

Triangle Bandage

- ◆ Material: 100% cotton or non woven fabric
- ◆ With or without safety pin
- ◆ Size: 36"x36"x51", 40"x40"x56" etc
- ◆ Cotton year: 40x34, 50x30, 48x48 etc
- ◆ Color: Unbleached or bleached

Code no	Size	Packing	Carton size
SUTGB127	90x90x127cm	1pc/polybag, 500pcs/ctn	52*32*42cm
SUTGB140	100x100x140cm	1pc/polybag, 500pcs/ctn	52*36*42cm



First Aid Kit (Code No.:SUFAB001)

We provide wide range of first aid kits available for you different needs, below is a sample list:

Item name	Specification	Quantity	Item name	Specification	Quantity
Adhesive bandage	72x19mm	12	First aid blanket	204x140cm	1
Iodine cotton bar	1pc/bag	24	Triangular bandage	90x90x129cm	1
Absorbent adhesive dressing	6x7cm	5	PBT elastic bandage	10x450cm	1
Absorbent adhesive dressing	10x10cm	5	Adhesive tape	1cmx10m	1
Dressing pad	5x5cm	5	Safety pin		4
Dressing pad	7.5x7.5cm	5	Month-to-mouth mask	20x20cm	1
Dressing pad	10x10cm	4	Instant ice bag	100g	1
Scissor	13.5cm	1	thermometer		1
Tweezer	12.5cm	1	First aid booklet		1
Iodine cotton ball	5pcs/bag	1	First aid instruction		1
Alcohol pad	5x5cm	4	First aid bag	21x14.5x6.5cm	1



First Aid Blanket

- ◆ Material: PET also named emergency blanket
- ◆ Color: gold silver/Silver Silver
- ◆ Size: 160*210cm, 140*210cm or custom size
- ◆ Feature: windproof, waterproof & against coldness

Code no	Item	Size	Packing	Carton size
SUFAB160210G	Gold/Silver blanket	160x210cm	1pc/PE bag, 200pcs/ctn	50*30*30cm
SUFAB140210G	Gold/Silver blanket	140x210cm	1pc/PE bag, 200pcs/ctn	50*30*30cm
SUFAB160210S	Silver/Silver blanket	160x210cm	1pc/PE bag, 200pcs/ctn	50*30*30cm
SUFAB160210S	Silver/Silver blanket	140x210cm	1pc/PE bag, 200pcs/ctn	50*30*30cm



First Aid Bandage

- ◆ Material: Bandage and pad; Bandage: Conforming bandage or gauze bandage
- ◆ Pad: Low adherent dressing, cotton covered with gauze or non woven
- ◆ Basic packing: Individually packed in cellophane

Code no	Spec.	Packing	Carton size
SUFAB064M	6cm*4m	1roll/bag, 600rolls/ctn	62*24*40cm
SUFAB084M	8cm*4m	1roll/bag, 480rolls/ctn	66*24*40cm
SUFAB104M	10cm*4m	1roll/bag, 360rolls/ctn	62*24*40cm



Aperture Zinc Oxide Plaster

- ◆ Material: 100% cotton fabric
- ◆ Colors: White/skin
- ◆ Glue: natural zinc oxide glue
- ◆ Package: 1roll/box
- ◆ Width: 18cm, 10cm etc
- ◆ Length: 10m, 10yards, 5m, 5yards etc

Code no	Size	Packing	Carton size
SUAZOP185M	18cmx5m	1roll/box, 30rolls/ctn	37.5*31.5*21cm
SUAZOP185Y	18cmx5yds	1roll/box, 30rolls/ctn	37.5*31.5*21cm
SUAZOP105M	10cmx5m	1roll/box, 60rolls/ctn	37.5*31.5*24.5cm
SUAZOP105Y	10cmx5yds	1roll/box, 60rolls/ctn	37.5*31.5*24.5cm



Zinc Oxide Adhesive Tape

- ◆ Material: 100% cotton fabric
- ◆ Colors: White/skin
- ◆ Glue: natural zinc oxide glue / hot-melt glue
- ◆ Package: Simple package / Plastic tin / Iron tin
- ◆ Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
- ◆ Length: 10m, 10yards, 5m, 5yards etc

Code no	Size	Packing	Carton size
SUZOT1255M	1.25cmx5m	48rolls/box, 12boxes/ctn	39*37*39cm
SUZOT255M	2.5cmx5m	30rolls/box, 12boxes/ctn	39*37*39cm
SUZOT055M	5cmx5m	18rolls/box, 12boxes/ctn	39*37*39cm
SUZOT755M	7.5cmx5m	12rolls/box, 12boxes/ctn	39*37*39cm
SUZOT105M	10cmx5m	9rolls/box, 12boxes/ctn	39*37*39cm
SUZOT12510Y	1.25cmx9.14m	48rolls/box, 12boxes/ctn	39*37*39cm
SUZOT2510Y	2.5cmx9.14m	30rolls/box, 12boxes/ctn	39*37*39cm
SUZOT0510Y	5cmx9.14m	18rolls/box, 12boxes/ctn	39*37*39cm
SUZOT7510Y	7.5cmx9.14m	12rolls/box, 12boxes/ctn	39*37*39cm
SUZOT1010Y	10cmx9.14m	9rolls/box, 12boxes/ctn	39*37*39cm



Non-woven Tape

- ◆ Material: non woven
- ◆ Color: Skin or white color
- ◆ Glue: hot-melt or acrylic acid glue
- ◆ Package: Simple package / Plastic tin / Iron tin
- ◆ Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
- ◆ Length: 10m, 10yards, 5m, 5yards etc

Code no	Size	Packing	Carton size
SUNT1255Y	1.25cmx5yds	48rolls/box, 12boxes/ctn	24*23.5*28.5cm
SUNT255Y	2.5cmx5yds	30rolls/box, 12boxes/ctn	24*23.5*28.5cm
SUNT055Y	5cmx5yds	18rolls/box, 12boxes/ctn	24*23.5*28.5cm
SUNT765Y	7.5cmx5yds	12rolls/box, 12boxes/ctn	24*23.5*41cm
SUNT105Y	10cmx5yds	9rolls/box, 12boxes/ctn	38.5*23.5*33.5cm
SUNT12510M	1.25cmx10m	48rolls/box, 12boxes/ctn	24*23.5*28.5cm
SUNT2510M	2.5cmx10m	30rolls/box, 12boxes/ctn	24*23.5*28.5cm
SUNT0510M	5cmx10m	18rolls/box, 12boxes/ctn	24*23.5*28.5cm
SUNT7510M	7.5cmx10m	12rolls/box, 12boxes/ctn	24*23.5*41cm
SUNT1010M	10cmx10m	9rolls/box, 12boxes/ctn	38.5*23.5*33.5cm



Silk Tape

- ◆ Material: Acetate silk fabric
- ◆ Color: Skin / white
- ◆ Glue: Hot melt glue
- ◆ Package: Simple packing/Plastic tin/Blister packing
- ◆ Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
- ◆ Length: 10m, 10yards, 5m, 5yards etc

Code no	Size	Packing	Carton size
SUSKT1255Y	1.25cmx4.5m	24rolls/box, 30boxes/ctn	39*18*29cm
SUSKT255Y	2.5cmx4.5m	12rolls/box, 30boxes/ctn	39*18*29cm
SUSKT055Y	5cmx4.5m	6rolls/box, 30boxes/ctn	39*18*29cm
SUSKT755Y	7.5cmx4.5m	6rolls/box, 20boxes/ctn	43*26.5*26cm
SUSKT105Y	10cmx4.5m	6rolls/box, 20boxes/ctn	43*26.5*26cm



PE Tape

- ◆ Material: PE film
- ◆ Color: transparent
- ◆ Glue: Hot melt or acrylic acid glue
- ◆ Package: Simple packing/Plastic tin/with dispenser
- ◆ Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
- ◆ Length: 10m, 10yards, 5m, 5yards etc

Code no	Size	Packing	Carton size
SUPET1255Y	1.25cmx5yards	24rolls/box, 30boxes/ctn	39*18.5*29cm
SUPET255Y	2.5cmx5yards	12rolls/box, 30boxes/ctn	39*18.5*29cm
SUPET055Y	5cmx5yards	6rolls/box, 30boxes/ctn	39*18.5*29cm
SUPET755Y	7.5cmx5yards	6rolls/box, 30boxes/ctn	44*26.5*26cm
SUPET105Y	10cmx5yards	6rolls/box, 30boxes/ctn	44*26.5*33.5cm
SUPET1255M	1.25cmx5m	24rolls/box, 30boxes/ctn	39*18.5*29cm
SUPET255M	2.5cmx5m	12rolls/box, 30boxes/ctn	39*18.5*29cm
SUPET055M	5cmx5m	6rolls/box, 30boxes/ctn	39*18.5*29cm
SUPET755M	7.5cmx5m	6rolls/box, 30boxes/ctn	44*26.5*26cm
SUPET105M	10cmx5m	6rolls/box, 30boxes/ctn	44*26.5*33.5cm



Kinesiology Tape

- ◆ Material: cotton+spandex
- ◆ Color: red, blue, orange, green, black, skin, etc
- ◆ Adhesive: Rubber
- ◆ Width: 2.5cm, 3.8cm, 5cm, 7.5cm, 10cm etc
- ◆ Length: 5m, 5yards etc
- ◆ Package: Simple packing, colorful box packing

Code no	Size	Packing	Carton size
SUKT1255Y	1.25cmx4.5m	24rolls/box, 30boxes/ctn	39*18*29cm
SUKT255Y	2.5cmx4.5m	12rolls/box, 30boxes/ctn	39*18*29cm
SUKT055Y	5cmx4.5m	6rolls/box, 30boxes/ctn	39*18*29cm
SUKT755Y	7.5cmx4.5m	6rolls/box, 20boxes/ctn	43*26.5*26cm
SUKT105Y	10cmx4.5m	6rolls/box, 20boxes/ctn	43*26.5*26cm



Sport Tape

- ◆ Material: 100% cotton
- ◆ Color: Beige, black, red, pink, blue, green etc
- ◆ Width: 2.5cm, 3.8cm, 5cm, 7.5cm etc
- ◆ Length: 5m, 5yards, 4m, 6m etc
- ◆ Feature: Latex free, waterproof

Code no	Size	Packing	Carton size
SUSPT1255Y	1.25cmx4.5m	24rolls/box, 30boxes/ctn	39*18*29cm
SUSPT255Y	2.5cmx4.5m	12rolls/box, 30boxes/ctn	39*18*29cm
SUSPT055Y	5cmx4.5m	6rolls/box, 30boxes/ctn	39*18*29cm
SUSPT755Y	7.5cmx4.5m	6rolls/box, 20boxes/ctn	43*26.5*26cm
SUSPT105Y	10cmx4.5m	6rolls/box, 20boxes/ctn	43*26.5*26cm



Capsicum Plaster

- ◆ Material: cotton+glue
- ◆ Color: skin or white
- ◆ Packing: 24pcs/box
- ◆ Size: 10*14cm, also 8*12cm/10*20cm on applying

Code no	Size	Packing	Carton size
SUCSP1218	12x18cm	24pcs/box, 1200pcs/ctn	62*35*25cm
SUCSP0610	6x10cm	24pcs/box, 2400pcs/ctn	62*35*25cm



Skin Closure Strip

- ◆ Material: Non woven fabric reinforcement
- ◆ Color: transparent
- ◆ Size: 3x75mm, 6x38mm, 6x75mm, 6x100mm, 12x100mm
- ◆ Packaging: 100pcs/box 10boxes/carton for Skin Closure Strip

Code no	Size	Packing
SUSCS001	1/8in.x3in./3x75mm	5 strips/pouch, 50pouches/box, 10box/ctn
SUSCS002	1/4in.x1-1/2in./6x38mm	6 strips/pouch, 50pouches/box, 10box/ctn
SUSCS003	1/4in.x3in./6x75mm	3 strips/pouch, 50pouches/box, 10box/ctn
SUSCS004	1/4in.x4in./6x100mm	10 strips/pouch, 50pouches/box, 10box/ctn
SUSCS005	1/2in.x4in./12x100mm	6 strips/pouch, 50pouches/box, 10box/ctn



Cooling Patch

- ◆ Material: Glycerin, water, bean dilute acid sodium etc
- ◆ Size: 4x9cm, 4x11cm, 5x10cm, 5x12cm etc
- ◆ Color: blue, green or other colors
- ◆ Packing: 1pc/bag, 3pcs/box or 4pcs/bag etc
- ◆ Application: Physical cooling

Code no	Size	Packing
SUCLP0512	5x12cm	1pc/foil bag, 3pcs/box, 144boxes/ctn
SUCLP0411	4x11cm	1pc/foil bag, 4pcs/box, 120boxes/ctn



Acne Plaster

Code No.:SUAP001
Material: transparent PE film+glue
Size: Dia 12mm/8mm Thickness:0.4mm
Package: 1pc, 8pcs, 12pcs/sheet, 36pcs, 50pcs/box 60boxes, 144boxes/ctn

- ◆ Suitable for all kinds of small wounds.
- ◆ Excellent absorbency, it can clean secretion from follicle.
- ◆ Reduce wound inflammation, make the skin flat and cured gradually.
- ◆ It offers a moist environment, accelerate wound healing and avoid scars



Burn Dressing

- ◆ Sterile trauma hydrogel burn dressing, it's soothing, Cooling Gel for buns, safe for use on children&facial burns, protect from extreme temperatures.

Code no	Burncare	Specification	Pieces/Box	Box/Carton
SUBD001	5.030101	3.5g	6	120
SUBD002	5.030102	5x5	20	16
SUBD003	5.030104	5x15	12	16
SUBD004	5.030103	10x10	20	8
SUBD005	5.030105	20x20	6	8
SUBD006	5.030106	30x30	6	3
SUBD007	5.030107	30x40	6	16
SUBD008	5.030108	40x60	6	16
SUBD009	5.030151	120x160cm	6	2
SUBD010	5.030109	50spray	6	8
SUBD011	5.030112	59bottle	6	8
SUBD012	5.030113	118bottle	4	3
SUBD013	5.030110	120spray	1	50
SUBD014	5.030121	25tube	1	50
SUBD015	5.030122	50tube	1	50



Forehead Thermometer Strip

- ◆ Code No.:SUFTS001
- ◆ Material: PET Sheet + Liquid Crystal
- ◆ Size: 16*90mm or Die Cut
- ◆ Temperature range: 35, 36, 37, 38, 39, 40°C
- ◆ Display Resolution: 1°C
- ◆ Packing: 1pc/oppbag, 1pc/paper bag, or 1pc/box



Non Woven Wound Dressing

- ◆ Material: Made of spunlace non woven
- ◆ Package: 1pc/pouch, 50pouches/box
- ◆ Sterile way: EO sterile

Code no	Specification	Packing	Carton size
SUNWD0505	5x5cm	50pcs/box 2500pcs/ctn	50x20x45cm
SUNWD0507	5x7cm	50pcs/box 2500pcs/ctn	52x24x45cm
SUNWD0607	6x7cm	50pcs/box 2500pcs/ctn	52x24x50cm
SUNWD0608	6x8cm	50pcs/box 1200pcs/ctn	50x21x31cm
SUNWD0510	5x10cm	50pcs/box 1200pcs/ctn	42x35x31cm
SUNWD0610	6x10cm	50pcs/box 1200pcs/ctn	42x34x31cm
SUNWD7510	10x7.5cm	50pcs/box 1200pcs/ctn	42x34x37cm
SUNWD1010	10x10cm	50pcs/box 1200pcs/ctn	58x35x35cm
SUNWD1012	10x12cm	50pcs/box 1200pcs/ctn	57x42x29cm
SUNWD1015	10x15cm	50pcs/box 1200pcs/ctn	58x44x38cm
SUNWD1020	10x20cm	50pcs/box 600pcs/ctn	55x25x43cm
SUNWD1025	10x25cm	50pcs/box 600pcs/ctn	58x33x38cm
SUNWD1030	10x30cm	50pcs/box 600pcs/ctn	58x38x38cm



Wound Dressing Roll


- ◆ Code No.:SUWDR001
- ◆ Material: Made of spunlace non woven
- ◆ Size: 5cmx10m, 10cmx10m, 15cmx10m, 20cmx10m etc
- ◆ Color: white(mostly), green, blue
- ◆ Package: non sterile, 1roll/box



Transparent Dressing Film

- ◆ Material: Made of transparent PU film
- ◆ Package: 1pc/pouch, 50pouches/box
- ◆ Sterile way: EO sterile

Code no	Specification	Packing	Carton size
SUTDF0505	5x5cm	50pcs/box 2500pcs/ctn	50x20x45cm
SUTDF0507	5x7cm	50pcs/box 2500pcs/ctn	52x24x45cm
SUTDF0607	6x7cm	50pcs/box 2500pcs/ctn	52x24x50cm
SUTDF0608	6x8cm	50pcs/box 1200pcs/ctn	50x21x31cm
SUTDF0510	5x10cm	50pcs/box 1200pcs/ctn	42x35x31cm
SUTDF0610	6x10cm	50pcs/box 1200pcs/ctn	42x34x31cm
SUTDF7510	10x7.5cm	50pcs/box 1200pcs/ctn	42x34x37cm
SUTDF1010	10x10cm	50pcs/box 1200pcs/ctn	58x35x35cm
SUTDF1012	10x12cm	50pcs/box 1200pcs/ctn	57x42x29cm
SUTDF1015	10x15cm	50pcs/box 1200pcs/ctn	58x44x38cm
SUTDF1020	10x20cm	50pcs/box 600pcs/ctn	55x25x43cm
SUTDF1025	10x25cm	50pcs/box 600pcs/ctn	58x33x38cm
SUTDF1030	10x30cm	50pcs/box 600pcs/ctn	58x38x38cm



Adhesive Eye Pad

- ◆ Code No.:SUAEP001
- ◆ Material: Made of spunlace non woven
- ◆ Size: Adult 6.5x9.5cm, child 4.5x6.7cm
- ◆ Packing: 1pc/pouch, 50pouches/box
- ◆ Sterile way: EO sterile



Alcohol Prep Pad

- ◆ Code No.:SUAPP001
- ◆ Material: Non woven, 70% isopropyl alcohol
- ◆ Size: 3x6.5cm, 4x6cm, 5x5cm, 7.5x7.5cm etc
- ◆ Package: 1pc/pouch, 100,200pouches/box
- ◆ Sterile way: EO



Povidone-Iodine Prep Pad

- ◆ Code No.:SUPPP002
- ◆ Material: 1%Povidone iodine + non woven swab
- ◆ Size: 3x6cm
- ◆ Packing: 100pcs/box, 100boxes/ctn
- ◆ Tear the package directly to wipe the part to be disinfected




Wound Plaster(Band Aid)

- ◆ Code No.: SUWP003
- ◆ Material: PE, PVC, Cotton/Nylon Fabric
- ◆ Color: skin, transparent or cartoon etc
- ◆ Size: 72*19mm or other
- ◆ Packing: individual pack in color box
- ◆ Sterilized: EO
- ◆ Shapes: available in various sizes



Wound Plaster Strip

- ◆ Code No.:SUWPS004
- ◆ Material: PVC, PE, PU, or Fabric, Elastic Fabric available
- ◆ Width: 4cm, 6cm, 8cm
- ◆ Length: 10cm, 100cm, 500cm, other sizes available on request
- ◆ Color: Skin color
- ◆ Package: 10pcs/paper box, or 6pcs/box, or 1pc/box



IV Wound Dressing

- ◆ Code No.:SUIWD005
- ◆ Material: Waterproof PU Film & Medical acrylate adhesive
- ◆ Color: White transparent iv dressing
- ◆ Size: 10*15cm with absorbent pad
- ◆ Package: 100pcs/box,2000pcs/ctn

Specification	Packing	Carton size
6x8cm	100pcs/box 2000pcs/ctn	49x29x29cm



Face Mask

NON WOVEN FACE MASK (Code No.:SUNFM001)

Material	Non-woven PP material
Style	With elastic ear loop or tie on
Color	Blue, green, white, pink, etc.
Layer	Usually 3ply, 1ply 2ply and 4ply also available.
Weight	18gsm+20gsm+25gsm etc
Size	17.5x9.5cm, 14.5x9cm, 12.5x8cm
BFE	≥ 99% & 99.9%
Packaging	50pcs/box 40boxes/ctn



PAPER FACE MASK (Code No.:SUPFM002)

Material	Paper material
Style	With elastic ear loop
Color	white
Layer	1ply, 2ply
Size	20x7cm
Packaging	100pcs/box, 100boxes/ctn

N95 FACE MASK (Code No.:SUNFM003)

Material	Made of multi-layered nonpoisonous non-allergic, non-stimulating materials
Color	White
Valve	With or without exhalation valve
Style	Earloop
Size	Standard 132x115x47mm; Large 140x125x52mm
Standard	NIOSH N95
Shape	Cup



CONIC FACE MASK (Code No.:SUCFM004)

Material	Material Acupuncture cotton
Style	Style with elastic ear loop
Color	White
Valve	Valve Without
Packaging	50pcs/box, 20boxes/ctn



Face mask with shield

- ◆ Code No.:SUFMS005
- ◆ Material: PP non woven, Filter Paper, PVC
- ◆ Size: 17.5x9.5cm for adult, 14.5x9cm for child
- ◆ Color: Green, Blue, Pink, White etc
- ◆ Style: Ear loop or Tie on; 1ply, 2ply, 3ply, 4ply
- ◆ Usage: Dust-proof, anti-virus, anti-bacteria, blood-proof, food-protection when cooking



Cap

1.DOCTOR CAP (Code No.:SUDTC001)

- ◆ Material: PP non woven/SMS
- ◆ Weight: 20gsm, 25gsm, 30gsm etc
- ◆ Type: with tie or elastic
- ◆ Size: 62*12.5cm/63*13.5cm
- ◆ Weight: blue, green, yellow etc
- ◆ Packing: 100pcs/bag, 1000pcs/ctn



3.BOUFFANT CAP (Code No.:SUBFC003)

- ◆ Material: PP non woven
- ◆ Weight: 10gsm, 12gsm, 15gsm etc
- ◆ Size: 18", 19", 20", 21" etc
- ◆ Color: white, blue, green, yellow etc
- ◆ Packing: 100pcs/bag, 1000pcs/ctn



5.SHOWER CAP (Code No.:SUSWC005)

- ◆ Material: PE material
- ◆ Weight: 3.0g, 3.5g, 4.0g etc
- ◆ Size: 19", 20", 21" etc
- ◆ Color: transparent
- ◆ Packing: 100pcs/bag, 1000pcs/ctn



2.CLIP CAP (Code No.:SUCPC002)

- ◆ Material: PP non-woven
- ◆ Weight: 10gsm, 12gsm, 15gsm etc
- ◆ Type: double or single elastic
- ◆ Size: 18", 19", 20", 21" etc
- ◆ Color: white, blue, green, yellow etc
- ◆ Packing: 100pcs/bag, 1000pcs/ctn



4.SPACE CAP (Code No.:SUSPC004)

- ◆ Material: PP non-woven / SMS
- ◆ Weight: 20gsm, 25gsm, 30gsm etc
- ◆ Size: 36*42cm etc
- ◆ Color: white, blue, green etc
- ◆ Packing: 100pcs/bag, 1000pcs/ctn



6.NET CAP (Code No.:SUNTC006)

- ◆ Material: Nylon material
- ◆ Size: 18", 19", 20", 21", 24", 28" etc
- ◆ Color: white, blue, green, yellow etc
- ◆ Packing: 100pcs/bag, 1000pcs/ctn



Shoe Cover

- ◆ Material: PP non-woven, PE, CPE
- ◆ Weight: 20gsm, 25gsm, 30gsm, 35gsm etc
- ◆ Type: Normal or anti-skid treatment, hand or machine
- ◆ Color: White, green, yellow, red, blue etc
- ◆ Size: 15x40cm, 17x40cm, 17x41cm etc

Code no	Size	Packing	Carton size
SUSC251540	25gsm, 15x40cm	100pcs/bag, 2000pcs/ctn	60x31x30cm
SUSC301540	30gsm, 15x40cm	100pcs/bag, 2000pcs/ctn	60x31x30cm



Anti-skid treatment

PE

PP non woven

CPE

Coverall

- ◆ Code No.: SUCVA007
- ◆ Material: PP/SMS/SF/MP
- ◆ Size: S, M, L, XL, XXL, XXXL
- ◆ Weight: 35gsm, 40gsm, 50gsm, 60gsm etc
- ◆ With or without hood/shoe cover
- ◆ Elastic cuff or knitted cuff
- ◆ Color: White, blue, yellow etc
- ◆ Packing: 1pc/pouch, 50pcs/ctn (Sterile)
5pcs/bag, 100pcs/ctn (Non sterile)



Surgical Gown

- ◆ Code No.: SUSGG008
- ◆ Material: PP/SMS/Reinforced
- ◆ Type: Standard or Reinforced
- ◆ Weight: 30gsm, 35gsm, 40gsm, 45gsm, 50gsm, 60gsm etc
- ◆ Cuff: Elastic cuff or knitted cuff
- ◆ Size: 115x137cm/120x140cm/125x150cm/130x160cm etc
- ◆ Color: Blue, light blue, green, yellow, dark blue etc
- ◆ Packing: 10pcs/bag, 10bags/ctn (Non-sterile)
1pc/pouch, 50pcs/ctn (Sterile)



Scrub Suit

- ◆ Code No.: SUSBS009
- ◆ Material: PP non-woven or SMS
- ◆ Weight: 30gsm, 35gsm, 40gsm, 45gsm etc
- ◆ Size: S, M, L, XL, XXL etc
- ◆ Short coat and trousers
- ◆ Color: Dark blue, light blue etc
- ◆ Packing: 1set/pack, 50sets/ctn



Patient Gown

- ◆ Code No.: SUPTG010
- ◆ Material: PP non-woven/SMS
- ◆ Weight: 14gsm, 18gsm, 25gsm, 30gsm etc
- ◆ Color: White, green, blue, yellow etc
- ◆ Size: S, M, L, XL, XXL etc
- ◆ Long sleeve, short sleeve, without sleeve
- ◆ Packing: 10pcs/bag, 10bags/ctn



Isolation Gown

- ◆ Code No.: SUILG011
- ◆ Material: PP/PP+PE film/SMS/SF
- ◆ Weight: 15gsm, 20gsm, 25gsm, 30gsm, 35gsm etc
- ◆ Size: S, M, L, XL, XXL, XXXL etc
- ◆ Elastic cuff and knitted cuff
- ◆ Color: White, green, blue, yellow etc
- ◆ Packing: 10pcs/bag, 10bags/ctn



CPE Clean Gown

- ◆ Code No.: SUCCG012
- ◆ Material: 100% polyethylene
- ◆ Style: apron, long sleeves, back empty, thumbs up/elastic wrists,
2 ties at the waist
- ◆ Size: S, M, L, XL, XXL etc
- ◆ Color: white, blue, green etc
- ◆ Weight: 50g/pc, 40g/pc, 35g/pc etc
- ◆ Packing: 1pc/bag, 20pcs/medium bag, 100pcs/ctn



AAMI Surgical Gown

- ◆ Code No.: SUASG013
- ◆ Material: AAMI SMS
- ◆ Level: AAMI1, AAMI 2, AAMI 3, AAMI 4
- ◆ Weight: 45gsm, 50gsm, 55gsm, 60gsm etc
- ◆ Size: S, M, L, XL, XXL, XXXL etc
- ◆ Elastic cuff and knitted cuff
- ◆ Color: blue, light blue etc
- ◆ Packing: 1pc/blister, 25pcs/ctn (Sterile)



Bed Cover

- ◆ Code No.: SUBDC014
- ◆ Material: PP,CPE,SMS,Laminated PE film
- ◆ Weight: 30g,35g,40g etc
- ◆ Color: white, blue etc
- ◆ Size: 102x190cm,140x240cm
- ◆ Packing: 10pcs/bag, 100pcs/ctn



Bed Sheet

- ◆ Code No.: SUBDS016
- ◆ Material: SPP/PP+PE/SMS
- ◆ Weight: 30g,35g,40g,45g etc
- ◆ Color: white, green, blue, yellow etc
- ◆ Size: 170x230cm, 120x220cm, 100x180cm etc
- ◆ Packing: 10pcs/bag, 100pcs/ctn(Non-sterile)
1pc/sterile bag, 50pcs/ctn(Sterile)



Sleeves Cover

- ◆ Code no.: SUSVC018
- ◆ Material: LDPE, HDPE
- ◆ Size: 20x40cm,20x46cm,22x46cm
- ◆ Thickness: 0.02mm, 0.025mm, 0.03mm etc
- ◆ Weight: 2.0g, 2.5g, 3.0g, 3.5g, 4.0g etc
- ◆ Color: transparent, white, blue, yellow etc
- ◆ Packing: 10pcs/roll, 10rolls/bag, 10bags/ctn




Pillow Cover

Code No	SUPLC015
Material	PP non woven
Color	white, blue etc
Size	60x60cm+10cm flap or customized
Style	with elastic ends / square ends or plain
Application	hotel, hospital, beauty salon, household etc
Packing	10pcs/bag, 100pcs/ctn



PE Apron

- ◆ Code No.: SUPEA017
- ◆ Material: LDPE, HDPE
- ◆ Size: 24"x42", 28"x46", 31.5"x49" etc
- ◆ Type: smooth or embossed surface
- ◆ Thickness: 0.022mm, 0.016mm, 0.02mm etc
- ◆ Weight: 10g, 14g, 16g, 19g, 20g etc
- ◆ Color: white, green, yellow, blue, red etc
- ◆ Packing: 100pcs/bag, 10bags/ctn



Cleaning Wipe

- ◆ Code no.: SUCLW019
- ◆ Material: viscose & polyester
- ◆ Style: dyed, wave, grid etc
- ◆ Tech: cross lapped & parallel laying
- ◆ Type: sheet, 1/4 folded, perforated roll
- ◆ Size: 35x60cm, 40cmx50cm, 38x40cm etc
- ◆ Weight: 40g - 100g as per customer request
- ◆ Color: red, blue, yellow, green, pink, white etc
- ◆ Usage: Cleaning wipes for kitchen, vehicle, computer etc



Bed Kit

- ◆ Code no.: SUBDK020
- ◆ Material: PP or SMS
- ◆ Color: white, blue, green etc
- ◆ Size: Pillow Case: 50x70cm
Bed Sheet: 200x130cm
Bed Cover: 240x145cm
- ◆ Packing: 1set/bag, 50sets/ctn
- ◆ Application: Ideal product for hospital, clinic, hotel etc



Mattress Protector

- ◆ Code no.: SUMTP021
- ◆ Material: 70% bamboo fiber 30% polyester
- ◆ Type: Plain white, zippered style
- ◆ Size: Queen or customized
- ◆ Packing: 1pc/PVC bag + Insert Card + Cardboard, 10pcs/ctn
- ◆ Feature: Waterproof, dustmite proof, bed bug proof, breathable, machine wash, tumble dry.



Anti-fog face shield

- ◆ Code no.: SUAFS022
- ◆ Material: Anti fog pet film
- ◆ Foam: Soft PU foam, latex free
- ◆ Elastic strap: White elastic strap, latex free
- ◆ Size: 33x22cm(+/-0.01mm) Color: Transparent
- ◆ Thickness: 0.18mm(+/-0.01mm)
- ◆ Function: Anti flog, anti glare, prevent liquid splashing



Protective Head Hood

- ◆ Code no.: SUPHH023
- ◆ Material: 50g SMS + anti-fog PET eye shield
- ◆ Color: blue or white + transparent shield
- ◆ Full protection for eyes, nose, mouth, ears might exposed to the radiation contaminate




Body bag (set)

- ◆ Code No.: SUBBS024
- ◆ Material: PVC or Nylon etc
- ◆ Thickness: 150micron, 0.35mm etc
- ◆ Size: 216x90cm, 230x90cm 220x80cm etc
- ◆ Color: white, black, green etc
- ◆ Center zipper, with 2/4/6 handles



Drape

Code No	SUDRP001
Material	27gsm blue or green film + 28gsm blue or green viscose
Size	35x50cm, 50x50cm, 50x75cm, 75x90cm etc
Color	Blue or Green
Packing	1pc/pouch, 50pcs/ctn(Sterile)



Laparotomy Pack

Code No.:SULTP002

Accessories	Material	Size	Quantity
Instrument cover	55g film + 28g PP	140x190cm	1pc
Standard surgical gowns	35gsm SMS	XL: 130x150cm	3pcs
Hand towels	Flat pattern	30x40cm	3pcs
Plain sheets	35gsm SMS	140x160cm	2pcs
Utility drapes with adhesive	35gsm SMS	40x60cm	4pcs
Laparotomy drape horizontal	35gsm SMS	190x240cm	1pc
Mayo cover	35gsm SMS	58x138cm	1pc



General Pack

Code No.:SUGNP003

Accessories	Material	Size	Quantity
Wrapping	Blue, 35gsm SMMS	100x100cm	1pc
Table cover	55g PE+30g Hydrophilic PP	160x190cm	1pc
Hand towels	60g white spunlace	30x40cm	6pcs
Standard surgical gown	Blue, 35gsm SMMS	L/120x150cm	1pc
Reinforced surgical gowns	Blue, 45gsm SMMS	XL/130x155cm	2pcs
Drape sheets	Blue, 40gsm SMMS	40x60cm	4pcs
Suture bag	80g paper	15x30cm	1pc
Mayo stand cover	Blue, 43g PE	80x145cm	1pc
Side drapes	Blue, 40g SMMS	120x200cm	2pcs
Head drape	Blue, 40g SMMS	160x240cm	1pc
Foot drape	Blue, 40g SMMS	190x200cm	1pc



Delivery Pack

Code No.:SUDLP004

Accessories	Material	Size	Quantity
Side drape with adhesive tape	Blue, 40g, SMS	75*150cm	1pc
Baby drape	White, 60g, spunlace	75*75cm	1pc
Table cover	55g PE film+30g PP	100*150cm	1pc
Drape	Blue, 40g, SMS	75*100cm	1pc
Leg cover	Blue, 40g, SMS	60*120cm	2pcs
Reinforced surgical gowns	Blue, 40g, SMS	XL/130*150cm	2pcs
Umbilical clamp	Blue or white	/	1pc
Hand towels	White, 60g, spunlace	40*40cm	2pcs



Disposable Underpad

- ◆ Material: Hot rolling fabric+Wood paper+PE
- ◆ Weight: 30g,40g,45g,50g,60g,70g,75g,80g,90g etc
- ◆ Size: 40x60cm, 60x90cm, 80x180cm, 100x120cm etc
- ◆ Color: Blue, pink, white, green etc, with or without SAP
- ◆ Application: For hospital use and homecare



Code no	Size	Packing	Carton size
SUDUP01-001	30x45cm-20g	20pcs/bag,20bags/ctn	42x25x42cm
SUDUP01-002	30x45cm-30g	20pcs/bag,20bags/ctn	42x25x52cm
SUDUP01-003	45x45cm-20g	20pcs/bag,20bags/ctn	46x31x32cm
SUDUP01-004	45x45cm-30g	20pcs/bag,20bags/ctn	46x31x45cm
SUDUP01-005	40x60cm-40g	15pcs/bag,20bags/ctn	42x30x48cm
SUDUP01-006	40x60cm-50g	15pcs/bag,20bags/ctn	42x36x48cm
SUDUP01-007	60x60cm-40g	15pcs/bag,10bags/ctn	50x22x28cm
SUDUP01-008	60x60cm-50g	15pcs/bag,10bags/ctn	50x30x32cm
SUDUP01-009	60x60cm-60g	15pcs/bag,10bags/ctn	52x36x40cm
SUDUP01-010	60x90cm-50g	15pcs/bag,10bags/ctn	48x30x40cm
SUDUP01-011	60x90cm-60g	15pcs/bag,10bags/ctn	52x32x46cm
SUDUP01-012	60x90cm-70g	15pcs/bag,10bags/ctn	58x36x48cm
SUDUP01-013	60x90cm-80g	15pcs/bag,10bags/ctn	62x40x52cm
SUDUP01-014	80x180cm-100g	10pcs/bag,10bags/ctn	57x30x46cm
SUDUP01-015	100x120cm-90g	10pcs/bag,10bags/ctn	52x30x42cm



Diaper

BABY DIAPER (Code No.:SUDIPB01)

Code no	Model	Length*Width(mm)	SAP Weight	Weight/pc	Packing
SUDIPB01-S	S	380*280	4g	23g±1g	52pcs/bag, 4bags/ctn
SUDIPB01-M	M	450*320	5g	28g±1g	52pcs/bag, 4bags/ctn
SUDIPB01-L	L	500*320	6g	33g±1g	52pcs/bag, 4bags/ctn
SUDIPB01-XL	XL	540*320	7g	38g±1g	52pcs/bag, 4bags/ctn



ADULT DIAPER (Code No.:SUDIPA02)

standard type: anti-leakage PE film, leg elastics, left/right tapes, frontal tape, leg cuffs

Code no	Model	Length*Width(mm)	SAP Weight	Weight/pc	Packing
SUDIPA02-M	M	800*650	7.5g	85g	10pcs/bag, 10bags/ctn
SUDIPA02-L	L	900*750	9g	95g	10pcs/bag, 10bags/ctn
SUDIPA02-XL	XL	980*800	10g	105g	10pcs/bag, 10bags/ctn

standard type: anti-leakage PE film, leg elastics, left/right tapes, frontal tape, leg cuffs, wetness indicator

Code no	Model	Length*Width(mm)	SAP Weight	Weight/pc	Packing
SUDIPA02-M-WI	M	800*650	7.5g	85g	10pcs/bag, 10bags/ctn
SUDIPA02-L-WI	L	900*750	9g	95g	10pcs/bag, 10bags/ctn
SUDIPA02-XL-WI	XL	980*800	10g	105g	10pcs/bag, 10bags/ctn



ABD Pad

- ◆ Disposable ABD pad is constructed of highly absorbent rayon and cellulose inner layers, the non-woven facing and highly absorbent cotton filler with rayon fibers diffuse liquid throughout the pad.
- ◆ All edges are sealed to prevent linting.
- ◆ Multiple dressings of this soft, patient-comfortable pad may be used for additional protection against strike trough.
- ◆ Available in wide range of sizes 5.5"x9", 8"x10" etc
- ◆ Both sterile and non-sterile.

NON STERILE ABD PAD, 25PCS/PACK

Code no	Size	Packing	Carton size
SUAP001-1010N	10x10cm	1250pcs/ctn	52x42x42cm
SUAP001-1020N	10x20cm	1250pcs/ctn	52x42x42cm
SUAP001-1520N	15x20cm	375pcs/ctn	34x32x26cm
SUAP001-2020N	20x20cm	375pcs/ctn	44x32x26cm
SUAP001-2040N	20x40cm	375pcs/ctn	44x32x52cm

STERILE ABD PAD, 1PC/POUCH

Code no	Size	Packing	Carton size
SUAP001-1010S	10x10cm	800pouches/ctn	34x30x48cm
SUAP001-1020S	10x20cm	800pouches/ctn	54x50x51cm
SUAP001-1520S	15x20cm	400pcs/ctn	48x30x32cm
SUAP001-2020S	20x20cm	400pcs/ctn	52x30x37cm
SUAP001-2040S	20x40cm	200pcs/ctn	52x30x37cm



Examination Paper Roll

- ◆ Code No.: SUEPR001
- ◆ Material: 1ply paper+1ply film or 2ply paper
- ◆ Weight: 10gsm-35gsm etc
- ◆ Width: 50cm, 60cm, 70cm, 100cm etc
- ◆ Length: 50m, 100m, 150m, 200m etc
- ◆ Pre-cut: 50cm, 60cm etc
- ◆ Usually white, blue or yellow color
- ◆ With perforated lines, easy for tearing



Medical Crepe Paper

- ◆ Code No.: SUMCP002
- ◆ Material: All pure wood pulp
- ◆ Full color: Blue, white, green etc
- ◆ Sterilization method: high temperature
- ◆ High pressure, plasma
- Size: 30x30cm, 60x60cm, 50x50cm, 90x90cm, 100x100cm etc

Code no	Size	Packing	Carton size
SUMCP002-100100	100x100cm	250pcs/ctn	103x39x12cm
SUMCP002-120120	120x120cm	200pcs/ctn	123x45x12cm
SUMCP002-120180	120x180cm	200pcs/ctn	123x92x16cm
SUMCP002-3030	30x30cm	1000pcs/ctn	35x33x15cm
SUMCP002-6060	60x60cm	500pcs/ctn	63x35x15cm
SUMCP002-9090	90x90cm	250pcs/ctn	93x35x12cm
SUMCP002-7575	75x75cm	500pcs/ctn	77x35x10cm
SUMCP002-4040	40x40cm	1000pcs/ctn	42x33x15cm



Latex Surgical Gloves

- ◆ Code no.: SULSG001
- ◆ 100% latex
- ◆ 6.5# 7# 7.5# 8# 8.5#(7 5# 17g/pair)
- ◆ Powder and powder free
- ◆ 1pair/pouch, 50pairs/box, 10boxes/ctn



Latex Examination Gloves

- ◆ Code no.: SULEG002
- ◆ 100% latex
- ◆ S:5g/M:5.5g/L:6.0g/XL:6.0g
- ◆ Powder and powder free
- ◆ 100pcs/box, 10boxes/ctn



Nitrile Gloves

- ◆ Code no.: SUNG003
- ◆ 100% nitrile
- ◆ XL/L/M/S
- ◆ Purple/Dark Blue/Light Blue ect
- ◆ Powder free, durable and stretch
- ◆ 100pcs/box, 1000pcs/ctn



Vinyl Gloves

- ◆ Code no.: SUVG004
- ◆ 100% vinyl
- ◆ S/M/L
- ◆ S 4.0g/M 4.5g/L 5.0g
- ◆ Powder and powder free
- ◆ 100pcs/box, 10boxes/ctn



PE Gloves

- ◆ Code no.: SUPG005
- ◆ HDPE/LDPE
- ◆ 5 fingers/2 fingers
- ◆ 0.7g, 1g, 1.5g etc
- ◆ S(27*24cm)/M(28*26cm)/L(29*26cm)
- ◆ 100pcs/bag, 100bags/ctn



Surgical Suture


Code No.:SUSST001

Suture classification	Synthetic absorbable; Nature absorbable; Non absorbable
Suture type	Polyglycolic acid; Polyglactine; Polyglactine rapid; Polydioxanone; Chromic catgut; Plain catgut; Silk; Nylon; Polyester; Polypropylene
Suture Length	45cm, 75cm, 90cm, 150cm or customized
Thread NO	0, 1, 2, 2/0, 3/0, 4/0, 5/0, 6/0, 7/0, 8/0
Needle type	1/2, 3/8, 1/4, 5/8 circle; J shape; Compound curve; 1/2 curved; straight
Needle Length	6mm, 8mm, 10mm, 12mm, 13mm, 15mm, 16mm, 18mm, 20mm, 22mm, 25mm, 30mm, 38mm, 40mm or customized
Needle Shape	Round Bodied (Taper Point), Curved Cutting (Regular Cutting), Reverse Cutting, Tapercut, Micro-point Reverse Cutting, Micro-point Spatula, Curved, Blunt Taper Point, Taper Point Fishhook



Symbol	Point Type of Needle
TAPER POINT NEEDLE	
BLUNT TAPER POINT NEEDLE	
TAPERCUT NEEDLE	
CUTTING NEEDLE	
REVERSE CUTTING NEEDLE	
MICRO-PPINT SPATULA NEEDLE	

Shapes



End Type

ROILLED-END	
ROILLED-END	
REGULAR RYE	
SPRING EYE	
SPRING DOUBLE EYE	



Item	Spe.	Packing	Carton size
Polyglycolic acid, length 75cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0, 6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Polyglactine, length 75cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0, 6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Polyglactine rapid, length 75cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0, 6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Polydioxanone, length 75cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0, 6/0, 7/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Chromic catgut, length 75cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0, 3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Plain catgut, length 75cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0, 3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Silk, length 75cm	0, 1, 2, 2/0, 3/0, 4/0, 5/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	7/0, 8/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Nylon, length 75cm	0, 1, 2, 2/0, 3/0, 4/0, 5/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	7/0, 8/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Polyester, length 75cm	0, 1, 2, 2/0, 3/0, 4/0, 5/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	7/0, 8/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Polypropylene, length 75cm	0, 1, 2, 2/0, 3/0, 4/0, 5/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	6/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	7/0, 8/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
Polypglycolic acid, length 90cm	0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	1	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	2/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	3/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm
	4/0, 5/0, 6/0, 7/0, 8/0	12pcs/box, 100boxes/ctn	52 * 28 * 33cm

Disposable Syringe

- ◆ Code no.: SUDSS002
- ◆ Material: medical grade PVC(latex or latex free)
- ◆ Size: 1cc,2cc,3cc,5cc,10cc,20cc,30cc,50cc,60cc,100cc
- ◆ With or without needle, luer lock or luer slip
- ◆ Needle type: mounted and unmounted
- ◆ Needle size: 16G-30G
- ◆ Needle Hub: Made of medical grade PP
- ◆ Needle protector: Made of medical grade transparent PP
- ◆ Silicone oil: ensure the sharpness and smooth of needle
- ◆ Adhesive: Epoxy resin is used to fix the needle and hub.
- ◆ Type: normal type; auto disable type; safety type
- ◆ Packing: PE packing, blister packing.



Insulin syringe



Auto-destruct safety syringe



Tipped syringe

PE UNIT PACKING WITH MOUNTED NEEDLE

Specs	Pcs/Box	Boxes/Carton	Pcs/Carton	Carton size	Cartons/20gp	Cartons/40hq
1ml	100	36	3600	58x43x44.5cm	260	620
2ml	100	27	2700	58x43x44.5cm	260	620
3ml	100	27	2700	58x43x44.5cm	260	620
5ml	100	24	2400	58x43x44.5cm	260	620
10ml	100	16	1600	67x44x42cm	240	560
20ml	50	18	900	57x50x43cm	240	560

BLISTER UNIT PACKING WITH UNMOUNTED NEEDLE

Specs	Pcs/Box	Boxes/Carton	Pcs/Carton	Carton size	Cartons/20gp	Cartons/40hq
1ml	100	32	3200	61.5x38x48.5cm	260	600
2ml	100	30	3000	61.5x38x48.5cm	260	600
3ml	100	30	3000	61.5x38x48.5cm	260	600
5ml	100	24	2400	61.5x38x48.5cm	260	600
10ml	100	16	1600	61.5x38x48.5cm	260	600
20ml	50	16	800	61.5x38x48.5cm	260	600
30ml	40	16	640	61.5x38x48.5cm	260	600
50ml	25	16	400	61.5x38x48.5cm	260	600



PE pack



Blister pack

Infusion Set

- ◆ Code no.: SUIFS003
- ◆ Medical grade PVC, EO sterile
- ◆ Luer lock or luer slip
- ◆ With or without needle
- ◆ Drip chamber with or without filter
- ◆ Spike with or without air vent
- ◆ Y type, flash bulb injection, ordinary type, burette type
- ◆ Tube length: 1.35m, 1.55m, 2m or customized
- ◆ PE or Blister packing



Burette Infusion Set

- ◆ Code no.: SUBIS005
- ◆ Made of medical grade PVC
- ◆ Child, infant
- ◆ 100ml, 150ml
- ◆ 20 drops or 60drops
- ◆ PE packing or Blister packing



Three Way Stopcock

- ◆ Code no.: SUTWS007
- ◆ Material: PC, PE
- ◆ With or without extension tube
- ◆ Tube length: 10cm, 15cm, 20cm, 25cm, 50cm etc
- ◆ PE Packing or blister packing
- ◆ Also available with Y injection port(latex and latex free)
- ◆ Packing: 50pcs/box 1000pcs/ctn
- ◆ Carton size: 57.5x38x38cm



Blood Transfusion Set

- ◆ Code no.: SUBTS004
- ◆ Medical grade PVC, EO sterile
- ◆ Luer lock or luer slip
- ◆ With or without needle
- ◆ Flow rate: 20drops/15drops per min
- ◆ Drip chamber with filter to prevent passage of any clot
- ◆ Flow regulator: PE or ABS
- ◆ Tube length: 1.5m/1.8m/2.0m
- ◆ PE or Blister packing



Flow Controuers

- ◆ Code no.: SUFCT006
- ◆ Material: Medical grade PV pr polystyrene
- ◆ Main accessories: Vented spike, dripping Chamber, Solution Filter, flow regulator, latex tube, Luer Lock connector
- ◆ Approximately 15 drops/ml, 20 drops/ml



Hypodermic Needle

- ◆ Medical grade PP or PVC
- ◆ Luer Slip or Luer Lock
- ◆ Size: 14G, 16G, 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G, 30G
- ◆ Transparent barrel and plunger, Latex or latex-free
- ◆ EO sterile and non-sterile
- ◆ PE pack or Blister pack

Code no.	Specification	Colour code
SUHN30G	30Gx1/2"	Yellow
SUHN29G	29Gx1/2"	Red
SUHN27G	27Gx1/2", 27Gx5/8"	Gray
SUHN26G	26Gx1/2", 26Gx5/8", 26Gx3/4"	Brown
SUHN25G	25Gx5/8", 25Gx3/4", 25Gx1"	Orange
SUHN24G	24Gx3/4", 24Gx1"	Purple
SUHN23G	23Gx1", 23Gx1 1/4", 23Gx1 1/2"	Deep blue
SUHN22G	22Gx1", 22Gx1 1/4", 22Gx1 1/2"	Black
SUHN21G	21Gx1 1/4", 21Gx1 1/2"	Deep green
SUHN20G	20Gx1 1/4", 20Gx1 1/2"	Yellow
SUHN19G	19Gx1 1/4", 19Gx1 1/2"	Crean
SUHN18G	18Gx1 1/4", 18Gx1 1/2"	Pink
SUHN16G	16Gx1 1/4", 16Gx1 1/2"	White



Blood Collection Needle

- ◆ Medical grade PVC and stainless steel
- ◆ Pen type, SVS, Holder
- ◆ 18G, 20G, 21G, 22G, 23G
- ◆ Cap color: Yellow, green, black, pink, purple
- ◆ Sterilized by EO gas, non-toxic and non-pyrogenic
- ◆ Components: Needle, cap and sleeve
- ◆ Packing: 100pcs/box, 50boxes/carton

Code no.	Gauge	Colour
SUBCN001	23Gx1"	Deep blue
SUBCN002	22Gx1"	Black
SUBCN003	22Gx1 1/2"	Black
SUBCN004	21Gx1"	Deep green
SUBCN005	21Gx1 1/2"	Deep green
SUBCN006	20Gx1"	Yellow
SUBCN007	20Gx1 1/2"	Yellow
SUBCN008	18Gx1 1/2"	Pink
SUBCN009	16Gx1 1/2"	White



IV Cannula

- ◆ Medical grade PVC+stainless steel
- ◆ Pen type, with wings, with injection port, Y type
- ◆ 14G, 16G, 18G, 20G, 22G, 24G, 26G
- ◆ Orange, Grey, Green, Pink, Blue, Yellow, Violet
- ◆ Sterilized by EO gas, non-toxic, non-pyrogenic
- ◆ Blister pack

I.V. CANNULA: PEN-LIKE, WITH WINGS, WITH INJECTION PORT TYPE

Code no.	Gauge	Product Specification	Flow	Color Code
SUIC14GP	14G	2.00*(2.1mmx50mm)	300mL/min	Orange
SUIC16GP	16G	2.00*(1.7mmx50mm)	200mL/min	Medium Gray
SUIC18GP	18G	1.77*(1.3mmx45mm)	90mL/min	Dark Green
SUIC20GP	20G	1.30*(1.1mmx33mm)	61mL/min	Pink
SUIC22GP	22G	1.00*(0.9mmx25mm)	36mL/min	Dark Blue
SUIC24GP	24G	0.75*(0.7mmx19mm)	18mL/min	Yellow
SUIC26GP	26G	0.63*(0.6mmx16mm)	12mL/min	Purple

I.V. CANNULA: Y TYPE

Code no.	Gauge	Product Specification	Flow	Color Code
SUIC18GY	18G	1.16*(1.3mmx30mm)	80mL/min	Dark Green
SUIC20GY	20G	1.16*(1.1mmx30mm)	50mL/min	Pink
SUIC22GY	22G	1.00*(0.9mmx25mm)	33mL/min	Dark Blue
SUIC24GY	24G	0.75*(0.7mmx19mm)	24mL/min	Yellow
SUIC26GY	26G	0.63*(0.6mmx16mm)	12mL/min	Purple



Injection port type



Butterfly type



Pen like type



Y type

Scalp Vein Set

- ◆ Medical grade PVC(Butterfly Needle)
- ◆ Luer slip or luer lock
- ◆ 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
- ◆ Tube length: 30cm standard
- ◆ EO sterile and non-sterile
- ◆ PE, blister, paper-plastic bag

Code no.	Specification	Color Code
SUSVS27G	27Gx1/2", 27Gx3/4"	Grey
SUSVS26G	26Gx5/8", 26Gx3/4"	Brown
SUSVS25G	25Gx3/4"	Orange
SUSVS24G	24Gx3/4"	Purple
SUSVS23G	23Gx3/4", 23Gx1", 23Gx1 1/8"	Blue
SUSVS22G	22Gx3/4", 22Gx1", 22Gx1 1/8"	Black
SUSVS21G	21Gx3/4", 21Gx1 1/8"	Green
SUSVS20G	20Gx3/4", 20Gx1 1/8"	Yellow
SUSVS19G	19Gx3/4", 19Gx1 1/8"	Beige
SUSVS18G	18Gx3/4", 18Gx1 1/8"	Pink



Needle For Anesthesia

1.Spinal Needle

- ◆ Medical grade PVC+stainless steel
- ◆ Type: quinke tip or pencil-point
- ◆ 16G, 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G
- ◆ Length:45mm, 60mm, 75mm, 90mm, 120mm, 150mm etc
- ◆ With or without intr
- ◆ Sterilized by EO gas, sterile, non-toxic, non-pyrogenic

QUINCKE TIP

Code no.:	Size	Introducer	Color	
SUANQT18	18Gx90mm	none	Pink	Pink
SUANQT19	19Gx90mm	none	Milk white	Cream
SUANQT20	20Gx90mm	none	Yellow	Yellow
SUANQT21	21Gx90mm	none	Green	Green
SUANQT22	22Gx90mm	none	Black	Black
SUANQT23	23Gx90mm	19Gx38mm	Blue	Blue
SUANQT24	24Gx90mm	20Gx38mm	Purple	Purple
SUANQT25	25Gx90mm	20Gx38mm	Orange	Orange
SUANQT26	26Gx90mm	20Gx38mm	Brown	Brown
SUANQT27	27Gx90mm	22Gx38mm	Medium gray	Grey

Length can be customized(38-220mm)

PENCIL POINT

Code no.:	Size	Introducer	Color	
SUANPP22	22Gx90mm	none	Black	Black
SUANPP23	23Gx90mm	19Gx38mm	Blue	Blue
SUANPP24	24Gx90mm	20Gx38mm	Purple	Purple
SUANPP25	25Gx90mm	20Gx38mm	Orange	Orange
SUANPP26	26Gx90mm	20Gx38mm	Brown	Brown
SUANPP27	27Gx90mm	22Gx38mm	Medium gray	Grey

Length can be customized(38-220mm)



Spinal Needle



Epidural Needle



Dental Anesthesia Needle

2.Epidural Needle

PENCIL POINT

Code no.:	Size	Size	Length(mm)
SUENPP16	16Gx90mm	16G x 3 1/2	80/90mm
SUENPP17	17Gx90mm	16G x 3 1/2	80/90mm
SUENPP18	18Gx90mm	16G x 3 1/2	80/90mm

80mm and 90mm two length available

3.Dental Needle

- ◆ British and metric system
- ◆ Length:8mm, 12mm, 14mm, 16mm, 21mm, 25mm, 35mm, 38mm
- ◆ Size: 27G, 30G

Code no.	Specification		Length(mm)
	Gauge	Color code	
SUDN27G	27G	Green	8mm, 12mm, 14mm, 16mm, 21mm, 25mm, 35mm, 38mm
SUDN30G	30G	Yellow	8mm, 12mm, 14mm, 16mm, 21mm, 25mm, 35mm, 38mm

Vacuum Blood Collection Tube

Code No.:SUYBCT001

Item	Size	Volume	Tube material	Additive	Colour	Main Intended Use	Packing
Plain Tube	13x75mm	1-5ml	Glass or PET	No additive	Red	Determinations in serum for clinical chemistry, Microbiological serology, immunology	13x75mm 100pcs/pallet, 1800pcs/CTN Carton size:51x31x29mm
	13x100mm	5-7ml					
	16x100mm	8-10ml					
Pro-Coagulation Activator Tube	13x75mm	1-5ml	Glass or PET	Clot activator	Red	Determinations in serum for clinical chemistry, Microbiological serology, immunology	13x75mm 100pcs/pallet, 1800pcs/CTN Carton size:51x31x29mm
	13x100mm	5-7ml					
	16x100mm	8-10ml					
Gel&Clot Activator Tube	13x75mm	1-4ml	Glass or PET	Gel&Clot Activator	Yellow	Determinations in serum for clinical chemistry, Microbiological serology, immunology	13x75mm 100pcs/pallet, 1800pcs/CTN Carton size:51x31x29mm
	13x100mm	4-6ml					
	16x100mm	7-8.5ml					
K:EDTA Tube	13x75mm	1-5ml	Glass or PET	K:EDTA	Purple	Determinations in EDTA whole blood for hematology	13x100mm 100pcs/pallet, 1800pcs/CTN Carton size:51x36x36mm
	13x100mm	5-7ml					
K:EDTA Tube	13x75mm	1-5ml	Glass or PET	K:EDTA	Purple	Determinations in EDTA whole blood for hematology	13x100mm 100pcs/pallet, 1800pcs/CTN Carton size:51x36x36mm
	13x100mm	5-7ml					
PT Tube	13x75mm	1.8-4.5ml	Glass	Sodium Citrate	Blue	Determinations in citrated plasma for coagulation testing	16x100mm 100pcs/pallet, 1800pcs/CTN
		1.8-3.6ml					
ESR Tube	13x75mm	1.6-2.4ml	Glass or PET	Sodium Citrate	Black	Test for the blood cell sedimentation rate	16x100mm 100pcs/pallet, 1800pcs/CTN
	8x120mm	1.6ml					
Heparin Tube	13x75mm	1-5ml	Glass or PET	Sodium Heparin or Lithium Heparin	Green	Determinations in heparinized plasma for clinical chemistry	16x100mm 100pcs/pallet, 1800pcs/CTN
	13x100mm	5-7ml					
Glucose Tube	13x75mm	1-5ml	Glass or PET	Inhibitor for Glucolysis Sodium fluoride +EDTA Sodium fluoride +Potassium oxalate	Grey	Determinations in stabilized anticoagulated whole blood or plasma for glucost and lactate testing	8x120mm 100pcs/pallet, 1800pcs/CTN Carton size:46x36x42mm
	13x100mm	5-7ml					
Gel&Heparin Tube	13x75mm	1-4ml	Glass or PET	Gel&Sodium or Lithium Heparin	Green	For plasma determination in chemistry	8x120mm 100pcs/pallet, 1800pcs/CTN Carton size:46x36x42mm
	13x100mm	4-6ml					
	16x100mm	7-8.5ml					

Capillary Tube

- ◆ Code No.: SUCT002
- ◆ Specially designed for infants, dying and other patients with blood collection problems.
- ◆ The groove design ensures that the venous blood flows into the tube smoothly and quickly, and easy to mix the additive inside which prevent blood clotting.
- ◆ The toothed bumps are two side. It makes the cap and tube combine tightly. It effectively prevent blood samples leaking.




Capillary Tube

Vacuum Blood Collection Tube

Silicone Foley Catheter

- ◆ Code No.: SUSFC001
- ◆ Made from 100% medical grade silicone
- ◆ Good for long term placement
- ◆ Size: 2-way pediatric; length: 270mm, 8Fr-10Fr, 3/5cc (balloon)
2-way standard; length: 400mm, 12Fr, 14Fr, 5/10cc (balloon)
16Fr-24Fr, 5/10/30cc (balloon)
- ◆ 3-way standard; length: 400mm, 16Fr-26Fr, 30cc balloon
- ◆ Color-coded for visualization of size
- ◆ Length: 310mm (Pediatric); 400mm (Standard)
- ◆ Sterilized by EO Gas, 1pc/pouch
- ◆ Single use only.

Size	Packing	Carton size
2way, F8-F10	10pcs/box, 500pcs/ctn	52.5x41x43cm
2way, F12-F22	10pcs/box, 500pcs/ctn	52.5x41x43cm
2way, F24-F26	10pcs/box, 500pcs/ctn	52.5x41x43cm
3way, F14-F22	10pcs/box, 500pcs/ctn	52.5x41x43cm
3way, F24-F26	10pcs/box, 500pcs/ctn	52.5x41x43cm



Male External Catheter

- ◆ Medical grade latex
- ◆ Size: S, M, L, XL
- ◆ Diameter: 20mm, 25mm, 30mm, 35mm
- ◆ For male's urine drainage use only
- ◆ Sterilized by EO Gas

Code no.	Size	Packing	Carton size
SUMEC003-S	S, 20mm	1pc/blister, 50pcs/box, 20boxes/ctn	55x28x57cm
SUMEC003-M	M, 25mm	1pc/blister, 50pcs/box, 20boxes/ctn	55x28x57cm
SUMEC003-L	L, 30mm	1pc/blister, 50pcs/box, 20boxes/ctn	55x28x57cm
SUMEC003-XL	XL, 35mm	1pc/blister, 50pcs/box, 20boxes/ctn	55x28x57cm



Male Condoms

- ◆ Code no.: SUMCM004
- ◆ Material: 100% latex
- ◆ Width: (49/52/53/55) ±2mm
- ◆ Length: ≥170mm
- ◆ Type: Smooth type, Float type, Thread type, Floating point thread combination type
- ◆ Packing: 12 pcs/box, 12 boxes/medium box, 36 medium boxes/carton



Latex Foley Catheter

- ◆ Code No.: SULFC002
- ◆ Made of natural latex
- ◆ Size: 1 way, 6Fr-24Fr
2-way, pediatric, 6Fr-10Fr, 3-5ml
2-way, standard, 12Fr-20Fr, 5ml-15ml/30ml/cc
2-way, standard, 22Fr-24Fr, 5ml-15ml/30ml/cc
3-way, standard, 16Fr-24Fr, 5ml-15ml/cc 30ml-50ml/cc
- ◆ Length 25-27cm (Children and Female) 40cm (Male)
- ◆ Sterilized by EO Gas, 1pc/pouch
- ◆ Single use only.

Size	Packing	Carton size
2way, F8-F10	10pcs/box, 500pcs/ctn	52.5x41x43cm
2way, F12-F22	10pcs/box, 500pcs/ctn	52.5x41x43cm
2way, F24-F26	10pcs/box, 500pcs/ctn	52.5x41x43cm
3way, F14-F22	10pcs/box, 500pcs/ctn	52.5x41x43cm
3way, F24-F26	10pcs/box, 500pcs/ctn	52.5x41x43cm



Feeding Tube

- ◆ Code No.: SUFDT005
- ◆ Medical grade PVC
- ◆ Transparent or frosted
- ◆ Tube length: 40cm±2cm (Infant), 120cm±2cm (Adult)
- ◆ Size: F4, F5, F6, F8, F10, F12, F14, F16, F18, F20, F22
- ◆ PE packing or blister packing
- ◆ Sterilized by EO GAS



Suction Catheter

- ◆ Code No.: SUSTC007
- ◆ Medical grade PVC
- ◆ Transparent or frosted
- ◆ Tube length: 30cm, 45cm, 60cm etc
- ◆ Size: F5, F8, F10, F12, F14, F16, F18, F20, F22, F24
- ◆ PE packing or blister packing
- ◆ Sterilized by EO GAS



Rectal Tube

- ◆ Made of non-toxic medical-grade PVC, clear and soft
- ◆ Length: 38cm
- ◆ Size: 14-36fr
- ◆ PE packing or blister packing

Code No.	Size(Fr)	Color	Size(Fr)	Color
SURCT009-14	14	Green	26	Pink
SURCT009-16	16	Orange	28	Grass green
SURCT009-18	18	Red	30	Black
SURCT009-20	20	Yellow	32	Dark blue
SURCT009-22	22	Purple	34	Jade green
SURCT009-24	24	Dark blue	36	Light blue

Stomach Tube

- ◆ Code No.: SUSMT006
- ◆ Medical grade PVC
- ◆ Transparent or frosted
- ◆ Tube length: 90cm, 110cm, 125cm etc
- ◆ Size: F6, F8, F10, F12, F14, F16, F18, F20, F22
- ◆ PE packing or blister packing
- ◆ Sterilized by EO GAS



Nelaton Catheter

- ◆ Code No.: SUNTC008
- ◆ Medical grade PVC
- ◆ Transparent or frosted
- ◆ Tube length: 20cm, 30cm, 40cm, 60cm etc
- ◆ Size: F6, F8, F10, F12, F14, F16, F18, F20, F22
- ◆ PE packing or blister packing
- ◆ Sterilized by EO GAS



Nasal Oxygen Cannula

- ◆ Code No.: SUNOC010
- ◆ Made of Medical grade PVC.
- ◆ Suitable for Adult, Pediatric, Infant and New born.
- ◆ Color: Green or white.
- ◆ Technique: injection & deep liquid.
- ◆ Component: with 2m main tube, 550mm.
- ◆ Branch tube, total length is 2.5m, 5mm OD.
- ◆ Packing: 1pc/PE bag, 250pcs/ctn

Oxygen Mask

- ◆ Made of Medical grade PVC.
- ◆ Suitable for Adult, Pediatric and Infant.
- ◆ Size: S/M/L/XL.
- ◆ Color: Green or white.
- ◆ Tube length 2m, 6mm OD

Code No.	Specification	Packing	Carton size
SUOGM011-S	S-New born	1pc/PE bag, 50pcs/ctn	49x28x24cm
SUOGM011-M	M-Child	1pc/PE bag, 50pcs/ctn	49x28x24cm
SUOGM011-L	L-Adult	1pc/PE bag, 50pcs/ctn	49x28x24cm
SUOGM011-XL	XL-Adult	1pc/PE bag, 50pcs/ctn	49x28x24cm

Nebulizer Mask

- ◆ Made of Medical grade PVC.
- ◆ Adjustable nasal clip.
- ◆ Size: S/M/L/XL.
- ◆ Suitable for Adult, Pediatric and Infant.
- ◆ Color: Green or white.
- ◆ Tube length 2m, 6mm OD

Code No.	Specification	Packing	Carton size
SUNBM012-S	S-New born	1pc/PE bag, 50pcs/ctn	49x28x24cm
SUNBM012-M	M-Child	1pc/PE bag, 50pcs/ctn	49x28x24cm
SUNBM012-L	L-Adult	1pc/PE bag, 50pcs/ctn	49x28x24cm
SUNBM012-XL	XL-Adult	1pc/PE bag, 50pcs/ctn	49x28x24cm

PVC Anesthesia Mask

- ◆ Mask for anesthesia, resuscitator and other, applications involving oxygen treatment.
- ◆ Material: Medical PVC materials.
- ◆ Available in 7 sizes with color-coded rings.
- ◆ Single patient use

Code No.	Specification	Packing	Carton size
SUPAM013-0	0#	1pc/PE bag, 350pcs/ctn	50x40x40cm
SUPAM013-1	1#	1pc/PE bag, 250pcs/ctn	50x40x40cm
SUPAM013-2	2#	1pc/PE bag, 170pcs/ctn	50x40x40cm
SUPAM013-3	3#	1pc/PE bag, 350pcs/ctn	50x40x40cm
SUPAM013-4	4#	1pc/PE bag, 125pcs/ctn	50x40x40cm
SUPAM013-5	5#	1pc/PE bag, 80pcs/ctn	50x40x40cm

CPR Mask

- ◆ Code No.: SUCM014
- ◆ Made of medical grade PVC, soft and smooth.
- ◆ Suitable for Adult, Pediatric and Infant
- ◆ Equip with a replaceable one-way valve with biological filter in one housing.



Manual Resuscitator

- ◆ Code No.: SUMRT016
- ◆ Material: PVC or silicone.
- ◆ Accessories: PVC or silicone mask, PVC oxygen tubing, EVA reservoir bag.
- ◆ Size: Adult, Pediatric and Infant.
- ◆ Packing: PE bag, PP box, Paper box.
- ◆ Latex free.



Oxygen Mask With Reservoir Bag

- ◆ Code No.: SUOMRB017
- ◆ Material: Medical grade non-toxic PVC.
- ◆ Component: Mask, oxygen tubing, connector, reservoir bag.
- ◆ Mask size: L/XL(Adult), M(Pediatric), S(Infant).
- ◆ Tube size: With or without 2m anti-crush tube(Customized).
- ◆ Reservoir bag: 1000ml.
- ◆ Color: Green/transparent



Colostomy Bags

- ◆ Code No.: SUCTB018
- ◆ Material: hydrocolloid adhesive, film
- ◆ Liner: non woven
- ◆ Color: transparent and skin
- ◆ Closure: Closed system and open system
- ◆ Type: one-piece, two pieces
- ◆ Packing: 20pcs/box, 40boxes/ctn

Product	Barrier Type	Size
One Piece Colostomy Bag, Closed	Hydrocolloid Adhesive	40mm,50mm,60mm,70mm
One Piece Colostomy Bag, Drainable	Hydrocolloid Adhesive	40mm,50mm,60mm,70mm
One Piece Colostomy Bag, Drainable	Hydrocolloid & Non-woven Border	40mm,50mm,60mm,70mm
Two Pieces Colostomy Bag, Closed		40mm,50mm,60mm,70mm
Two Pieces Colostomy Bag, Drainable		40mm,50mm,60mm,70mm
Two Pieces Colostomy Bag		40mm,50mm,60mm,70mm
Two Pieces Colostomy Bag, skin Barrier	Hydrocolloid Adhesive	40mm,50mm,60mm,70mm
Two Pieces Colostomy Bag, skin Barrier	Hydrocolloid & Non-woven Border	40mm,50mm,60mm,70mm



Yankauer Handle

- ◆ Code no.: SUYKH019
- ◆ Made of K-resin
- ◆ Available with suction connection tube
- ◆ Size: crown tip, with control vent.
crown tip, without control vent.
standard tip, with control vent.
standard tip, without control vent.

Suction Connection Tube

- ◆ Code no.: SUSCT020
- ◆ Made of non-toxic PVC, medical grade.
- ◆ Flexible PVC resist well the negative pressure from sucking flat.
- ◆ Available with three different connectors.
- ◆ Size: 3/16" 1.8m; 3/16" 3.6m; 1/4" 1.8m; 1/4" 3.6m

Breathing Circuit

- ◆ Code no.: SUBTC021
- ◆ Material: medical grade PVC
- ◆ Normal length: 2.0m, 2.4m, can be customized
- ◆ Model: Adult, Child
- ◆ Color: Transparent, green, white
- ◆ Diameter: 22mm, 15mm
- ◆ Length: 1.6m/1.8m/2m/3m
- ◆ Application: Anesthesia operation

Code No.	Specification	Packing	Carton size
SUBTC021-22	22mm	100pcs/ctn	74*39*35.5cm
SUBTC021-15	15mm	100pcs/ctn	74*39*35.5cm

Breathing Exerciser

- ◆ Code no.: SUBTE022
- ◆ Material: Medical PVC, ABS, PP, PE
- ◆ Sterilization: EO sterilization
- ◆ Size: 600cc/900cc/1200cc

Item	Packing	Carton size
Single ball breathing Exerciser	1pc/box, 50pcs/ctn	57*37*47cm
Three ball breathing Exerciser	1pc/box, 50pcs/ctn	57*37*47cm

Mucus Extractor

- ◆ Code no.: SUMET023
- ◆ Material: non-toxic PVC, medical grade
- ◆ Sterilization: EO GAS
- ◆ Feature: Soft and clear

Item	Model Size	Packing (Individual pack)	Carton size
Mucus Extractor	25ml	20pcs/box, 10boxes/ctn	57*37*47cm
Infant Mucus Extractor	F6 - F14	20pcs/box, 10boxes/ctn	57*37*47cm



Umbilical Catheter

- ◆ Code No: SUUBC024
- ◆ Material: medical grade PVC
- ◆ Tube length: 40cm
- ◆ Size: F4, F5, F6, F8, etc
- ◆ PVC tube: transparent, frosted
- ◆ PE packing or blister packing
- ◆ Sterilized by EO Gas



Closed Suction System

- ◆ Code No: SUCSS025
- ◆ Material: 100% medical grade PVC
- ◆ Size: Fr6, Fr8, Fr10, Fr12, Fr14, Fr16
- ◆ Component: marked tube, swivel female 15mm (standard connector) and etc.
- ◆ Color: light green, blue, black white, red etc
- ◆ Blister packing, Sterilized by EO GAS.



Saliva Ejector

- ◆ Code No: SUSET026
- ◆ Material: Transparent PVC tube, copper interior wiring with chrome plated
- ◆ Size: length 150±2mm Diameter 6.8±0.5mm
- ◆ Body color: Green, red, blue, yellow, clear optional as client's request
- ◆ Sterilized: It is sterilized by EO Gas and safety to use.
- ◆ Saliva ejector features a securely fastened, smooth tip for patient comfort.



Chest tube with trocar

- ◆ Code No: SUCTT027
- ◆ Material: medical grade PVC/TPU
- ◆ Size: 7F, 8F, 10F, 12F, 16F, 20F, 24F, 28F, 32F
- ◆ Length: 250mm/280mm/310mm/350mm
- ◆ With trocar or without trocar ◆ Packing: 1pc/pouch, 10pcs/box, 100pcs/ctn
- ◆ Carton size: 46 5x34 5x28cm



Bubble humidifier bottle

- ◆ Code No: SUBHB028
- ◆ Material: PP
- ◆ Color: transparent, white, green lid, blue lid etc.
- ◆ Size: 40ml, 200ml, 250ml, 500ml
- ◆ Packaging: 1pc/PE bag, 50pcs/carton
- ◆ Carton size: 45*45*40cm
- ◆ Sterilized by EO Gas



Aerochamber

- ◆ Code No: SUACB029

Spacer body size	52x131mm (bottom diameter x height)
Volume	175 ml
Mouthpiece	22 .5x15 5mm (top of mouthpiece dimension)
Spacer body weight	40g
Mask weight	Small: 19.7g Medium: 23.7g Large: 43.8g
Warranty	One year for defects in materials and workmanship
Materials	100% latex free, silicone rubber anti-static polypropylene and TPR
Shipping carton size	63x37x52cm 25200pcs per 20t container (general quantity for reference)
Other: OEM accepted	Body printing available, color alterable



Blood Bag

- ◆ Code No.: SUBDB030
- ◆ Main material : Medical grade PVC
- ◆ Bag size choice: 200, 250, 350 400, 450 500ml
- ◆ Type: Single, double, triple
- ◆ Material origin: China or Netherlands
- ◆ Function: Blood container
- ◆ Color: Transparent

Code No. X means with NP,SP,CT	Ant-Coagulant	Capacity	Needle
SU250C(X)	CPDA1	250/150/150 ml	16G
SU350C(X)	CPDA1	350/300/300 ml	16G
SU450C(X)	CPDA1	450/400/400 ml (450/300/300) ml	16G
SU500C(X)	CPDA1	500/400/400 ml	16G




Central Venous Catheter Kit

- ◆ Code No.: SUCVCK031
- ◆ Catheters made of high quality TPU, latex-free
- ◆ Type: Single lumen, double lumen, triple lumen and Quad-lumen
- ◆ Lumen size: 14G, 16G, 18G, 20G, 4F, 5F, 7F, 8F, 5.5F, 8.5F etc.
- ◆ Catheter length: 15cm, 16cm, 20cm, 30cm etc.
- ◆ Packing: 1set/Blister packing



Hemodialysis Catheter Kit

- ◆ Code No.: SUHCK032
- ◆ Catheters made of high quality TPU and transparent silicone extension tube
- ◆ Type: Single lumen, double lumen, triple lumen
- ◆ Lumen size: 8F, 6.5F, 8.5F, 11.5F, 12F etc.
- ◆ Catheter length: 10cm, 13cm, 16cm, 20cm, 30cm etc.
- ◆ Packing: 1set/Blister packing



Epidural Anaesthesia Kit

- ◆ Code No.: SUEAK033
- ◆ Epidural needle Size: 16/18G
- ◆ Epidural needle length: 80mm/90mm
- ◆ Epidural filter: 0.22 micron
- ◆ LOR Syringe: 7ml/10ml
- ◆ Epidural catheter: 900/1000mm
- ◆ Packing: 1set/hard blister
- ◆ Sterilized by EO GAS
- Content:**
- ◆ 1pc LOR syringe, 7ml, 10ml available
- ◆ 1pc Epidural needle (tuohy needle), 16/18Gx 80mm
- ◆ 1pc Epidural catheter: multiport with mark 0.8/1.0mm±850mm
- ◆ 1pc Epidural filter: 0.22um, hydrophobic membrane, luer lock
- ◆ 1pc Catheter adaptor: TBA. secure the catheter complete
- ◆ 1pc Non woven adhesive dressing



Suction Bag

- ◆ Code No.: SUSTB034
- ◆ Material: PE
- ◆ Cover Color: Pink/ Blue
- ◆ Type: Gel(with and without), Filter(with and without)

Code no.	Model	Filter	Gel	Packing	Carton size
SUSTB034-001	1200ml	No	No	60pcs/box	80*36*38cm
SUSTB034-002	1200ml	No	Yes	60pcs/box	80*36*38cm
SUSTB034-003	1200ml	Yes	Yes	60pcs/box	80*36*38cm
SUSTB034-004	1200ml	Yes	No	60pcs/box	80*36*38cm
SUSTB034-005	2000ml	No	No	60pcs/box	80*36*38cm
SUSTB034-006	2000ml	No	Yes	60pcs/box	80*36*38cm
SUSTB034-007	2000ml	Yes	Yes	60pcs/box	80*36*38cm
SUSTB034-008	2000ml	Yes	No	60pcs/box	80*36*38cm
SUSTB034-009	2800ml	No	No	60pcs/box	80*36*38cm
SUSTB034-010	2800ml	No	Yes	60pcs/box	80*36*38cm
SUSTB034-011	2800ml	Yes	Yes	60pcs/box	80*36*38cm
SUSTB034-012	2800ml	Yes	No	60pcs/box	80*36*38cm

- ◆ Application: Operating room, ICU, Ambulance and other areas with suction source
- ◆ Volume: 1200ml, 2000ml, 2800ml



Suction Canister

- ◆ Code No.: SUSTC035
- ◆ Material: PC
- ◆ Volume: 1200ml, 2000ml, 2800ml
- ◆ Application: Operating room, ICU, Ambulance and other areas with suction source

Code no.	Model	Packing	Carton size
SUSTC035-1200	1200ml	60pcs/box	80*36*38cm
SUSTC035-2000	2000ml	60pcs/box	80*36*38cm
SUSTC035-2800	2800ml	60pcs/box	80*36*38cm



Closed Wound Drainage System

- ◆ Code No.: SUCWDS036
- ◆ Made of non-toxic PVC.
- ◆ Component: 100/200ml silicone reservoir or 200mlV400ml/800ml PVC reservoir
- ◆ Flat or round fluted/perforated drain(Fr10-Fr18)
- ◆ Stainless steel trocar, Fr10, 15, 19 etc.
- ◆ Used to remove blood and excess fluid from wound.
- ◆ Packing: 1pc/Tyvek pouch, 50pcs/ctn.
- ◆ Sterilized by EO gas.



Hollow

Drainage Bottle

- ◆ Code No.: SUDNB037
- ◆ Material medical grade PVC
- ◆ Chamber single, double and triple
- ◆ Capacity: 1600ml, 2000ml, 2500ml
- ◆ Y connector for 2 drainage points, can connect 2 thoracic catheters
- ◆ Sterilized by EO gas



Standard Endotracheal Tube


- ◆ Made of medical grade PVC.
- ◆ With cuff or without cuff
- ◆ Radio opaque line provided.
- ◆ In blister packing or peelable pouch as per request
- ◆ Sterile by EO gas, single use

WITHOUT CUFFED

Code no.	Size(mm)	Code no.	Size(mm)
SUSET038-20	2.0	SUSET038-60	6.0
SUSET038-25	2.5	SUSET038-65	6.5
SUSET038-30	3.0	SUSET038-70	7.0
SUSET038-35	3.5	SUSET038-75	7.5
SUSET038-40	4.0	SUSET038-80	8.0
SUSET038-45	4.5	SUSET038-85	8.5
SUSET038-50	5.0	SUSET038-90	9.0
SUSET038-55	5.5	SUSET038-95	9.5

WITH CUFFED

Code no.	Size(mm)	Code no.	Size(mm)
SUSET038-25C	2.5	SUSET038-65C	6.5
SUSET038-30C	3.0	SUSET038-70C	7.0
SUSET038-35C	3.5	SUSET038-75C	7.5
SUSET038-40C	4.0	SUSET038-80C	8.0
SUSET038-45C	4.5	SUSET038-85C	8.5
SUSET038-50C	5.0	SUSET038-90C	9.0
SUSET038-55C	5.5	SUSET038-95C	9.5
SUSET038-60C	6.0		



Reinforced Endotracheal Tube

- ◆ Made of medical grade PVC.
- ◆ With cuff and without cuff.
- ◆ In blister packing or peelable pouch as per request
- ◆ Sterile by EO gas, single use
- ◆ 1pc/pouch, 100pcs/ctn

WITHOUT CUFFED

Code no.	Size(mm)	Code no.	Size(mm)
SURET039-20	2.0	SURET039-60	6.0
SURET039-25	2.5	SURET039-65	6.5
SURET039-30	3.0	SURET039-70	7.0
SURET039-35	3.5	SURET039-75	7.5
SURET039-40	4.0	SURET039-80	8.0
SURET039-45	4.5	SURET039-85	8.5
SURET039-50	5.0	SURET039-90	9.0
SURET039-55	5.5	SURET039-95	9.5

WITH CUFFED

Code no.	Size(mm)	Code no.	Size(mm)
SURET039-25C	2.5	SURET039-65C	6.5
SURET039-30C	3.0	SURET039-70C	7.0
SURET039-35C	3.5	SURET039-75C	7.5
SURET039-40C	4.0	SURET039-80C	8.0
SURET039-45C	4.5	SURET039-85C	8.5
SURET039-50C	5.0	SURET039-90C	9.0
SURET039-55C	5.5	SURET039-95C	9.5
SURET039-60C	6.0		



Guedel Tube Airway

- ◆ Material: semi-rigid non-toxic polyethylene
- ◆ Size: 000#/00#/0#/1#/2#/3#/4#/5#/6#
- ◆ Sterile: EO sterilization

Code no.	Size(Fr)	Color	Length(mm)	Packing
SUGTA040-000	000		40	50pcs/box, 10boxes/ctn
SUGTA040-00	00		50	50pcs/box, 10boxes/ctn
SUGTA040-0	0		60	50pcs/box, 10boxes/ctn
SUGTA040-1	1		70	50pcs/box, 10boxes/ctn
SUGTA040-2	2		80	50pcs/box, 10boxes/ctn
SUGTA040-3	3		90	50pcs/box, 10boxes/ctn
SUGTA040-4	4		100	50pcs/box, 10boxes/ctn
SUGTA040-5	5		110	50pcs/box, 10boxes/ctn
SUGTA040-6	6		120	50pcs/box, 10boxes/ctn



Silicone Laryngeal Mask Airway

- ◆ Made of imported medical-grade silicone
- ◆ Silicone smooth cuff reducing stimulation to the patient's body
- ◆ Single-use and reusable 2 types available
- ◆ Cuff can be with bar or without
- ◆ Sizes: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#
- ◆ EO gas sterilized

Code no.	Size	Maximum Cuff inflation Volumes (ml)	Patient weight (KG)	Remark
SUSLMA001	1	4	0-5	Disposable
SUSLMA002	1.5	7	5-10	
SUSLMA003	2	10	10-20	
SUSLMA004	2.5	14	20-30	
SUSLMA005	3	20	30-50	
SUSLMA006	4	30	50-70	
SUSLMA007	5	40	70-100	



PVC Laryngeal Mask Airway

- ◆ Made of medical grade PVC
- ◆ Size : Suitable for neonate , infant and children and adult
- ◆ Seven sizes: 1#, 1.5#, 2#, 2.5#, 3#, 4#, 5#
- ◆ Cuff can be with bar or without
- ◆ EO gas sterilized

Code no.	Size	Maximum Cuff inflation Volumes (ml)	Patient weight (KG)	Remark
SUPLMA001	1	4	0-5	Disposable
SUPLMA002	1.5	7	5-10	
SUPLMA003	2	10	10-20	
SUPLMA004	2.5	14	20-30	
SUPLMA005	3	20	30-50	
SUPLMA006	4	30	50-70	
SUPLMA007	5	40	70-100	



Tracheotomy Tube

- ◆ Made of clear, non-toxic PVC
- ◆ Elongated tube is used for specific patients' case
- ◆ Full length Radio-opaque line
- ◆ EO gas sterilized

WITHOUT CUFFED

Code no.	Size(mm)	Code no.	Size(mm)
SUTTT-30	3.0	SUTTT-65	6.5
SUTTT-35	3.5	SUTTT-70	7.0
SUTTT-40	4.0	SUTTT-75	7.5
SUTTT-45	4.5	SUTTT-80	8.0
SUTTT-50	5.0	SUTTT-85	8.5
SUTTT-55	5.5	SUTTT-90	9.0
SUTTT-60	6.0	SUTTT-95	9.5

WITH CUFFED

Code no.	Size(mm)	Code no.	Size(mm)
SUTTT-50C	5.0	SUTTT-75C	7.5
SUTTT-55C	5.5	SUTTT-80C	8.0
SUTTT-60C	6.0	SUTTT-85C	8.5
SUTTT-65C	6.5	SUTTT-90C	9.0
SUTTT-70C	7.0	SUTTT-95C	9.5



Adult Urine Bag

- ◆ Material: Made of PVC, non-toxic
- ◆ Sizes: 750ml, 1000ml, 1500ml, 2000ml
- ◆ Thickness: 1.2mm
- ◆ Medical PVC tube: $\phi 0.65 \times 90\text{cm}$
- ◆ Milky or transparent
- ◆ Supplied in individual peelable polybag or blister pack sterile
- ◆ Type bottom with T valve/push-pull valve/screw valve/without valve (single tube)
- ◆ EO gas sterilized



Code No.	Specification	Packing	Carton size
SUAUB001	2000ml without bottom outlet	1pc/individual pack, 10pcs/PE bag, 25bags/ctn	55x41x29cm
SUAUB002	2000ml with push-pull valve	1pc/individual pack, 10pcs/PE bag, 25bags/ctn	55x41x29cm
SUAUB003	2000ml with screw valve	1pc/individual pack, 10pcs/PE bag, 25bags/ctn	55x41x31cm
SUAUB004	2000ml with T valve	1pc/individual pack, 10pcs/PE bag, 25bags/ctn	55x41x33cm

Pediatric Urine Bag

- ◆ Material: Made of PVC, non-toxic
- ◆ Type : For infant
- ◆ Sizes: 100ml, 200ml
- ◆ Should be stored in a dry, ventilated environment
- ◆ EO gas sterilized

Code No.	Specification	Packing	Carton size
SUPUB-100	100ml	1pc/individual pack, 10pcs/PE bag, 25bags/ctn	54x28x43cm
SUPUB-200	200ml	1pc/individual pack, 10pcs/PE bag, 25bags/ctn	54x28x48cm



Luxury Urine Bag

- ◆ Material: Made of PVC, non-toxic
- ◆ Anti-reflux valve is available, to prevent urine reflux
- ◆ Volume 2000ml is the most common used volume while 4000ml
- ◆ 5000ml volumes are available for Larger urine patient
- ◆ Drainage Tube Length 90mm or 120mm
- ◆ T Valve
- ◆ Color : Frosted, transparent, translucence
- ◆ Used for 7 days collection of urine and other body fluids
- ◆ EO gas sterilized

Code No.	Specification	Packing	Carton size
SULUB2000	2000ml	1pc/individual pack, 50pcs/ctn	40x38x32cm



Surgical Blades

- ◆ Code No.: SUSBD041
- ◆ Material: Carbon steel / stainless steel
- ◆ Size: 10#, 10A, 11#, 12#, 15#, 20#, 21#, 22#, 23#, 24#
- ◆ Packing: 100pcs/box, 50boxes/ctn
- ◆ Carton size: 36*32*10.5cm
- ◆ Sterilized by Gamma



Surgical Scalpel

- ◆ Code No.: SUSSP042
- ◆ Material: Carbon steel / stainless steel
- ◆ Size: 10#, 10A, 11#, 12#, 15#, 20#, 21#, 22#, 23#, 24#
- ◆ Packing: 100pcs/box, 50boxes/ctn
- ◆ Carton size: 40*32*22cm
- ◆ Sterilized by Gamma



Electrical Surgical Pencil

- ◆ Code No.: SUESP043
- ◆ PVC or Silicone material
- ◆ Stainless steel: 4cm, diameter 2.3mm
- ◆ ABS pencil: 170cm Cable length: 3m
- ◆ Frequency confine: 0-1.0MHz
- ◆ Function: cut and coagulation
- ◆ EO gas sterilized

Spec.	Packing	Carton size	G.W./N.W.
3m silicone or PVC cable 3-bananas plug	1pc/bag, 100pcs/ctn	38*30*23cm	8KGS/7KGS



Electrical Surgical Plate

- ◆ Code No.: SUESP044
- ◆ Foam, Aluminum-foil, Hydro-gel
- ◆ Type: Bipolar and Unipolar
- ◆ Apply for child and adult
- ◆ EO gas sterilized
- ◆ Shelf life: 2 years
- ◆ Packing: 1pc/pack, 100pcs/ctn
- ◆ 50x30x25cm

Type	Jack diameter	Lead length	Conductive area	Size
Bipolar	2.3mm	3m	122cm ²	81x165mm
Unipolar	2.3mm	3m	132cm ²	105x203mm



Skin Stapler

- ◆ Code No.: SUSSP045
- ◆ Model: 15W, 25W, 35W, 15R, 25R, 35R, 35H, 25H, 15H
- ◆ Type : W: 7.0mm x 4.0mm, 0.6mm Wire Diameter Type
- ◆ R: 5.4mm x 3.6mm, 0.55mm Wire Diameter Type
- ◆ H: 7.0mm x 4.5mm, 0.6mm wire Diameter Staples
- ◆ Packing: 80pcs/ctn
- ◆ Carton Size: 46.5x45x22cm



Vaginal Speculum

- ◆ Code No.: SUVSC046
- ◆ It can be adjusted at will, with smooth edges from plastic resin-less, apprehension from patients.
- ◆ Material: medical PS ,or PC, non-toxic, non-irritant .
- ◆ Sizes: large, medium, small.
- ◆ Vaginal speculum with side screw: often used in European, Southeast Asia and some Africa countries and Latin America.
- ◆ Vaginal speculum(hook button type): used in Australia.
- ◆ Vaginal speculum with medium screw: used in Southeast Asia.
- ◆ Vaginal speculum press up trigger: used in America and Canada.
- ◆ Sterilized by EO gas, sterile, non- toxic, pyrogen free.
- ◆ Packing: PE, PP bag or blister packing, 100pcs/ctn.



Vaginal Douche

- ◆ Code No.: SUVDC047
- ◆ Material: medical grade PVC and natural rubber
- ◆ Color: blue, black, red, white or pink etc
- ◆ Capacity: 310ml, 180ml, 120ml
- ◆ Usage: enema for gynecology
- ◆ Packing: Single OPP bag or Individual box

Pkg(pcs/ctn)	Ctn size(cm)	G.W.(kgs)
1pc/opp bag or 1pc/box, 100pcs/ctn	56x46x32cm	10kgs




Cervical Scraper

- ◆ Code No.: SUCVS048
- ◆ Material: plastic or wooden
- ◆ Sterilized by EO Gas
- ◆ Packing: 1pc/pouch, 100pcs/box, 2000pcs/ctn
- ◆ Carton size: 54x23x45cm
- ◆ Application: for gynecological examination diseases.



Cervical Brush

- ◆ Code No.: SUCVB049
- ◆ Material: plastic with bristle
- ◆ Color: blue, pink, white etc
- ◆ Usage: gynecological examination
- ◆ Packing: 1pc/pouch, 100pcs/box, 2000pcs/ctn
- ◆ Carton size: 54x23x45cm



Umbilical Cord Clamp

- ◆ Code No.: SUUCC050
- ◆ Medical grade PVC
- ◆ For new born use
- ◆ Blue, white, pink
- ◆ 5cm length
- ◆ Sterilized by EO GAS
- ◆ Packing: 1pc/PE bag



I.D Bracelet

- ◆ Code No.: SUIDB052
- ◆ Medical grade PVC
- ◆ Adult, child, infant
- ◆ Pocket insert card type & write-on type
- ◆ Color: Blue, red, green, yellow, orange etc



Hot Water Bag

- ◆ Code No.: SUHWB054
- ◆ Material: PVC or natural rubber
- ◆ Capacity: 500ml, 750ml, 1000ml, 1500ml, 2000ml
- ◆ Packing: 1pc/OPP bag .50pcs/ctn
- ◆ Carton size: 63x45x58cm or other



Tourniquet

- ◆ Code No.: SUTQT051
- ◆ Disposable or reusable
- ◆ Adult and Child
- ◆ Size: 2.5*45cm, 2.5*40cm, 2.5*35cm
- ◆ Silk Printing or Heat transfer printing
- ◆ Color: Blue, red, green, yellow, orange
- ◆ Packing: 1pc/PE bag or 100pcs/box



Hot Cold Pack

- ◆ Code No.: SUHCP053
- ◆ Material: Fabric outside, Gel inside
- ◆ Size: 23*13cm, 250g
- ◆ Packing: 56pcs/ctn
- ◆ Carton: 42*25*24cm
- ◆ Temperature range: 5°C to 70°C



Ice Bag

- ◆ Code No.: SUICB055
- ◆ Material: Polyester and PVC
- ◆ Size: 7", 9", 11"
- ◆ Packing :1pc/bag.50pcs/ctn
- ◆ Carton size: 63x45x58 or other



Stool Container

Code No.	Spec.	Material	Packing	Carton Size
SUSTC056-001	30ml	PS	100pcs/bag, 1000pcs/carton	56*42*42cm
SUSTC056-002	30ml	PP	100pcs/bag, 1000pcs/carton	57*42*52cm
SUSTC056-003	60ml	PP	100pcs/bag, 500pcs/carton	57*38*40cm
SUSTC056-004	100ml/120ml	PP	100pcs/bag, 500pcs/carton	59*38*65cm



Urine Container

Code No.	Spec.	Material	Packing	Carton Size
SUURC057-001	30ml	PS	100pcs/bag, 1000pcs/carton	56*42*42cm
SUURC057-002	30ml	PP	100pcs/bag, 1000pcs/carton	57*42*52cm
SUURC057-003	60ml	PP	100pcs/bag, 500pcs/carton	57*38*40cm
SUURC057-004	100ml/120ml	PP	100pcs/bag, 500pcs/carton	59*38*65cm



Microscope Slide

Code No.	Specification	Packing	Carton size
SUMS7101	25.4*76.2mm	50pcs/box, 50boxes/ctn	36*18*15cm
SUMS7102	25.4*76.2mm	72pcs/box, 50boxes/ctn	44*20*15cm
SUMS7103	25.4*76.2mm	50pcs/box, 50boxes/ctn	36*18*15cm
SUMS7104	25.4*76.2mm	72pcs/box, 50boxes/ctn	44*20*15cm
SUMS7105-1	25.4*76.2mm	50pcs/box, 50boxes/ctn	36*18*15cm
SUMS7107	25.4*76.2mm	72pcs/box, 50boxes/ctn	44*20*15cm
SUMS7107-1	25.4*76.2mm	50pcs/box, 50boxes/ctn	36*18*15cm

Cover Glass


Code No.	Specification	Packing	Carton size
7201	18*18mm	100pcs/boxes, 500boxes/carton	36*21*16cm
	20*20mm	100pcs/boxes, 500boxes/carton	36*21*16cm
	22*22mm	100pcs/boxes, 500boxes/carton	37*25*19cm
	22*50mm	100pcs/boxes, 250boxes/carton	41*25*17cm
	24*24mm	100pcs/boxes, 500boxes/carton	37*25*17cm
	24*32mm	100pcs/boxes, 400boxes/carton	44*27*19cm
	24*40mm	100pcs/boxes, 250boxes/carton	41*25*17cm
	24*50mm	100pcs/boxes, 250boxes/carton	41*25*17cm
	24*60mm	100pcs/boxes, 250boxes/carton	46*27*20cm



Pipette

- ◆ Made of PP material
- ◆ Normal colors: yellow, blue, clear
- ◆ With or without scale
- ◆ Tip for Eppendorf, Gilson, Finland, Dalong etc
- ◆ Volume from 10ul to 500ul

Code No.	Description	Specification	Material	Pack quantity	Carton Size
SUPPT001	Tips for Eppendorf pipette	10ul white	PP	50000	56x42x30cm
SUPPT002	Tips for Eppendorf pipette	200ul yellow	PP	50000	56x42x42cm
SUPPT004	Tips for Eppendorf pipette	1000ul blue	PP	15000	56x42x42cm
SUPPT005	Tips for Gilson pipette	10ul white	PP	140000	56x42x42cm
SUPPT006	Tips for Gilson pipette	200ul yellow	PP	50000	56x42x42cm
SUPPT007	Tips for Gilson pipette	1000ul blue	PP	15000	56x42x42cm
SUPPT008	Tips for Finland pipette	10ul white	PP	140000	56x42x42cm
SUPPT009	Tips for Finland pipette	200ul white	PP	50000	56x42x42cm
SUPPT010	Tips for Finland pipette	1000ul white	PP	15000	56x42x42cm
SUPPT011	Tips for Gilson pipette	5000ul white	PP	7500	56x42x42cm
SUPPT012	Filter tips for Gilson	10ul white	PP	140000	56x42x42cm
SUPPT013	Filter tips for Gilson	200ul yellow	PP	50000	56x42x42cm
SUPPT014	Filter tips for Gilson	200ul white	PP	50000	56x42x42cm
SUPPT015	Filter tips for Gilson	1000ul blue	PP	15000	56x42x42cm



Bulb Irrigation Syringe

- ◆ Code No: SUBIS058
- ◆ 50ml, 60ml, 80ml, 100ml available
- ◆ Non-toxic PVC
- ◆ PE bag, pouch or blister
- ◆ Sterilized by EO



Ear Wash bulb

- ◆ Code No: SUEWB059
- ◆ Material: PVC
- ◆ Feature: Flush syringe with water before using
- ◆ Specification: 30ml, 60ml, 90ml



Urine Dip Strips

- ◆ Types: 10A, 11A, 12A, 14A
- ◆ Application: used in clinical for urine leukocyte, nitrite, protein, PH, blood, specific gravity, glucose, Ca, Creatinine, Ascorbic acid, Microalbumin for qualitative and semi quantitative detection

Code No.	Item	Packing	Carton size
SUUDS060	Urine dip strips	100strips/barrel, 100barrels/bottle	59*39*31cm



Flocking swab

- ◆ Code No: SUFKS061
- ◆ Material: sponge(foam head)&handle(plastic/wood)
- ◆ Color: white,blue etc
- ◆ Outer dimension: 150mm/190mm/13x100mm
- ◆ Usage: for gynecological examination
- ◆ Feature: disposable, for single use only



Penrose drainage tube


- ◆ Code No: SUPDT062
- ◆ Material: natural latex
- ◆ Size: 1/8"1/4", 3/8", 1/2", 5/8", 3/4", 7/8", 1"
- ◆ Length: 12-17
- ◆ Usage: for surgical wound drainage
- ◆ Packed: 1pc in an individual blister bag, 100pcs/ctn



Petri Dish

- ◆ Material: PS
- ◆ Type: one room,two rooms,three rooms
- ◆ Color: transparent
- ◆ Size: 35x15mm,60x15mm,70x15mm,90x15mm,150x15mm etc
- ◆ Shape: round,square

Code No.	Size	Material	Shape	Packing
SUPTD001	35mm*15mm	PS	Round	2000pcs/ctn
SUPTD002	55mm*15mm			1000pcs/ctn
SUPTD003	60mm*15mm			1000pcs/ctn
SUPTD004	70mm*15mm			960pcs/ctn
SUPTD005	90mm*15mm(1room)			600pcs/ctn
SUPTD006	90mm*15mm(2 rooms)			600pcs/ctn
SUPTD007	90mm*15mm(3rooms)			600pcs/ctn
SUPTD008	150mm*15mm			110pcs/ctn
SUPTD009	130mm*13mm		Square	150pcs/ctn



Transport swab

- ◆ Code No: SUTPS063
 - ◆ Stick material: wooden or plastic
 - ◆ Size: 12x75mm/13x150mm
 - ◆ Type: for female or for male
 - ◆ Packing: 1000pcs/pack,2000pcs/ctn
 - ◆ Used for sampling from the parts of vagina,urogenital canal etc
- Transport swab with medium
- Classified : amies medium, amies with charcoal, stuart medium, stuart charcoal, cary-blair medium and cary-blair with charcoal adapted to sampling from the parts of pharyngeal, nose, ear, eye, vagina, urogenital canal, skin,intestinal tract etc



Washing Bottle

- ◆ Material: PE
- ◆ Size: 125/250/500/1000ml
- ◆ Type: red mouth and white mouth
- ◆ Packing: 600, 400, 200, 100pcs/ctn
- ◆ Used in a chemical laboratory

Code no.	Capacity (ml)	Total weight (g)	Internal diameter (±0.2mm)	Bottle diameter (±0.2mm)	Bottle height (±0.2mm)	Total height (±0.2mm)
SUWBT250	250	34.2	22.5	62.5	130	175



Single Channel Pipette

- ◆ Code No.: SUSCP064
- ◆ Material: ABS+PA+PC
- ◆ Single channel,8 channel,12 channel
- ◆ Type: Fixed volume and variable volume
- ◆ Fixed volume: 5ul,10ul,20ul,25ul,50ul,100ul,200ul,250ul,500ul,1000ul,2000ul,5000ul,1ml etc.
- ◆ Variable volume: 0.1-0.5ul,0.5-10ul,2-20ul,5-50ul,10-100ul,20-200ul,50-500ul,100-1000ul, 200-1000ul,1000-5000ul,2-10ml etc.
- ◆ Packing: 1pc/box,50pcs/ctn
- ◆ Increment: 0.1vl,0.5vl,5vl etc



Pasteur Pipette

- ◆ Material: PE
- ◆ Capacity: 1ml,2ml,3ml,5ml,10 etc
- ◆ Length: 145mm/148mm/150mm/155mm/158mm/160mm/178mm/185mm/200mm etc
- ◆ Type: Sterile and Non sterile

Code No.	Description	Volume/Length	Material	Qty/pk	Qty/Ctn
SUPPT-001	Graduated pipette, Non-sterile	1ml,145mm	PE	500	5000
SUPPT-002	Graduated pipette, Individually, wrapped, EO sterile	1ml,145mm	PE	250	5000
SUPPT-003	Graduated pipette, Non-sterile	1ml,148mm	PE	500	5000
SUPPT-004	Graduated pipette, Non-sterile	2ml,150mm	PE	500	5000
SUPPT-005	Graduated pipette, Non-sterile	3ml,155mm	PE	500	5000
SUPPT-006	Graduated pipette, Non-sterile	3ml,155mm	PE	500	5000
SUPPT-007	Graduated pipette, Individually, wrapped, EO sterile	3ml,155mm	PE	250	5000
SUPPT-008	Graduated pipette, EO sterile	3ml,155mm	PE	500	5000
SUPPT-009	Graduated pipette, Non-sterile	3ml,158mm	PE	500	4000
SUPPT-010	Graduated pipette, Non-sterile	3ml,178mm	PE	500	5000
SUPPT-011	Graduated pipette, Non-sterile	3ml,185mm	PE	500	5000
SUPPT-012	Graduated pipette, Non-sterile	5ml,200mm	PE	500	5000
SUPPT-013	Graduated pipette, Non-sterile	10ml,295mm	PE	500	5000



Centrifuge Tube

- ◆ Code No.: SUCFT065
- ◆ Material: PP
- ◆ Capacity: 0.2ml, 0.5ml, 1.5ml, 2ml
- ◆ With screw cap, round or conical bottom
- ◆ With printed graduation
- ◆ Packing: 500pcs/polybag, 10000pcs/ctn



Plastic Test Tube

- ◆ Material: PS, PP, Glass
- ◆ Capacity: 3ml, 5ml, 7ml, 10ml
- ◆ Color: Transparent
- ◆ All of our test tube can be printed (Graduated)

Code No.	Size(mm)	Volume	Material	Packing(bag/ctn)	Carton Size
SUPTT-001	12*55, U-bottom	3ml	PS	500/5000pcs	46*32*39mm
SUPTT-002	12*60	4ml	PS	500/5000pcs	46*32*39mm
SUPTT-003	12*75	5ml	PS	250/5000pcs	46*32*39mm
SUPTT-004	12*75	5ml	PP	250/5000pcs	46*32*39mm
SUPTT-005	13*78	6ml	PS	200/5000pcs	51*32*42mm
SUPTT-006	12*100	7ml	PS	100/4000pcs	51*32*42mm
SUPTT-007	12*100	7ml	PE	100/4000pcs	51*32*42mm
SUPTT-008	13*100	8ml	PS	100/4000pcs	51*32*42mm
SUPTT-009	13*100	8ml	PP	100/4000pcs	51*32*42mm
SUPTT-010	16*100	14ml	PS	100/2400pcs	51*32*42mm
SUPTT-011	16*100	14ml	PP	100/2400pcs	51*32*42mm



Measuring Cup

- ◆ Material: PP, Glass
- ◆ Capacity: 10ml, 20ml, 30ml, 100ml, 500ml etc
- ◆ Color: Solid/transparent

Code No.	Volume	Size/mm
SUMSC-001	10ml	Φ32*H35
SUMSC-002	10ml	Φ31*H31
SUMSC-003	20ml	Φ36*H35
SUMSC-004	20ml	Φ37*H40
SUMSC-005	30ml	Φ43*H35
SUMSC-006	50ml	Φ45*H48
SUMSC-007	100ml	Φ55*H67
SUMSC-008	150ml	Φ62*H70
SUMSC-009	500ml	90*131*125
SUMSC-010	1200ml	121*170*185
SUMSC-011	2000ml	151*208*214
SUMSC-012	3000ml	168*228*240
SUMSC-013	5000ml	200*250*270



Autoclave Bag

- ◆ Code No.: SUACB066
- ◆ Material: LDPE, HDPE, PP/Biodegradable material.
- ◆ Size: 19*23inch, 23*35inch, 38*48inch, 50*75cm, 60*90cm, 100*120cm
- ◆ Color: Yellow, Red, Black etc
- ◆ Usage: Collecting medical rubbish
- ◆ Durable, Waterproof, Eco-friendly



Heparin Cap

- ◆ Code No.: SUHPC068
- ◆ Material: Medical grade PVC
- ◆ Color: Transparent or yellow
- ◆ Packing: 1pc/blister packing, 100pcs/box, 2000pcs/ctn
- ◆ Carton size: 31* 28*24cm



Lubricating Jelly

- ◆ Code No.: SULCJ070
- ◆ Greaseless, transparent, water soluble and non-irritating
- ◆ Non-irritating, Non-staining, Non-greasy.
- ◆ Size: 3g, 5g, 10g, 30g, 35g, 42g, 45g, 82g, 85g, 100g and 150g, Etc.
- ◆ Packing: Tubes, Sachets and Syringe packing .
- ◆ It's recommended for gynecological and surgical lubrication and when additional vaginal lubrication is needed.



Vomit Bag

- ◆ Code No.: SUVMB067
- ◆ Size : 17x37.5cm, custom required
- ◆ Material: virgin PE
- ◆ Color: blue/white etc.
- ◆ Thickness: 0.03mm per layer or customer required
- ◆ Weight: around 17g per piece
- ◆ Usage: for patients, pregnant woman, carsickness and seasickness people



Hand Scrub Brush

- ◆ Code No.: SUHSB069
- ◆ Component: brush 8*5*2.1cm, sponge and PVC
- ◆ Color: Yellow or White
- ◆ Type : With Chlorhexidine Digluconate, Providine-iodine, parachlorometaxyleno, Isopropyl alcohol available.
- ◆ Packing: 250 sets/ctn 54x43x33cm



Alcohol Lamp

- ◆ Code No.: SUACL071
- ◆ Material: Boro 3.3 Glass
- ◆ Capacity: 150ml, 250ml etc
- ◆ Packing: 90pcs/ctn
- ◆ Carton size : 53x42x34cm



Medicine Bag

◆ Code No.: SUMCB072

Item	Pill bag, Medicine bag
Material	Laminated material/Multi payer: LDPE
Size	9x9cm, 10x9cm, 10x12cm, 12x15cm, 12x18cm etc
Shape	Zip-lock bag, back sealing bag
Printing	Tablet name, instruction, patient name, date, 4-pictures (time) etc
Advantage	Waterproof, airproof, dampproof
Packing	50, 100pcs/pack, 10, 20packs/ctn



Dropper Glass Bottle

- ◆ Code No.: SUDGB073
- ◆ Material: Glass
- ◆ Color: Amber
- ◆ Volume: 5ml-100ml



Pills Box

- ◆ Code No.: SUPLB075
- ◆ Material: ABS stand PP
- ◆ Product Size: 28.4 *10.4 *2.8cm
- ◆ Function: Pills storage for 31 days max



Pop Up Vial

Code no.	Spec.	Packing	Carton Size	Material
SUPUV13DR	13DR	325	33*31*38cm	PP
SUPUV19DR	19DR	230	33*31*38cm	PP
SUPUV30DR	30DR	155	33*31*38cm	PP
SUPUV60DR	60DR	75	33*31*38cm	PP



Plastic Jars

- ◆ Code no.: SUPTJ074
- ◆ Material: PET
- ◆ Color: clear
- ◆ Capacity: 50g, 100g



Pill Cutter

- ◆ Code No.: SUPCT076
- ◆ Material: ABS
- ◆ Function: cut pills, store pills, crush pills
- ◆ Size: 4.1*5*6.4cm
- ◆ Packing: 1pc/box, 360boxes/ctn, 45*30*34cm



Autoclave Sterilization Indicator Tape

- ◆ Innovative Designs as per your request
- ◆ Sterilization method: STEAM and ETO
- ◆ Sizes: 19mm*50m
- ◆ Material: masking tape

Sterilization Indicator Tape					
Code no.	Size	Roll/Carton	Carton Size L*W*H(mm)	N.W(kg)	G.W(kg)
	W(mm)*L(m)				
SUSIT-001	12*50	76	27.5*27.5*24	5.5	6
SUSIT-002	19*50	48	27.5*27.5*24	5.5	6
SUSIT-003	25*50	36	27.5*27.5*24	5.5	6



Tyvek Sterilization Reel

Code no	Size	Roll/Carton	Carton Size
SUTSR001	50mm*100m	10rolls/ctn	400*200*315mm
SUTSR002	100mm*100m	6rolls/ctn	400*200*315mm
SUTSR003	150mm*100m	4rolls/ctn	400*200*315mm
SUTSR004	200mm*100m	2rolls/ctn	400*200*215mm
SUTSR005	250mm*100m	2rolls/ctn	400*200*265mm
SUTSR006	300mm*100m	2rolls/ctn	400*200*315mm
SUTSR007	350mm*100m	2rolls/ctn	400*200*375mm
SUTSR008	400mm*100m	2rolls/ctn	400*200*415mm
SUTSR009	500mm*100m	2rolls/ctn	400*200*515mm

Auto-sealing(self-sealing)Pouches

Code no	Size	Roll/Carton	Carton Size
SUASP-001	70*260mm	200pcs/box, 15boxes/ctn	320*280*245mm
SUASP-002	90*260mm	200pcs/box, 10boxes/ctn	320*280*210mm
SUASP-003	135*280mm	200pcs/box, 10boxes/ctn	320*300*300mm
SUASP-004	190*360mm	200pcs/box, 5boxes/ctn	380*320*215mm
SUASP-005	200*320mm	200pcs/box, 5boxes/ctn	340*320*225mm
SUASP-006	200*360mm	200pcs/box, 5boxes/ctn	380*320*225mm
SUASP-007	230*395mm	200pcs/box, 5boxes/ctn	420*320*255mm
SUASP-008	305*430mm	200pcs/box, 5boxes/ctn	450*320*330mm
SUASP-009	400*535mm	200pcs/box, 5boxes/ctn	550*320*425mm

Heat-Sealing Flat Reel Pouch

Code no	Size	Roll/Carton	Carton Size
SUHRFP001	50mmx200m	12rolls/ctn	430*210*315mm
SUHRFP002	75mmx200m	8rolls/ctn	430*210*315mm
SUHRFP003	100mmx200m	6rolls/ctn	430*210*315mm
SUHRFP004	150mmx200m	4rollg/ctn	430*210*315mm
SUHRFP005	200mmx200m	2rolls/ctn	430*210*215mm
SUHRFP006	250mmx200m	2rolls/ctn	430*210*265mm
SUHRFP007	300mmx200m	2rolls/ctn	430*210*315mm
SUHRFP008	400mmx200m	2rolls/ctn	430*210*415mm
SUHRFP009	500mmx200m	2rolls/ctn	430*210*515mm
SUHRFP010	600mmx200m	2rollg/ctn	430*210*615mm

Heat-Sealing Gusseted Pouches

Code no	Size	Roll/Carton	Carton Size
SUHSGP001	5mmx100m	10rolls/ctn	400*200*315mm
SUHSGP002	75mmx100m	8rollg/ctn	400*200*315mm
SUHSGP003	100mmx100m	6rolls/ctn	400*200*315mm
SUHSGP004	150mmx100m	4rollg/ctn	400*200*315mm
SUHSGP005	200mmx100m	2rollg/ctn	400*200*215mm
SUHSGP006	250mmx100m	2rollg/ctn	400*200*265mm
SUHSGP007	300mmx100m	2rollg/ctn	400*200*315mm
SUHSGP008	400mmx100m	2rollg/ctn	400*200*415mm



Sport Support

Code No.	Type	Material	Size	Packing	Carton size
SUSST001	Arm support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm
SUSST002	Ankle support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm
SUSST003	Knee support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm
SUSST004	Waist support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm
SUSST005	Shoulder support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm
SUSST006	Calf support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm
SUSST007	Palm support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100 pcs/ctn	36x43x45cm
SUSST008	Thigh support	Nylon,Neoprene,Rubber, Polyester,Elastic	S/M/L/XL	1pc/box,100pcs/ctn	36x43x45cm



Skin Traction Kit

- ◆ Style: child and adult
- ◆ Material: bandage and foam,belt and tape, Bandage has two style, adhesive or non-adhesive
- ◆ Package:1pc/box, 100pcs/ctn

Code No.	Size	Packing	Carton size
SUSTK077A	Adult	1pc/box,100boxes/ctn	50×37×58cm
SUSTK077C	Child	1pc/box,100boxes/ctn	50×37×58cm



Sportwatch

Code No.: SUSTW077
This mode counts in 1/100th a second interval for the first 30 minutes of timing(29min 59.999 Sec.)
--1 second increment up to maximum of 24 hours
--Split time,total time and lap time can be displayed
--It can record 2 laps memories
--Display hour,minute,second,year, month,day,week and Auto Calendar
--Able to set the Alarm,snooze function
--12 hours or 24 hours time display,Hourly Chime
--Any color is available
--Battery included(Lithium cell*1pcs)



Qty: 100pcs
Ctn size: 44*29.5*28.5cm
G.W: 9.8KG
N.W: 8.8KG

Sharp Container/Safety box

- ◆ Material: plastic or paper
- ◆ Color: yellow, red etc
- ◆ Capacity: 1L,3L,5L,10L,12L etc
- ◆ Application: medical waste container

Code No	Products size	Capacity(L)	Pcs/ctn	Carton size
SUSPC-001	9.5*7.5*13.5cm	0.7	100	35*35*46cm
SUSPC-002	12*9*11.5cm	0.8	90	38*38*39.5cm
SUSPC-003	11*8.5*15cm	1	100	36*33*46cm
SUSPC-004	16*14*15cm	2.8	30	36*36*37cm
SUSPC-005	22*18*24.5cm	8	24	53*42*45cm
SUSPC-006	14*11.5*20cm	2.0	40	46*32*41cm
SUSPC-007	14*12*24cm	3.0	30	48*32*48cm
SUSPC-008	22.5*20*19cm	6.5	24	44*44*40.5cm
SUSPC-009	16.5*14*19cm	3.6	100	52*52*57cm
SUSPC-010	22*19*20cm	6.0	24	44*44*40.5cm



Rubbish Bin With Lid And Wheels

- ◆ Material: VIRGIN HDPE
- ◆ Capacity: 100L
- ◆ Color: Red, dark gray, blue, green, yellow, green army
- ◆ Uses: Square, green space, commercial street, residential, home or medical

Code No	TYPE	A(mm)	B(mm)	C(mm)
SURBB-001	GT-240D	720±5	570±5	1050±5
SURBB-002	GT-240H	715±5	560±5	1000±5
SURBB-003	GT-120A	550±5	470±5	900±5
SURBB-004	GT-100A	520±5	470±5	810±5
SURBB-005	GT-50A	460±5	420±5	640±5
SURBB-006	GT-25A	365±5	310±5	455±5
SURBB-007	GT-15A	265±5	255±5	345±5



Urine Bottle

- ◆ Code No.: SUUBT078
- ◆ Material: PE
- ◆ Type: Male,Female
- ◆ Size: 800ml,1000ml
- ◆ Packing: 66pcs/ctn,60pcs/ctn



Bedpan

- ◆ Code No.: SUBDP079
- ◆ Material: PP
- ◆ Color: White or blue.
- ◆ Size: 510*290*90mm
- ◆ Feature: with handle, with or without cover
- ◆ Packing: 1pc/bag,20pcs/carton Carton size: 52*30*60cm



Dental Dressing Set

- ◆ Code No.: SUDDS080
- ◆ Material: ABS+PP+Paper
- ◆ Color: Blue, white
- ◆ Items: Probe, Mouth mirror, Forceps, Dental bib, Saliva ejector, Tray
- ◆ Probe: Single ended dental probe
- ◆ Sterilized: Sterilized by EO gas

Components	Description	Quantity
Dental probe	Single ended or double ended	1 piece
	PVC+ stainless steel	
Dental tweezers	Iron coated zinc, Iron coated nickel, or plastic	1 piece
	PVC+ stainless steel	
Mouth mirror	Plastic or glass lens	1 piece
Dental bib	Dental bib with tie	1 piece
Saliva ejector	PVC tube, copper interior wiring	1 piece
	with chrom plated 150*6.5mm	
Strong straw	straight	1 piece
Strong straw	bending	1 piece
operation pipe	PVC	1 piece
Cotton roll	cotton	2 pieces
Face mask	3-ply Ear-loop	1 piece
Tray	PVC	1 piece



Gynecological Examination Set

- ◆ Code No.: SUGES081
- ◆ Name: Disposable gynecological examination kit/packs/set
- ◆ Type: Lateral screw type, push type, middle screw type, fastener type
- ◆ Component of the kits: vaginal speculum/gloves/cervical scraper/non-woven napkin
- ◆ Sterilization: EO gas OEM service: available
- ◆ Package: paper and plastic bag, or PE bag, 100sets/ctn

Content	Quantity
Disposable vaginal speculum	1pc
Natural latex examination gloves with powdered (S/M/L)	1pair
Disposable cervical scraper	1pc
Non-woven napkin (50*50cm)	1pc
PE bag/plastic and paper bag	1pc



Wound Dressing Set

- ◆ Code No.: SUWDS082

Item	Description (wound dressing set)
SUM-15080	Basic dressing set
Contents	quantity
Dressing forceps	1
Tray	1
Cotton ball	6
Cotton gauze swabs 7.5x7.5cm 8ply	6
Hand towel	1
Patient drapes	2
Yellow Clinical bag	1
LDPE apron 125cm	1



Tongue Depressor(Spatula)

- ◆ Code No.: SUTDP083
- ◆ Sterile or Non sterile
- ◆ Material: Wooden / Plastic
- ◆ Size: Adult / Child
- ◆ Packing: 100pcs/box, 50boxes/carton
- ◆ Carton size: 41.5*31.5*20cm



Mouth Mirror

- ◆ Code No.: SUMMR085
- ◆ Material Polymer
- ◆ Length: 15cm, 15.8cm, 16cm
- ◆ Diameter: 2.2cm, 2.4cm
- ◆ Carton: 57*45*29cm
- ◆ Packaging: 150pcs/box, 2700pcs/ctn



Surgical Tip

- ◆ Code No.: SUSGT087
- ◆ Material of tubes: PVC (latex free, non-toxic)
- ◆ Color: Green, white, blue
- ◆ Size: Large: 5.8mm, Medium: 3.8mm, Small: 2.5mm



Skin marker pen

- ◆ Code No.: SUSMP089
- ◆ Material: PP
- ◆ Market size: 166mmx10mm
- ◆ Nib size: 0.5mm
- ◆ Ink color: gentian violet & blue
- ◆ Usage: surgical skin marking, striated eyebrow marking
- ◆ Packing: 1pc/pouch, 25pcs/box 40boxes/ctn
- ◆ Carton size: 50*39*46cm



Denture Box

- ◆ Code No.: SUDTB084
- ◆ Material: PP
- ◆ Color: Yellow, Blue, Green, White, Pink
- ◆ Size: 95*73*56mm, 95*74*39mm, 84*75*78mm, 86*78*78mm,
- ◆ 100*98*86mm
- ◆ Different shape available



Dental Tweezers

- ◆ Code No.: SUDTT086
- ◆ Material: ABS/Metal
- ◆ Color: blue, white etc
- ◆ Size: L,M,S
- ◆ Packaging: 100pcs/box, 1800pcs/ctn
- ◆ Carton: 57*45*29cm



Disposable Razor

- ◆ Code No.: SUDSR088
- ◆ Model: single blade, double blade
- ◆ Material: plastic handle with stainless steel
- ◆ Color: blue, green, customized color



Dental Bib

- ◆ Code No.: SUDTB090
- ◆ Material: 2-ply paper+ 1-ply poly film or 1-ply paper+ 1-ply film
- ◆ Color: Red, green, yellow, blue, white etc
- ◆ Size: 33cmx45cm(Standard)
- ◆ Weight: Paper: 16-19gsm, Film:12-14gsm



Dental Probe

- ◆ Code no.: SUDTP092
- ◆ Material: ABS
- ◆ Color: White .Blue
- ◆ Size: S,M,L
- ◆ Packing: one piece in one plastic bag ,1000 pcs in one carton



Dental Irrigation Syringe

- ◆ Code No.: SUDIS094
- ◆ Usage: Used in cleaning a certain part in the dental operation.
- ◆ Material: High transparent medical grade PP, synthetic rubber(latex free)or natural rubber. Medical grade stainless steel SUS304, medical grade silicone
- ◆ Type: with needle, without needle



Endoscopy Mouth piece

- ◆ Code No.: SUEMP096
- ◆ Material: medical polymer materials.
- ◆ Color: blue, white, green etc
- ◆ Usage: for endoscopy and endoscopic surgery to prevent patients chewed endoscopy tube.
- ◆ Packing: 50pcs/box,10boxes/ctn Carton size: 59*29*35cm



Dental Protective Film

- ◆ Code No.: SUDPF091
- ◆ Material: PE
- ◆ Size: 4"x6"x1200pcs
- ◆ Thickness: 0.05mm
- ◆ Color: Blue,Clear
- ◆ Length: 180m
- ◆ G.N: 8/7.5kg



Dental Apron

- ◆ Code no.: SUDTA093
- ◆ Material: 1-ply paper+ 1-ply PE poly
- ◆ Color: White,Blue
- ◆ Size: 36cm*64cm
- ◆ Weight: Paper: 24gsm, Film: 14gsm



Dental Irrigation Needle

- ◆ Code No.: SUDIN095
- ◆ Usage: Used in cleaning the wound in dental treatment.
- ◆ Material: Medical grade stainless steel SUS304.
- ◆ Type: Double-Vented, Side-Vented, Notched
- ◆ Size: 23G, 27G, 30G



X Ray Film

- ◆ Type: blue sensitive or green sensitive
- ◆ Storage Humidity: 40-60%
- ◆ Storage Temperature : 10-20°C/50-70°F
- ◆ Size: 8x10in; 11x14in; 14x14in etc

Code No.	Specification(inch)	Quantity(Pcs)	Quantity/Carton
SUXRF-001	8*10	100/box	5boxes
SUXRF-002	10*12	100/box	5boxes
SUXRF-003	11*14	100/box	5boxes
SUXRF-004	12*15	100/box	5boxes
SUXRF-005	14*14	100/box	5boxes
SUXRF-006	14*17	100/box	5boxes
SUXRF-007	18*24cm	100/box	5boxes
SUXRF-008	30*340cm	100/box	5boxes



X-Ray Film Envelopes

- ◆ Material: Kraft paper
- ◆ Size: 14***14", 10***12", 8***10", 11***14", 14***17" etc
- ◆ Medical X-ray Film Blue, has the characteristics of low fog, Powerful contrast,quick sensitivity and high resolution ratio. It can supply photos with high -clarity and high-penetration. It can be mainly used in the common x-ray photos for bones or chest.

Code No.	Specification	Packing	Carton Size
SUXRFE001	14***14"	500pcs/ctn	38x38x32cm
SUXRFE002	10***12"	500pcs/ctn	38x38x32cm
SUXRFE003	8***10"	500pcs/ctn	38x38x32cm
SUXRFE004	11***14"	500pcs/ctn	38x38x32cm
SUXRFE005	14***17"	500pcs/ctn	38x38x32cm



X Ray Film Developer

- ◆ Code No.: SUXRFD097
- ◆ Material: Hydroquinone,Anhydrous Sodium Sulfate,Phenitane,Acetic Acid,Glutaric dialdehyde
- ◆ Contents: A, 5L-----1 bottle
B, 0.5L-----1 bottle
C, 0.5L-----1 bottle
- ◆ Packing: 2sets/ctn or customer OEM
- ◆ Store in an environment temperature below 45°C, keep dry and ventilated



X Ray Film Fixer

- ◆ Code No.: SUXRFF098
- ◆ Material: Ammonium thiosulfate, Anhydrous Sodium Sulfate, Orthoboric acid, Aluminium sulfate, etc.
- ◆ Contents: A, 5L-----1 bottle
B, 2.5L-----1 bottle
- ◆ Packing: 2sets/ctn or customer OEM
- ◆ Store in an environment temperature below 45°C, keep dry and ventilated



ECG Electrode Plate

- ◆ Code No.: SUEEP099
- ◆ Electrode Diameter: 48x48mm
- ◆ Wet gel, Snap electrode
- ◆ Silver/Silver Chloride(Ag/AgCl)Sensor
- ◆ Non-woven fabric
- ◆ Packing: 50pieces/Bag, 1500pcs/carton
- ◆ Carton size : 45*37*23.5cm



Electrical Conductive Gel

- ◆ Code No.: SUECG100
- ◆ Material: Water-soluble polymer conductive gel
- ◆ Color: white, blue, transparent
- ◆ Capacity: 60ml,100ml,250ml,1000ml etc
- ◆ Packing: 200pcs/ctn,100pcs/ctn,50pcs/ctn,40pcs/ctn,2 bottles/ctn
- ◆ Directly used in electrode pads can also be used on the skin



Ultrasound Thermal Paper

- ◆ Synthetic Thermal Paper Matt or Glossy(film).
- ◆ Each roll is packed in individual pack
- ◆ High resolution and preservation of printed video image
- ◆ 110S,110HG,110HD

Code No.	Specification	Size	Quantity	Package	Carton size
SUUTP-001	UTP-110S (Normal)	110mm*20m	Print 240 copies of each roll	10roll/box,10box/ctn	57x35x25cm
SUUTP-002	UTP-110HG (High Glossy)	110mm*18m	Print 220 copies of each roll	10roll/box,10box/ctn	57x35x25cm
SUUTP-003	UTP-110HD (High-density)	110mm*18m	Print 220 copies of each roll	10roll/box,10box/ctn	57x35x25cm



Ultrasound Gel

Features:

- ◆ Better sound transmission, distinctive image display.
- ◆ No toxic and irritant effect on the skin
- ◆ No damage and corrosion to the ultrasound probe.

Technical Parameters:

- ◆ Sound Speed 1520-1620m/s
- ◆ Sound impedance Frequency 1.5x 10⁶-1.7x 10⁶Pz
- ◆ Sound Attenuation <= 0.05dB/cm.
- ◆ MHz Viscosity >= 15Pa, S PH Value 5.5-8
- ◆ PH Value 5.5-8

Code No.	Specification	Package	Carton Size	G.W(KG)	N.W(KG)
SUUSG001	15ml/Bottle	60 Bottles/ctn	18*16*13cm	2KG	1.8KG
SUUSG002	250ml/Bottle	60 Bottles/ctn	46*32*20cm	16KG	15KG
SUUSG003	2.5L/Barrel	6 Barrels/ctn	34*32*25cm	17KG	15KG
SUUSG004	5L/Barrel	2 Barrels/ctn	37*20*25cm	11KG	10KG
SUUSG005	1Gallon/Barrel	2 Barrels/ctn	29*26*31cm	10.5KG	9.5KG



Rapid Tests (Code No.: SURPT101)

FERTILITY TESTS	FORMAT
HCG	strip/cassette/midstream
LH	strip/cassette/midstream
FSH	strip/cassette/midstream
TORCH TESTS	FORMAT
TOXO (Toxoplasma)	strip/cassette
CMV(Cytomegalovirus)	strip/cassette
HSV(Herpes Simplex Virus)	strip/cassette
RV(Rubella Virus)	strip/cassette
RESPIRATORY DISEASES TESTS	FORMAT
MP(Mycoplasma Pneumonia)	strip/cassette
CP(Chlamydia Pneumonia)	strip/cassette
TB(Tuberculosis)	strip/cassette
COX B(Coxsackievirus)	strip/cassette
RSV(Respiratory Syncytial Virus)	strip/cassette
Adenovirus	strip/cassette
Mumps	strip/cassette
HPV B19(Human Parvovirus)	strip/cassette
Measles	strip/cassette
FERTILITY TESTS	FORMAT
HAV	strip/cassette
HDV	strip/cassette
HEV	strip/cassette
HGV	strip/cassette
HBsAg	strip/cassette
HBsAb	strip/cassette
HBcAb	strip/cassette
HBV	strip/cassette
HCV	strip/cassette
Malaria	strip/cassette
Dengue	strip/cassette
SEXUALLY TRANSMITTED DISEASES TESTS	FORMAT
HIV	strip/cassette
TP	strip/cassette
NGH(Neisseria Gonorrhoeae)	strip/cassette
CT(Chlamydia trachomatis)	strip/cassette
GASTROINTESTINAL DISEASES TESTS	FORMAT
HP(Helicobacter Pylori)	strip/cassette
Rotavirus	strip/cassette
EV71	strip/cassette
Typhoid	strip/cassette
TUMOR MARKERS TESTS	FORMAT
PSA	strip/cassette
FOB	strip/cassette
CEA	strip/cassette
AFP	strip/cassette
CARDIAC MARKERS	FORMAT
Troponin I	strip/cassette
Myoglobin	strip/cassette
CK-MB	strip/cassette



Aneroid Sphygmomanometer

- ◆ Code No.: SUASM101
- ◆ Material: Gauge: zinc alloy
Bladder: PVC/LATEX
Bulb: PVC/LATEX
- ◆ Measuring range: 0-300mm HG or 20-300mmHG(+3.75mmHG)
- ◆ Packing: 1pc/box, 50pcs/ctn
- ◆ Carton size: 55*37*37cm
- ◆ Color: Black/Blue/Grey etc.



Aneroid Sphygmomanometer (Wall Type)

- ◆ Code No.: SUASM102
- ◆ Square shaped 14.5cm diameter dial
- ◆ Right and left adjustable for easy reading
- ◆ 2 tube adult latex bladder, standard latex bulb
- ◆ Inflation valve with spring, standard and value, PVC coil tubing
- ◆ 1pc/box, 20 boxes/ctn



Stethoscope

- ◆ Code No.: SUSTS103
- ◆ Type: normal or sprague rappaport type
- ◆ Specification: single head/double head, Y tube
- ◆ Color: black, pink, purple, blue
- ◆ Packing: 1pc/box, 50boxes/ctn
- ◆ Carton size: 66*31*43cm



Desk-top Sphygmomanometer

- ◆ Code No.: SUDTSM104
- ◆ Description: Aluminum shell
- ◆ Color: Black,yellow,red,gray,white,green,blue,light green
- ◆ Test scope: 0-300mmHg
- ◆ Precision: +/-3mmHg
- ◆ Mini scale division: 2 mmHg
- ◆ Packing quantity: 1pc/box, 10pcs/ctn
- ◆ Carton size: 30*22*37 cm



Baby Feeding Bottle

- ◆ Code No.: SUBFB105
- ◆ Material: PP or Silicone
- ◆ Type: With or without handle
- ◆ Size: 150ml with/without handle
250ml with/without handle
- ◆ Package: 1pc/box, 144pcs/ctn, 81*41*43cm(250ml)
1pc/box, 144pcs/ctn, 60*41*40cm(150ml)



Baby Pacifier

- ◆ Code No.: SUBBP106
- ◆ Material: High grade PP or silicone
- ◆ Color: Blue, Black, Pink, Purple
- ◆ Size: 5cm
- ◆ Packing: 1pc in one opp, 80pcs/ctn
- ◆ BPA free
- ◆ Steam sterilized



Breast Pump

- ◆ Code No.: SUBPP107
- ◆ Buld material: medical grade PVC, PP or Glass
- ◆ Cup material: medical grade PP or PC
- ◆ Volum: 41ml, 57ml etc
- ◆ Weight: 48g, 58g etc
- ◆ Color: red, blue, yellow, black etc
- ◆ Size: 6.8x13.5cm



Blood Glucose Meter

- ◆ Code No.: SUBGM108
- ◆ Result in only 8 seconds
- ◆ Automatic shut-off in 3 minutes
- ◆ 180 test memory
- ◆ 1.7.14.21.28-days average display
- ◆ Less than 0.8ul blood sample
- ◆ Code strip, no need hand setting
- ◆ Common type or auto code type



Oxymeter

- ◆ Code No.: SUOMT109
- ◆ Suite for children and adult
- ◆ OLED screen,
- ◆ product size: 55*33*28mm
- ◆ Pulse frequency range: 25-250BMP
- ◆ Pulse frequency accuracy: ±1%
- ◆ Blood oxygen range: 701-100%
- ◆ Blood oxygen precision: ±2%



Arm Blood Pressure Monitor

- ◆ Code No.: SUABPM110
- ◆ Product size: 119*107*52mm
- ◆ Display: LCD digital display
- ◆ Measurement method: pulse scanning
- ◆ Test Range: Hypertension/ Hypotension: 20mmHg~280mmHg
- ◆ Pulse: (40~200) times/minute
- ◆ Accuracy: $\pm 3\text{mmHg}(\pm 0.4\text{Kpa})$
- ◆ Pulse: Within $\pm 3\%$ of reading
- ◆ Voltage: 6V
- ◆ Packing: 1set/box, 30sets/ctn
- ◆ Carton size: 50*39*32cm



Wrist Blood Pressure Monitor

- ◆ Code No.: SUWBPM111
- ◆ Product size: 7.56*6.84*3.1cm
- ◆ Measurable wrist circumference: 13.5cm~19.5cm
- ◆ Display: LCD digital display
- ◆ Measurement method: pulse scanning
- ◆ Test Range: Hypertension/ Hypotension: 20mmHg~280mmHg
- ◆ Pulse: (40~200) times/minute
- ◆ Accuracy: $\pm 3\text{mmHg}(\pm 0.4\text{Kpa})$
- ◆ Pulse: Within $\pm 3\%$ of reading
- ◆ Voltage: 2.2V
- ◆ Power consumption: 0.01W
- ◆ Packing: 1set/box, 30sets/ctn
- ◆ Carton size: 50*39*32cm



Thermometer

- ◆ Code No.: SUTMT112
- ◆ Digital thermometer and mercury thermometer
- ◆ Digital thermometer: pen-like digital thermometer, infrared Thermometer, infrared ear thermometer, gun thermometer
- ◆ Mercury thermometer: for oral/axillary
- ◆ Packing for digital thermometer: 1pc/tube, 12pcs/box, 1-
- ◆ Packing for mercury thermometer: 1pc/box, 360pcs/ctn
- ◆ Feature: High sensitive and accurate, easy to use and no allergy, non-contact skin.



Pull Out Vacuum Apparatus

- ◆ Code No.: SUPOVA113
- ◆ Material: ABS Plastic/GPPS polystyrene resin
- ◆ Type: Set 6, Set 12, Set 24 etc
- ◆ Color: Transparent
- ◆ Air Pump: PP/ABS
- ◆ Gun: big/small
- ◆ Handle: Plastic Handle
- ◆ Extension tube: 60cm
- ◆ Air Pump: PP/ABS
- ◆ Function: magnetotherapy cupping therapy
- ◆ Cupping therapy is combining modern technology and traditional medicine, it is China treasure--a new health care product.



Medical Regulator

- ◆ Code No.: SUMRT114
- ◆ Input Pressure: 150bar;
- ◆ Output Pressure: 2bar-3bar;
- ◆ Control Pressure of Safety Valve: 3.5bar $\pm$ 0.5bar;
- ◆ Flowmeter Control Range: 0-10LPM/15LPM;
- ◆ Material: Heavy mainbody brass + aluminum alloy flow meter;
- ◆ Connecting thread Male G5/8";



IV Stand

- ◆ Code No.: SUIVS115
- ◆ Circular Ram Horn" (4 or 6 hooks)
- ◆ 18 inch diameter steel base provide the stand with excellent stability.
- ◆ Regulate the range: 100cm
- Technical Features:**
- ◆ Packing size for stand (6pcs/box, L x W x H): 56 x 56 x 43 cm
- ◆ Packing size for pipe (6pcs/box, L x W x H): 130 x 20 x 10cm
- ◆ N.W.: 15.5kg G.W.: 17.5kg (6 sets)



Crutch

- ◆ Code No.: SUCTC116
- ◆ Material: Aluminum alloy material
- ◆ High adjustable from 93 to 155cm rubber tip assembly.
- ◆ 9-position adjustable.
- ◆ TPR Grip. Five speed adjustable suitable for different people height.
- ◆ Maximum load-bearing 100kg.
- Technical Features:**
- ◆ Packing size: S/M/L 113*26*31CM 123*26*31CM 133*26*31CM
- ◆ 10Pair /ctn



Stretcher

This stretcher with outriggers and can be folded, It is made of high strength aluminum alloy material and waterproof leather materials, It has characterizes of light-weight-ed, small-sized, easy-carried, reliable quality, use-safely and easy for sterilization etc. It is commonly used in hospital, gymnasium, ambulance, and battle fields to carry patients and wounded person.

- ◆ Code No.: SUSTC117
- ◆ Product size(L*W*H): 208*55*13cm
- ◆ Folded size(L*W*H): 104*17*9cm
- ◆ Packing size(2pcs/carton):108*18.5*18cm
- ◆ N.W.: 4.7kg G.W.: 5.7kg
- ◆ Loading bearing: less than 159kg



Ambulance Stretcher

- ◆ Code No.: SUABS118
- ◆ It is made of high-strength aluminum alloy thick-wall pipes
- ◆ The bed surface is made of foamed cushion and adjustable back making the wounded feel comfortable
- ◆ The mechanical foldaway legs are controlled by the handles on both ends of the stretcher
- ◆ And the legs can folded under the stretcher automatically when loading into the vehicle
- ◆ Also the stretcher can be locked by the fixing devices when on the ambulance car
- ◆ This type of stretcher is characterized by its being light, durable antiseptic easy for sterilization



Spine Board

Uses the high strength engineering plastic, durable, not easy to aging, can the safe transfer of the sick and wounded, is a can be equipped to medical departments of emergency transfer equipment, can undertake X light perspective



Code no.	Size	Loading	N.W	G.W	Packing size
SUSPB119	183x45x6cm	≤159KG	7KG	10KG	1pc/ctn 185x46x7cm

Hospital Bed

- ◆ Code No.: SUHPB120
- ◆ Specification: L2150 x W1050 x H460-720mm.
- ◆ Standard Accessory:
- ◆ -Detachable ABS head/foot board
- ◆ -Collapsible ABS guardrails
- ◆ -Integrated stamping frame and mattress base are made of steel epoxy
- ◆ -I.V. Pole & Drainage hook with PU cover
- ◆ -Heavy duty caster with brakes 5"
- ◆ -Hand remote control
- ◆ -Optional accessory:
- ◆ -Mattress & over-bed table & bedside cabinet & Back-up battery & CPR control
- ◆ -Carton size: 2100 x 980 x 500mm
- ◆ -G.W.: 90kg N.W.: 85kg



Transfusion Chair

- ◆ Code No.: SUTFC121
- ◆ Fabric is pressed PVC artificial leather
- ◆ 3 angles adjustable
- ◆ With rotatable table
- ◆ Product size: 90x75x110 cm
- ◆ Packing size: 100x60x44 cm
- ◆ G.W.: 34kg N.W.: 31.5kg



Wheel Chair

- ◆ Code No.: SUWCR122
- ◆ Front wheel: 8 solid universal tires
- ◆ Rear wheel: 24 pneumatic tires
- ◆ T Front fork: Black powder coated
- ◆ Hand rim: PVC
- ◆ Frame: Silver powder coated steel frame
- ◆ Armrests: Fixed, PE
- ◆ Seat cushion: Oxford cloth
- ◆ Pedals: Silver powder coated iron
- ◆ Brake: Powder coated
- ◆ Cross Brace: double, powder coated
- ◆ Guard board: Powder coated iron
- ◆ Packing: 1pc/carton, 135pcs/20GP, 338pcs/40HQ
- ◆ N.W:14KG G.W:16KG
- ◆ Carton size: 90cm*26cm*87cm



Foot Stool Steel

- ◆ Code No.: SUFSS123
- ◆ Frame Material: 201 or 304 Stainless Steel
- ◆ Pedal Surface Material: TPR
- ◆ Stable structure and durable material
- ◆ With four anti-slipping rubber feet
- ◆ Product size: 490*290*205mm
- ◆ Safe bearing: 200kg
- ◆ Packing: 620*430*70mm



Otoscope & Ophthalmoscope Set

- ◆ Code No.: SUOOS124
- ◆ Specification: VM-T200 Otoscope & VM-P100 Ophthalmoscope
- ◆ Technology: LED/Without Optical fiber/3x Magnification/small view window
- ◆ Weight: 7.5kg
- ◆ Quantity: 20pcs
- ◆ Dimension: 40*24*23cm



Height Adjustable Stool

- ◆ Code No.: SUHAS125
- ◆ Electroplating steel frame
- ◆ Equipped with 5 castors
- ◆ The height can be adjusted by gas spring
- ◆ Size: 32*32*(46-59)cm
- ◆ Packing size: 42*42*100mm, 2pcs/ctn
- ◆ GW: 6kg, NW: 5kg



Medical Dressing Table

- ◆ Code No.: SUMDT126
- ◆ Product size: 190*60*75cm
- ◆ Bed frame made of stainless steel rectangular tube,
- ◆ Bed surface is consisted of six same pieces of stainless steel sheet, and one piece of them with a square basin.
- ◆ Load bearing: 150kg
- ◆ Knock down construction
- ◆ Packing size: 205*67*24cm, 2pcs/ctn
- ◆ Weight: 54KG



Examination Couch

- ◆ Code No.: SUEMC127
- ◆ Made of 20*20*1.5mm stainless steel square tube, 5cm
- ◆ Sponge surface with antibacterial PVC wrapped



Massage Bed

- ◆ Code No.: SUMSB128
- ◆ Description: 2-section wooden massage table
- ◆ Dimension: L186*W70*H50~78cm(not include headrest and armrest)
- ◆ Heavy Loading: 250kg
- ◆ Packing: 1pc/ctn
- ◆ Carton Dimension: 95*20.5*74cm
- ◆ Net weight: 15kgs
- ◆ Gross weight: 18kgs
- ◆ Quantity Per Container: 204pcs/20GP 480pcs/40HQ
- ◆ Legs system: German beech wood
- ◆ Material: PU leather, sponge of 5 cm thickness/high quality hardwood frame
- ◆ Upholstery Colour: available to any colors



Air Mattress

Air mattress with pump designed for the treatment and prevention of pressure sore stage in short-term and domiciliary therapies. The mattress in PVC is composed of 1307 cm-high bubble cells, particularly comfortable, it is fixed to the bed by extra flaps on both top and bottom side. Very easy to install, maintain and use thanks to the possibility of adjust comfort range. The pump can be hung to the bed by means of two hooks.

- ◆ Code No.: SUAMT129
- ◆ Dimensions without flaps 200*90*7cm
- ◆ Extension flaps length 50cm*50cm(head&foot side)
- ◆ Mattress material TwNan Ya EN71K80 Medical-Grade PVC
- ◆ Material thickness 0.32mm
- ◆ Cold-resistant -30℃
- ◆ Number of cells 130
- ◆ Weight support 135kg
- ◆ HF Welding one-time forming



X-ray Machine

- ◆ Code No.: SUXRM130
- ◆ Tube voltage: 70KV±10%
- ◆ Tube current: 8mA±20%
- ◆ Exposal time: 0.2~4.0s
- ◆ Input power: <600VA
- ◆ Output power: P=339VA
- ◆ Irradiation diameter: 6cm
- ◆ Fuse: RF1-20-5A
- ◆ Capacity: 70KV.8mA,4s
- ◆ The way of operation: brake load; continuous operation
- ◆ The weight of x-ray source: 7.5Kg
- ◆ The variation rate of high-voltage generator 1:226
- ◆ X-ray area: 3.14x302(mm²)
- ◆ The time of loading follows R'10
- ◆ Up-down range of X-ray head: 900mm
- ◆ Level turning angle of head: 360°

- ◆ Rounding angle of head: 270°
- ◆ Moving level angle: ±270°
- ◆ Life: 6years



Health scale/body scale

- ◆ Code No.: SUHSC131
- ◆ Type: mechanical or digital
- ◆ Rated load: 160KG/350LB
- ◆ Product division: 1kg(≤ 60kg), 1.5kg(> 60kg)
- ◆ Product size: 42.5x28.5x8.5cm, 45.8x34x10.2cm etc
- ◆ Warranty: one year quality warranty
- ◆ Packing: 3pcs/carton carton size: 49x35x33cm
- ◆ N.W/G.W: 11/12kg



Vertical Pressure Steam Sterilizer

- ◆ Code No.: SUVPS132

Model technical data	LS-35LD(automatic)	LS-50LD(automatic)	LS-75LD(automatic)	LS-100LD(automatic)
Chamber volume	35L(φ318X450)mm	50L(φ340X550)mm	75L(φ400X600)mm	100L(φ400X600)mm
Working pressure	0.22MPa			
Working temperature	134°C			
Max working pressure	0.23Mpa			
Heat average	±1°C			
Timer	0~99min or 0~99hour59min			
Adjustment of temperature	105~134			
Power	2.5kw/AC220V.50HZ	3kw/AC220V.50HZ	4.5KW AC220.50HZ	
Overall dimension	480x460x850(mm)	520x520x980(mm)	560x560x980(mm)	590x490x1080(mm)
Transport dimension	570x550x970(mm)	590x590x1110(mm)	650x630x1150(mm)	680x650x1220(mm)
G.W/N.W	56kg/42kg	68kg/50kg	90kg/70kg	105kg/85kg



Electric Suction Apparatus

- ◆ Code No.: SUESA133
- ◆ High vacuum, high flow
- ◆ Input power: 180VA
- ◆ Power supply: 220V,50HZ
- ◆ Max vacuum: ≥80kpa
- ◆ Adjustable vacuum range: 20kpa~Max vacuum
- ◆ Max airflow: ≥20L/Min(760mmHg)
- ◆ Sound Level: ≤60db(A)
- ◆ Suction bottle(glass): 2500ml/bottle, 2 bottles in a group
- ◆ N.W.: 14KG
- ◆ Dimension: 360*320*480(mm)



Oxygen Concentrator

- ◆ Code No.: SUOCT134
- ◆ Model: 5L/min single flow/dual flow
- ◆ Features:
 - ◆ PSA technology
 - ◆ Adjustable flow rate
 - ◆ Dehumidity function
- ◆ Timing: LCD show record the accumulated working time of the machine, accumulated time, set time
- ◆ Power failure alarm
- ◆ High and low pressure alarm
- ◆ Flow rate: 0-5LPM
- ◆ Purity: 93%±3%
- ◆ Outlet pressure(Mpa): 0.04-0.07
- ◆ Sound level(dB)≤50
- ◆ Power: AC220V±10%, 50/60±1Hz; AC110V±10%, 50/60±1Hz
- ◆ Powder consumption: ≤400W
- ◆ LCD display: switch times, operating pressure, present working time, accumulating time, presetting time from 10mins to 40hours



TD-4Z Benchtop low Speed Centrifuge

- ◆ Code no.: SUBLSC135
- ◆ Features & Advantages:
 - ◆ Variable frequency drive motor, digital display
 - ◆ All-steel body, stainless steel chamber
 - ◆ Electronic lid lock to ensure safety
 - ◆ 9 levels of acceleration and deceleration
- ◆ Aluminum alloy Rotor
- ◆ RCF can be set up directly
- ◆ Time, Speed, RCF, etc. can be modified during operation

Max speed	4200r/min
Max RCF	2680 xg
Max capacity	12x20ml or 4x50ml (Fixed-angle Rotor)
Speed Accuracy	±20r/min
Timer Range	1min~99min
Noise	≤55dB (A)
Power Supply	AC 220V 50Hz 5A
Dimension	440x340x290 (LxWxH)mm
Weight	23kg
Power	200W



ABS Medicine Trolley

- ◆ Code No.: SUAMT136
- ◆ Size: 625x475x930mm
- ◆ ABS table top with transparent dust-proof pad, convex design to avoid down, SUS304 guardrail
- ◆ With dust baskets, utility container
- ◆ 2 small drawers, 2 middle drawers, 1 big drawer-centralized lock
- ◆ Multi-bin container
- ◆ Inner with partitions can be well organized
- ◆ Luxurious Noiseless casters with brake
- ◆ Plastic-steel columns



Medical Privacy Screen

- ◆ Code No.: SUMPS137
- ◆ Patient screen on castors is mainly used for isolation and protecting patient privacy during examination in hospital
- ◆ Descriptions:
- ◆ Size: 2000x1870mm (4pieces)
- ◆ Epoxy powder coated steel frame
- ◆ Broad space casters design is for stability
- ◆ Each panel of screen is with a hinge, which can be turn into any angle
- ◆ Noiseless castors 2 inch
- ◆ Carton Size: 1850x520x120mm (1pc/ctn)



Binocular Microscope

- ◆ Code no.: SUBMC138
- ◆ Total Magnification: 40-1600x
- ◆ Eyepiece: wide field plane-scope eyepiece WF10X, WF16X
- ◆ Achromatic; 45mm DIN achromatic objective: 4X, 10X, 40X(S), 100X(S.Oil)
- ◆ Eyepiece tube: Articulated free binocular head, interpupillary distance: 55-75mm, 45° inclined, 360° rotatable
- ◆ Illumination: Built-in adjustable brightness halogen lamp 6V/20W and kohler illumination system
- ◆ Diaphragm: 2-30mm iris diaphragm and 32 filter
- ◆ Condenser: ABBE condensor.N.A.1.25
- ◆ Adjustable Mechanism: coaxial coarse and fine focusing adjustable mechanism: 20mm, precision: 0.002mm
- ◆ Working stage: double layer mechanical stage: 130x140mm, moving range X-Y: 70x50mm
- ◆ Net weight: 5.6kg
- ◆ Packing: 2sets/ctn, 57x37x49cm



Antibacterial Gel 500ml

- ◆ Code No.: SUABG139
- ◆ Icerik/Ingredients: Deionized water, Ethyl alcohol, Acrylates/ C-10-C30 Alkyl Acrylate Crosspolymer, Tri Ethanol Amine, Glycerin, Fragrance.
- ◆ Properties: It is used for having and keeping hygiene.
- ◆ Directions: Place enough product in your palm to thoroughly cover your hands. Rub hands together briskly until dry. No rinsing required, non towel needed.
- ◆ Warning: Flammable,keep away from fire of flame. For external use only. Keep out of eyes. In case of eye contact,immediately flush eyes with water.If swallowed call a doctor. Do not use directly to infected or open wounds. Do not swallow. Keep out of children. Keep in dry and away from direct sunlight, under 30°C



2% Cidex Solution

- ◆ Code No.: SUCST140

Active ingredient	Glutaraldehyde
Active content	2.0%—2.4%
Dosage form	Liquid
Brand	Retouch
Feature	Disinfection and sterilization of medical devices
Method of use manual	Automatic cleaning machine low-flavor and low-foam type
Time of action	1.Gastroscope,colonoscopy,duodenoscopy should be immersed for not less than 10 minutes; 2.immersion in bronchoscopy for not less than 20 minutes; 3.After infection of patients with special infections such as Mycobacterium tuberculosis and other mycobacteria, The endoscope is immersed for not less than 45 minutes.
Easy to use	Preservatives and activators are added before use and up to 14 days after activation.
Toxic rationality	This product contains high levels of active ingredients, strong volatility, strong odours and common allergies.



Povidone Iodine Solution

- ◆ Code no.: SUPLS141
- ◆ Dosage form Disinfectant
- ◆ Appearance: Red brown liquid
- ◆ Indication: Used for skin and mucosa disinfection, surgical equipment disinfection, aquatic animal body surface sterilization or fish eggs disinfection.
- ◆ Precaution: Don't use with alkalinity and QAC meanwhile.
- ◆ Storage: Avoid light, seal well and keep in dark place.

ANTIMICROBIAL		
No	Generic name	Dosage
1	Povidone Iodine	10%/200ml/bottle 10%/100ml/bottle 5%/100ml/bottle



Piston Nebulizer

- ◆ Code No.: SUPNZ142
- ◆ Type: Vibration Mesh Nebulizer
- ◆ Model: RJ-101A
- ◆ Size: 68mm (L) x 52mm (W) x 125 mm (H)
- ◆ Power: 3V DC (Two "AA" 1.5V Batteries/ optional AC adapter)
- ◆ Vibrating Frequency: Approx. 110kHz
- ◆ Nebulizer Speed: 0.2-0.7mL/min
- ◆ Mesh /Particle Size: 3 μm
- ◆ Medication capacity: 8 ml Maximum / 0.5 ml Minimum
- ◆ Battery Life: Up to 150 minutes used continuously with 2AA batteries
- ◆ Packing: 1pc/box, 20boxes/ctn
- ◆ Carton size: 68*37.5*33cm



Breathing Exerciser

- ◆ Code No.: SUBTE143
- ◆ Model: 5000ml
- ◆ Material: PP/PS or PE
- ◆ Product size: 165*75*170mm
- ◆ Packing: 1pc/box, 50pcs/ctn
- ◆ Carton size: 87*42*37cm



Micro Applicator

- ◆ Code No.: SUMAC144
- ◆ Applicator Rod Series: 1.0mm/1.5mm/2.0mm/3.0mm
- ◆ Material: Plastic
- ◆ Colour: Pink, red, blue, yellow, white etc.
- ◆ Packing: 100pcs/bottle, 4bottles/box, 10boxes/ctn
- ◆ Carton size: 49*40*26cm



Pill Counting Tray

- ◆ 1. Code no.: SUPCT145
- ◆ 2. Material: ABS+PP
- ◆ 3. Color: blue or customized
- ◆ 4. Product size: 20*14.5*3.8 cm



Back Stretcher

- ◆ Code No.: SUBST146
- ◆ Material: PP, ABS, NBR
- ◆ Color: Black with foam pad in blue, green, purple, red, pink
- ◆ Function: Relieve pains, Pains relief, Posture corrector
- ◆ Feature: Three gear adjustment



Enfant Warmer

- ◆ Code No.: SUEWM147
- ◆ Skin temp and manual controlled by micro-computer
- ◆ Over temperature alarm
- ◆ Available adjustment or radiator and mattress tilting
- ◆ With one illumination light
- ◆ APGAR timer
- ◆ RS-232 connector
- ◆ Packing size for pipe (1pc/ctn, L*W*H): 117*76*118cm



Portable Type Stainless Steel Sterilizer

- ◆ Code No.: SUPTSS148
- ◆ Chamber Volume: 18L (280*260mm)
- ◆ Working Pressure: 0.14-0.16MPa
- ◆ Working Temperature: 126°C
- ◆ Double scale indication pressure gauge. Electric heated
- ◆ Easy to operate, safe and reliable
- ◆ Carton size: 410*410*520mm



Mayo Table

Code No.: SUMTB149

This mayo trolley is made of stainless steel 201 material
The thickness of frame is 1mm, the thickness of tray is 1.2mm

- ◆ Product size (L*W*H): 60*40*82cm
- ◆ Packing size (10pcs/carton): 121*61*71cm
- ◆ N.W. : 40kg G.W. :50kg



Accompanying Chair

Code No.: SUAPC150

The chair is made of pressed PVC artificial leather fabric,
and powder coated steel tubes.

- ◆ Product size (L*W*H): 60*40*82cm
- ◆ Packing size (10pcs/carton): 121*61*71cm
- ◆ N.W. : 40kg G.W. :50kg



Baby Incubator

- ◆ Code No.: SUBIB151
- ◆ Microprocessor based servo controlled temperature
- ◆ System Control modes: Air mode
- ◆ Humidity is adjustable in two levels
- ◆ Set temperature, air temperature, heating power are displayed separately by LED
- ◆ Self - testing function, various failure alarms by audible and visual
- ◆ Alarm function
- ◆ Power failure, Over temperature, Temperature deviation, Fan failure, Temperature
- ◆ Sensor failure > 37° temperature set function
- ◆ Triple protection for over temperature with separate cut off device, more safety system.
- ◆ The inclination of infant bed is adjustable
- ◆ Single wall hood, 4 operating windows and 2 iris ports, 3 infusion seals locking device for front door
- ◆ RS232 connector, 1 tray
- ◆ Carton size: L 126*W 71*H 86cm



Dental Chair

- ◆ Code No.: SUDTC152
- ◆ All controlled by the electric valve
- ◆ Plastic pipe ◆ DC motor ◆ Air spring
- ◆ Air-lock balance arm system ◆ Handpiece tubes
- ◆ 3 way Syringe(one for hot and one for cold),2pcs
- ◆ Water Suction and Saliva Ejector,1set each
- ◆ The device for prevent reflux and the Built-in Water purified System
- ◆ Turnable Cuspidor ◆ Automatic Heating water system
- ◆ Operation Light ◆ Low Voltage DC X-film viewer
- ◆ Display System of the H.P.'s air pressure
- ◆ Dentist Stool,1pc(No.2 Old style)
- ◆ With low mounted instrument tray
- ◆ Carton size:155*105*125cm,1pc/ctn



Multi-functional chair

- ◆ Code No.: SUMFC153
- ◆ Length: 180cm, Width: 55cm
- ◆ Backrest: 78cm, Seat:52cm
- ◆ Legrest: 35+15cm, Armrest: 54cm
- ◆ Max. Height: 82cm, Min.Height: 56cm
- ◆ Max. Load: 250kgs; Packing size: 168*80*72cm
- ◆ Net weight: 78kgs; Gross Weight: 90kgs
- ◆ Motors: 4 or 5 Motors(CE), Leather: PVC/Fiber
- ◆ Upholstery: High density Foam, Base: Steel Structure reinforced plastic casting
- ◆ Control: Handset Control/Foot control optional
- ◆ Carton size: 168*80*72cm, 1pc/ctn



Diagnosis and treatment table

- ◆ Code No.: SUDTT154
- ◆ Dimension: 755*525*795mm
- ◆ Desk Material: Marble
- ◆ Power:AC 220V,50/60HZ; Power: 1000VA
- ◆ Positive pressure: 2.5 Kg/cm2
- ◆ Negative Pressure: 740mmHg(Max)
- ◆ Spray gun: 3pcs;Suction gun:1pc
- ◆ Laryngoscope Preheater: 1pc
- ◆ Lead spotlight: 1pc
- ◆ x-ray film viewer: 1pc
- ◆ Carton size:87*85*95cm,1pc/ctn



Hospital Laundry Waste Trolley

- ◆ Code No.: SUHLWT155
- ◆ Specification: 800*590*890mm
- ◆ Inner diameter of the cloth bag: 250*400*670mm
- ◆ Carton size: 890*570*980mm



Electric-heating Water Distiller

Code No.: SUEHWD156
Stainless steel material, adopt stamping and welding technology.
Corrosion and ageing resistance, easy to operate, safe and durable.

Model	WD-5	WD-10	WD-20
Water Output	=5 L/H	=10 L/H	=20 L/H
Consumption	5KW	7.5KW	15KW
Power Supply	AC110/220V, 50/60Hz	380V/50Hz	380V/50Hz
External Size (Dia.*H)mm	310*700mm	355*820mm	405*677mm
Package size(W*D*H)	360*360*820mm	360*350*900mm	540*390*740mm
Gross weight	7.5kg	8.5kg	11.5kg



Auto-control Electric-heating Water Distiller

Code No.: SUAEDW157
Stainless steel material, adopt stamping and welding technology.
Water shortage and power-off alarm; automatically make up water and heat again.
Power, heating indicator and show water shortage.

Model	WD-A5	WD-A10	WD-A20
Water Output	=5 L/H	=10 L/H	=20 L/H
Consumption	5KW	7.5KW	15KW
Power Supply	AC110/220V, 50/60Hz	380V/50Hz	380V/50Hz
External Size (Dia.*H)	310*740mm	235*855mm	415*720mm
Package Size(W*D*H)	360*360*830mm	360*360*920mm	535*380*770mm
Gross Weight	8kg	9.5kg	13kg



SOLUTION FOR INFUSION Glass bottle:-		
No	Generic name	Dosage
1	Ciprofloxacin 200 mg Solution for Infusion	100ml/bottle
2	Metronidazole 500mg Solution for Infusion	100ml/bottle
3	Paracetamol 1gm Solution for Infusion	100ml/bottle
4	Fluconazole 200mg Solution for Infusion	100ml/bottle

INTRAVENOUS INFUSION Glass bottle / PP bottle		
No	Generic name	Dosage
1	9% Sodium Chloride Solution for Infusion	500ml/bottle
2	9% Sodium Chloride & 5% Glucose Solution for Infusion	500ml/bottle
3	Dextrose 5% Solution for Infusion	500ml/bottle
4	Ringer Lactate Solution for Infusion	500ml/bottle
5	Sterile water for Injection	5ml-10ml

ANTIBIOTICS		
No	Generic name	Dosage
1	Azithromycin Tablets	250mg*100 tablets/box

ANTIMALARIAL		
No	Generic name	Dosage
1	Artemether + lumefantrine Tablets	20mg+120mg*24 tablets/box 40mg+240mg*12 tablets/box 80mg+480mg*6 tablets/box
2	Artesunate Tablets	50mg*10 tablets/box
3	Artesunate Injection	60mg /IM -IV
4	Quinine Injection	600mg/2ml*10ampoules/box 400mg/4ml*100ampoules/box

