



180011112242



(2018) 国认监认字 (244) 号



中国认可
国际互认
检测
TESTING
CNAS L6780

检 验 报 告

TEST REPORT



Scan is following us

Scan checks for authenticity



* W 2 0 2 0 1 0 8 8 5 E *

浙江省轻工业产品质量检验研究院

Zhejiang Light Industrial Products Inspection and Research Institute

国家纺织服装产品质量监督检验中心（浙江）

National Textiles and Garment Quality Supervision Inspection Center(Zhejiang)

Test report

Number:W202010885E

page1/3

Name of Customer	XIAMEN MIAO XING TECHNOLOGY CO.,LTD.	Address	Floor 5, WORKSHOP 4, 333 TONGFU ROAD, TONGAN DISTRICT, XIAMENG CITY
Manufacturer	XIAMEN MIAO XING TECHNOLOGY CO.,LTD.	Address	Floor 5, WORKSHOP 4, 333 TONGFU ROAD, TONGAN DISTRICT, XIAMENG CITY
Sample information	Name of sample: KN95 Face Masks Characteristics of sample: White Trademark of sample: --- Specification/model: --- Level: KN95 Category of safety specification: --- Art. No.: --- -----		
Above-mentioned information by Customer-supplied			
The sent way of sample	Courier	Sample quantity	30 pieces
Receiving Date of Sample	2020/04/25	Test Category	Entrusted inspection
Date of Testing	2020-04-26- 2020-04-29		
RatingRequirements	GB 2626-2006		
Test Summary: See the attached page for the results. <div style="text-align: right;"> [签章] Test Seal Date of Approval: 2020-04-29 </div>			
Remarks	The sample is not pretreated and the result is for reference.		

Approved by:



Test report

Number:W202010885E

page2/3

ITEM	STANDARD (KN95)	RESULT	RATING
1. FILTER EFFICIENCY(SALT MEDIUM) (GB 2626-2006 6.3)(%)			
-	≥ 95.0	99.3	PASS
2. RESPIRATORY RESISTANCE (GB 2626-2006 6.5 & 6.6)(Pa)			
-INSPIRATORY RESISTANCE	≤ 350	124	PASS
-EXPIRATORY RESISTANCE	≤ 250	110	PASS
3. HEADBAND (GB 2626-2006 6.11)			
-	$\geq 10N$,hold on10s	No breakage and slippage were observed	PASS
4. APPEARANCE (GB 2626-2006 6.1 & 6.2)			
-	Meet the requirements of GB 2626-2006 standard	The appearance quality of 2 masks conforms to the standard	PASS
5.VISUAL FIELD BELOW MASK (GB 2890-2009)(°)			
-	≥ 60	80	PASS

Test report

Number:W202010885E

page3/3

Picture(s) of sample



—End of report—



DECLARATION

1. Our organization guarantees impartiality, independence and honesty of inspection, and is responsible for the results of inspection, keeping the samples supplied by the entrusting party confidential and at the same time protecting the ownership of the samples supplied.
2. The test report will be deemed invalid without signatures of the inspector/reviewer and authorized personnel, and the red special inspection stamp of our organization.
3. The test report will be invalid if it is altered. Copies of the report are invalid without the red special inspection stamp of our organization.
4. The test results shown in this report are applicable only to the samples provided by customers.
5. All the pages of the report are integral parts of the report. Our organization will not be responsible for any misunderstanding or other results caused by using separate page(s) of the report.
6. If there is any dissent of the report, the entrusting party shall notify our organization timely. For the mandatory inspection given by governmental administration departments, and dissent about the sample being tested or test results on the report should be dealt with in accordance with national regulations.

浙江省轻工业品质量检验研究院及附设的检验中心 The Affiliated Inspection Centers

浙江省轻工业品质量检验研究院
Zhejiang Light Industrial Products Inspection and Research Institute
国家纺织服装产品质量监督检验中心（浙江）
National Textiles and Garment Quality Supervision Inspection Center (Zhejiang)
地址：浙江省杭州市江干区下沙路 300 号 6 号楼
Address: Building No.6, 300 XiaSha Road, Hangzhou, Zhejiang

联系电话：0571-85122669
Telephone: 0571-85122669

E-Mail: 1137907994@qq.com

国家家具产品质量监督检验中心（浙江）
National Center for Quality Supervision Inspection of Furniture (Zhejiang)
浙江省室内安全及家具产品质量检验中心
Zhejiang Center of Quality Test for Indoor Safety and Furniture Products
地址：浙江省杭州市余杭区良渚街道经一路 1 号良渚大学科技园 4 号楼
Address: Building 4 LiangZhu University Science and Technology Park, No.1,
JingYi Rd, Yuhang District, Hangzhou, Zhejiang

联系电话：0571-89009556
Telephone: 0571-89009556

E-Mail: 2047699564@qq.com

浙江省轻工及五金产品质量检验中心
Zhejiang Center of Quality Test for Light Industry and Hardware Products
浙江省体育用品质量检验中心
Zhejiang Center of Quality Test for Sports Products
地址：浙江省杭州市西湖区天目山路 222 号 3 号楼
浙江省杭州市余杭区良渚街道经一路 1 号良渚大学科技园 3 号楼
Address: No. 222 Tianmushan Rd., Hangzhou, Zhejiang
Building 3 LiangZhu University Science and Technology Park, No.1,
JingYi Rd., Hangzhou, Zhejiang

联系电话：0571-89001107
Telephone: 0571-89001107

E-Mail: wujtest@126.com

国家锁具产品质量监督检验中心（浙江）
National Center for Quality Supervision Inspection of Lock (Zhejiang)
浙江省锁具产品质量检验中心
Zhejiang Center of Quality Test for Lock Products
地址：浙江省杭州市塘苗路 24 号
Address: No. 24 Tangmiao Rd., Hangzhou, Zhejiang

联系电话：0571-85027738
Telephone: 0571-85027738

E-Mail: Locktest@126.com

网上业务受理/报告查询：<http://www.zjttj.cn:807>
Online Business-Reception/Report Inquires: <http://www.zjttj.cn:807>

投诉电话：0571-85023552
Complaint Tel: 0571-85023552